



# CARDIOVASCULAR & CRITICAL CARE SOLUTIONS

## PRODUCT CATALOGUE



CARDIOVASCULAR  
& CRITICAL CARE

# TABLE OF CONTENTS

4	<b>PRESSURE MONITORING</b>
31	<b>THERMODILUTION &amp; PULMONARY ARTERY CATHETERS</b>
35	<b>CENTRAL VENOUS CATHETERS</b>
45	<b>CATHETERS AND HAEMOSTASIS</b>
49	<b>FLUID MANAGEMENT</b>
65	<b>TUBING</b>
68	<b>SAFETY &amp; INFECTION CONTROL</b>
83	<b>COMPLETE PROCEDURAL SOLUTIONS</b>
91	<b>ACCESSORIES</b>

# IT ALL BEGAN with the coronary control syringe. The rest is history.

In 1987, Merit Medical pioneered a clear-as-glass polycarbonate syringe to reduce the risk of a glass syringe breaking in the hand of a healthcare provider. Since our beginning, we've consistently prioritised solutions to your procedural challenges, with the goal of making your job safer, more accurate and more efficient. We believe that when you have the right tools to do your job well, it results in better patient care.

Merit Medical's comprehensive Cardiovascular and Critical Care (CCC) product portfolio, presents these healthcare solutions. This includes custom kits that decrease set-up time and reduce waste, as well as products like blood pressure monitoring systems for the ICU and CICU that enhance the accuracy of patient care. Our products like the Futura® Safety Scalpels and Merit Advance® Safety Needles aim to prevent workplace injury, while our infection control solutions like the BackStop+™ Disposal System, the DualCap® Disinfection System, and the Safedraw® Closed Arterial Blood Sampling Kits reduce contaminant exposure for both you and your patients.

In addition to the benefits each of these products provide, partnering with Merit Medical will ensure that healthcare providers can rely on product consistency across the entire continuum of patient care.



Coronary Control  
Syringe  
CCS100



# PRESSURE MONITORING

---

## **MERITRANS DTXPLUS®, SAFEDRAW® AND RT2000™ TRANSDUCERS & TRANSDUCER SETS**

- 5 Meritans DTXPlus® Disposable Pressure Transducer Set Overview**
- 6 Meritans DTXPlus® Disposable Pressure Transducers**
- 7 Meritans DTXPlus® Gabarith™ Tested Disposable Transducer Sets**
- 8 Safedraw® Safe Arterial Sampling Disposable Transducer Sets**
- 9 Meritans DTXPlus® & Safedraw® Flexi Disposable Transducer Sets**
- 10 Safedraw®-P Safe Arterial Sampling Disposable Transducer Sets**
- 12 RT2000™ Reusable Transducers and Disposable Domes & Dome Systems**
- 14 Meritans DTXPlus® and Safedraw® Paediatric & Neonatal Transducer Sets**
- 15 Meritans DTXPlus® Cardiac Catheterisation Lab Transducer Sets**
- 16 Meritans DTXPlus® Accessories**
- 17 Meritans DTXPlus® Transducer Interface Cables**

## **MERITRANS® TRANSDUCERS, TRANSDUCER SETS & TRAM®/TRAM-P® MANIFOLDS WITH INTEGRATED TRANSDUCER**

- 22 Meritans® Disposable Transducer Sets**
- 26 Tram® Manifolds with Integrated Transducer**
- 27 Tram-P™ Manifolds with Integrated Transducer**
- 28 Meritans® Transducer Accessories**
- 29 Meritans® Interface Cables**
  
- 30 Merit Mentor™**

## Meritrans DTXPlus® Disposable Pressure Transducer Set Overview


Designed for ease of use and minimal wavelength distortion

Haemodynamic monitoring information is used to manage the care of critically ill patients. It is essential that one is confident about the accuracy of the information generated.

Accuracy in pressure monitoring is not just dependent on the transducer, it is dependent on all the components in the pressure monitoring system; catheter, pressure tubing, monitor filter and stopcocks - the "total chain." The Gabarith™ test is a scientific technique that provides information about the dynamic response performance of the total chain pressure monitoring system.

Meritrans DTXPlus® pressure monitoring sets confirmed by the Gabarith test have a known level of accuracy. The Gabarith™ test exceeds the demands of the AAMI standard (Association for the Advancement of Medical Instrumentation).

Billiet E, Colardyn F. (1998) Pressure Measurement Evaluation and Accuracy: The Gabarith Test. Intensive Care Medicine 24:1323-1326



### MERITRANS DTXPLUS® DISPOSABLE PRESSURE TRANSDUCER SETS SPECIFICATIONS

DETAIL	DTXPLUS® SPECIFICATIONS
<b>Environmental performance:</b> a) storage temperatures b) operating temperatures c) humidity d) atmospheric pressure	-25°C to 70°C 15°C to 40°C 10% to 90% 425 to 850 torr
<b>Pressure range:</b> a) operating pressure b) overpressure withstand	-30 to 300 mmHg -700 to 7,755 mmHg
<b>Mounting requirement:</b>	Will meet performance requirements when mounted on any axis
<b>Frequency response (Stand-alone Transducer):</b>	15% bandwidth shall be >200 Hz
<b>Transducer excitation:</b>	DC to 5,000 Hz 4 to 8 V RMS
<b>Phase shift:</b>	<5°
<b>Transducer excitation impedance (input impedance):</b>	>200 ohms
<b>Transducer signal impedance (output impedance):</b>	<3,000 ohms
<b>Transducer symmetry:</b>	±5%
<b>Nominal sensitivity (gain):</b>	5 μV/V/mmHg
<b>Unbalance (zero balance):</b>	±40 mmHg
<b>Accuracy:</b>	The maximum error due to the total effects of nonlinearity, hysteresis and sensitivity variations shall be no more than 2% of the reading or ±1 mmHg, whichever is greater, over the operating range.
<b>Risk current:</b>	<2μA
<b>Defibrillator withstand:</b>	Withstand five repeated discharges in 5 minutes of 400J sinusoidal energy delivered into 50 ohms when all transducer wires are shorted and saline is in fluid column
<b>Requirements for the cable:</b> a) cable length stated b) transducer venting to atm. c) cable will withstand five repeated discharges in 5 minutes of 360J sinusoidal energy delivered into 50 ohms when all wires are shorted and foil surrounds the cable	Cable length stated Venting to atmosphere Withstand five repeated discharges in 5 minutes of 400J sinusoidal energy delivered into 50 ohms when all wires are shorted and foil surrounds the cable
<b>Flow Rate, Meritrans DTXPlus® transducers:</b>	2 – 4 cc/hr or 20 – 40cc/hr, depending on model, with an infusion pressure (difference) of 200mmHg.
<b>Fast Flush Rate, Meritrans DTXPlus® transducers:</b>	Min 50cc/min using a microdrip administration set with 60 inch of IV tubing with an infusion pressure (difference) of 200mmHg.
<b>Meritrans DTXPlus® Transducer dimensions (inc Flush and Zero stopcock):</b>	Max Length: 8.4cm (Inc zero stopcock and male luer fitting); Max Width: 3.2cm; Max Height: 3cm (Inc Easy Vent® Deadender on Zero stopcock). Length of Transducer electric cable: 34.5cm

## DISPOSABLE PRESSURE TRANSDUCERS



## Meritans DTXPlus®

Disposable Pressure Transducers  
Sterile, single use

PRODUCT CODE	DESCRIPTION	BOX/CASE CONTENT (UNITS)
682018	Disposable transducer with Critiflex™ integral flush device and stopcock (3cc/hr nominal flow rate)	10/100
682019	Disposable transducer for neonatal/paediatric uses, used in conjunction with a mechanical infusion pump (max. 30cc/hr flow rate)	10/100
682020	Disposable transducers with integral flush device without stopcock (3cc/hr nominal flow rate)	10/100
682021	Disposable transducer without integral flush device or stopcock, suitable for cardiology/neurology applications	10/100

## STAND ALONE ADMINISTRATION SETS

PRODUCT CODE	DESCRIPTION	BOX/CASE CONTENT (UNITS)
682209	Single administration set with 15µ filter, angled cannula and roller clamp (red)	5/25
682210	Double administration set with 15µ filter, angled cannula and roller clamps (red/blue)	5/25
682211	Triple administration set with 15µ filter, angled cannula and roller clamps (red/blue/white)	5/25



## Meritrans DTXPlus® Gabarith™-Tested Disposable Transducer Sets

Sterile, single use



### DISPOSABLE TRANSDUCER SETS FOR THE INVASIVE MONITORING OF HAEMODYNAMIC PRESSURE

PRODUCT CODE	NO. OF DISPOSABLE TRANSDUCERS	FLUSH/FLOW RATE mL/HR	TOTAL PRESSURE LINE LENGTH	PRESSURE LINE(S) COLOUR	ADDITIONAL INFORMATION	BOX/CASE CONTENT (UNITS)	
682001	1		3	150cm	Clear	With R.O.S.E.	5/50
682002	1		3	150cm	Clear		5/50
682004/689006*	1		3	150cm	Red		5/50
682005	1		3	150cm	Blue	With R.O.S.E.	5/50
682006/689007*	1		3	150cm	Blue		5/50
682008	1		3	150cm	Yellow	Add-on set to set no. 682037, -38, -55	5/50
682009	1		3	175cm	Red	With R.O.S.E.	5/20
682010	1		3	175cm	Red		5/50
682011	1		3	25cm	Clear	Velcro strap for arm mount	5/50
682017	1		3	150cm	Clear	With non-shunting stopcock and extension set	5/25
682022	1		30	160cm small bore	Clear	For use with infusion pump	5/50
682025	1		30	35cm small bore	Clear	For use with infusion pump	5/50
682030	2		3	150/150cm	Red/Blue	With R.O.S.E. on blue line	5/25
682032	2		3	150/150cm	Red/Blue		5/25
682037	2		3	150/150cm	Red/Blue	With 3rd IV line for add-on set	5/25
682041	2		3	150/150cm	Blue/Yellow		5/25
682046	3		3	150/150/150cm	Red/Blue/Yellow		5/25

#### PRODUCT DESCRIPTION

- All sets have been verified by the unique Gabarith™ accuracy test
- A range of Single, Double and Triple Disposable Transducer Sets
- Safedraw® and Safedraw®-P closed loop arterial blood sampling system options available (please refer to next section)
- R.O.S.E.® (Resonance OverShoot Eliminator) included in selected sets
- Adult pressure sets include rotating male luer making it easier for the catheter connection
- \* DEHP-free Pressure Set

#### How Gabarith™ Sets can benefit you

- Unique Gabarith testing of sets, the total chain concept is designed for improved accuracy
- Merit sampling and zeroing stopcocks with on/off snap for confirmation of stopcock positioning
- Designed straight fluid pathway minimises the risk of air entrapment and the potential of subsequent damping of signal
- Robust monitor connection minimises water penetration and offers back pressure calibration facility
- Unique Easy Vent® zeroing deadender minimises the risk of touch contamination
- Transducer set configurations covering neonatal to adult pressure monitoring requirements

DISPOSABLE PRESSURE TRANSDUCER SETS



Safedraw® Safe Arterial Sampling  
Meritans DTXPlus®

Disposable Transducer Sets

Sterile, single use

Safedraw® safe arterial blood sampling systems are designed to reduce unnecessary blood discard and blood spillage. The system's easy-to-use, self-sealing, silicone septum facilitates easier blood sampling for clinical samples.


PRODUCT CODE	NO. OF DISPOSABLE TRANSDUCERS	FLUSH/ FLOW RATE (mL/HR)	TOTAL PRESSURE LINE LENGTH (CM)	PRESSURE LINE(S) COLOUR	ADDITIONAL INFORMATION	SAFEDRAW® SYRINGE VOLUME (mL)	BOX/CASE CONTENTS (UNITS)	
682051	1		155cm (125/30)	Clear	Single Admin Set. Single Transducer Set with Integrated Safedraw® Septum (Red)	2cc	5/25	
682053	2		A: 155cm (125/30) V: 150cm (125/25)	A: Clear V: Blue	Bifurcated Admin Set. Double Transducer Set with Integrated Safedraw® Septum (Red) and R.O.S.E on Venous Line	2cc	5/20	
682054	2		A: 155cm (125/30) V: 150cm (125/25)	A: Clear V: Blue	Bifurcated Admin Set. Double Transducer Set with Integrated Safedraw® Septum (Red)	2cc	5/25	
682055	2		A: 155cm (125/30) V: 150cm (125/25)	A: Clear V: Blue	Trifurcated Admin Set. Double Transducer Set with Integrated Safedraw® Septum (Red)	2cc	5/20	
682056	3		A: 155cm (125/30) V: 150cm (125/25) PA: 150cm (125/25)	A: Clear V: Blue PA: Yellow	Trifurcated Admin Set. Triple Transducer Set with Integrated Safedraw® Septum (Red)	2cc	5/20	
686337	1	3	190cm (160/30)	Red/Clear	Single Admin Set. Single 'Flexi' Transducer Set with Integrated Safedraw® Septum (Red)	2.5cc	5/25	
686585/ 689369*	1	3	190cm (75/85/30)	Red/Clear	Single Admin Set. Single 'Flexi' Transducer Set with 2 x Integrated Safedraw® Septums (Red)	2.5cc	5/20	
687714	1	3	190cm (160/30)	Red/Clear	No Admin Set. Single 'Flexi' Add-On Transducer Set with Integrated Safedraw® Septum (Red)	2.5cc	5/25	
682052	1		30 (Paed Set)	135cm (125/10)	Clear	Single Paediatric Transducer Set with Integrated Safedraw® Septum (Red) <b>For use with Infusion Pump</b>	1cc	5/25
688077	1	3 (Paed Set)	135cm (125/10)	Clear	Single Paediatric Transducer Set with Integrated Safedraw® Septum (Red) <b>For use with Infusion Pump</b>	2.5cc	5/25	
688849	1	3 (Paed Set)	29cm (22/7)	Clear	Single Paediatric Set with Integrated Safedraw® Septum (Red) <b>For use with Infusion Pump</b>	0.5cc	5/25	

\* Safedraw®-P DEHP-free sets

SAFEDRAW® ACCESSORIES

PRODUCT CODE	PRODUCT DESCRIPTION	MODEL	BOX/CASE CONTENT (UNITS)
682158	Safedraw® Vacutainer Transfer Device	TA-STV	25/125
682159	Safedraw® Luer-Lok Safe-Needle	TA-BPN	50/500





## Meritans DTXPlus® Flexi

### Disposable Transducer Sets

Sterile, single use

Meritans DTXPlus® Flexi Transducer Systems are a range of invasive pressure monitoring sets designed to maximise efficiency and ease of use by allowing the clinician to make line additions when needed. In the Flexi configuration, the Arterial disposable transducer system is used as the 'master' line. When clinically required, additional lines are connected to this 'master' line in a cascade. This is achieved via an additional stopcock positioned on the administration set side of the transducer, thereby having no impact on dynamic response.

This innovative Flexi solution, offers the adaptability of single or multiple transducer set configurations, and facilitates set selection, set-up time and minimises total cost.

Flexi sets are available in standard set format, Safedraw® or Safedraw®-P Closed-loop, arterial blood sampling formats.

PRODUCT CODE	NO. OF DISPOSABLE TRANSDUCERS	FLUSH/ FLOW RATE (mL/HR)	TOTAL PRESSURE LINE LENGTH (CM)	PRESSURE LINE(S) COLOUR	ADDITIONAL INFORMATION	SAFEDRAW® SYRINGE VOLUME (mL)	SEPTUM MODEL	BOX/CASE CONTENTS (UNITS)
686562	1	3	175cm (150/25)	Red	Single Admin Set. Single 'Flexi' Arterial Transducer Set	N/A	N/A	5/25
686495/ 689367*	1	3	175cm (160/15)	Blue	No Admin Set. Single 'Flexi' Add-On Venous Transducer Set	N/A	N/A	5/25
686496	1	3	175cm (175)	Yellow	No Admin Set. Single 'Flexi' Add-On PA Transducer Set	N/A	N/A	5/25
686572	1	3	175cm (160/15)	Clear	No Admin Set. Single 'Flexi' Add-On Universal Transducer Set	N/A	N/A	5/25
687715	1	3	175cm (160/15)	Blue	Single Admin Set. Single 'Flexi' Venous Transducer Set	N/A	N/A	5/25
686337	1	3	190cm (160/30)	Red/Clear	Single Admin Set. Single 'Flexi' Transducer Set with Integrated Safedraw® Septum (Red)	2.5cc	Safedraw® 1 x In-Line	5/25
686585/ 689369*	1	3	190cm (75/85/30)	Red/Clear	Single Admin Set. Single 'Flexi' Transducer Set with 2 x Integrated Safedraw® Septums (Red)	2.5cc	Safedraw® 2 x In-Line	5/20
687714	1	3	190cm (160/30)	Red/Clear	No Admin Set. Single 'Flexi' Add-On Transducer Set with Integrated Safedraw® Septum (Red)	2.5cc	Safedraw® 1 x In-Line	5/25
688396	1		193cm (163/30)	Red/Clear	Single Admin Set. Single 'Flexi' Transducer Set with non-shunting stopcock and Integrated Safedraw®-P Septum (Red)	2.5cc	Safedraw®-P 1 x In-Line	5/20
688389	1		193cm (163/30)	Red/Clear	Single Admin Set. Single 'Flexi' Transducer Set with Integrated Safedraw®-P Stopcock Septum (Red)	3cc	Safedraw®-P 1 x Stopcock	5/20
688420/ 689311*	1		213cm (183/30)	Red/Clear	Single Admin Set. Single 'Flexi' Transducer Set with Integrated Safedraw®-P Septum (Red)	2.5cc	Safedraw®-P 1 x In-Line	5/25
688406	1		203cm (53/120/30)	Red/Clear	Single Admin Set. Single 'Flexi' Transducer Set with 2 x Integrated Safedraw®-P Septums (Red)	2.5cc	Safedraw®-P 2 x In-Line	5/20
688513	1	3	193cm (73/90/30)	Red/Clear	Single Admin Set. Single 'Flexi' Transducer Set with 2 x Integrated Safedraw®-P Stopcock Septums (Red)	3cc	Safedraw®-P 2 x Stopcock	5/25

\* DTX Flexi DEHP-free set

### FLEXI SAFEDRAW® AND SAFEDRAW®-P ACCESSORIES

PRODUCT CODE	PRODUCT DESCRIPTION	BOX/CASE CONTENT (UNITS)
688446	Safedraw®-P Luer-lock Adapter System	50/500
682158	Safedraw® Vacutainer Transfer Device (TA-STV)	25/125
682159	Safedraw® Luer-Lok Safe-Needle (TA-BPN)	50/500

## DISPOSABLE PRESSURE TRANSDUCER SETS



## Safedraw®-P Safe Arterial Sampling Meritans DTXPlus®

Disposable Transducer Sets

Sterile, single use

Safedraw®-P arterial blood sampling systems are a unique range of invasive pressure monitoring sets designed for closed loop, safe arterial blood sampling procedures.

Whole blood samples are collected via direct luer-slip access into the unique Safedraw®-P split-septum blood sampling port, thereby negating the need for blood sampling adapters. The Safedraw®-P blood sampling port is either integrated directly into the tubing, or as an integral part of the special three-way sampling stopcock.

Safedraw®-P safe arterial blood sampling systems are designed to reduce unnecessary blood discard and blood spillage. The system's direct luer-slip access design facilitates easier blood sampling.

Available in Single, Double, Triple and 'Flexi' Disposable Transducer configurations, with Adult and Paediatric options.

**Product Codes are shown on next page**

### SAFEDRAW®-P ACCESSORIES

PRODUCT CODE	PRODUCT DESCRIPTION	BOX/CASE CONTENT (UNITS)
688446	Safedraw®-P Luer-slip Adapter System	50/500



## Safedraw®-P Safe Arterial Sampling Meritrans DTXPlus® Disposable Transducer Sets

Sterile, single use

PRODUCT CODE	NO. OF DISPOSABLE TRANSDUCERS	FLUSH/ FLOW RATE (mL/HR)	TOTAL PRESSURE LINE LENGTH (CM)	PRESSURE LINE(S) COLOUR	ADDITIONAL INFORMATION	SAFEDRAW SYRINGE VOLUME (mL)	SEPTUM MODEL	BOX/CASE CONTENTS (UNITS)
688413	1		183cm (153/30)	Clear	Single Admin Set. Single Transducer Set with Integrated Safedraw®-P Septum (Red)	2cc	1 x Inline	5/20
688410	1		213cm (183/30)	Red/Clear	Single Admin Set. Single Transducer Set with Integrated Safedraw®-P Septum (Red)	2cc	1 x Inline	5/25
688457	1		223cm (133/90)	Red/Clear	Single Admin Set. Single Transducer Set with Integrated Safedraw®-P Septum (Red)	3cc	1 x Inline	5/20
688538	1		78cm (63/15)	Red/Clear	Single Admin Set. Single Transducer Set with Integrated Safedraw®-P Septum (Red) with Velcro Strap for Arm mounting	1.5cc	1 x Inline	5/20
688403/ 689314*	1		A: 213cm (183/30) V: 150	A: Red/Clear V: Blue	Single Admin Set. Single Transducer Set with non-shunting stopcock and Integrated Safedraw®-P Septum (Red) (with Separately Packaged 150cm Blue pressure line)	2.5cc	1 x Inline	5/25
688407	1		A: 183cm (153/30) V: 125	A: Red/Clear V: Blue	Single Admin Set. Single Transducer Set with R.O.S.E., non-shunting stopcock, and Integrated Safedraw®-P Septum (Red) (with Separately Packaged 125cm Blue pressure line)	2cc	1 x Inline	5/20
688496	1		A: 193cm (163/30) V: 150	A: Red/Clear V: Blue	Single Admin Set. Single Transducer Set with non-shunting stopcock (with 150cm Blue pressure line) and Integrated Safedraw®-P Septum (Red)	2.5cc	1 x Inline (Red)	5/20
688416	2		A: 183cm (153/30) V: 150cm (125/25)	A: Clear V: Blue	Bifurcated Admin Set. Double Transducer Set with Integrated Safedraw®-P Septum (Red)	2cc	1 x Inline	5/20
688401	2		A: 213cm (183/30) V: 186cm (183/3)	A: Red/Clear V: Blue/Clear	Bifurcated Admin Set. Double Transducer Set with Integrated Safedraw®-P Septum (Red) and Integrated Safedraw®-P Stopcock Septum (Blue)	2.5cc	1 x Inline (Red) 1 x Stopcock (Blue)	5/25
688604	2		A: 223cm (133/90) V: 216cm (183/33)	A: Red/Clear V: Blue/Clear	Bifurcated Admin Set. Double Transducer Set with Integrated Safedraw®-P Septum (Red) and Integrated Safedraw®-P Stopcock Septum (Blue)	3cc	1 x Inline (Red) 1 x Stopcock (Blue)	5/25
688452	2		A: 223cm (133/90) V: 90cm	A: Red/Clear V: Blue	Trifurcated Admin Set. Double Transducer Set with Integrated Safedraw®-P Septum (Red) and Integrated Stopcock (Blue)	3cc	1 x Inline (Red)	5/20
688528	2		A: 223cm (133/90) V: 211cm (183/28)	A: Red/Clear V: Blue/Clear	Trifurcated Admin Set. Double Transducer Set with Integrated Safedraw®-P Septum (Red) and Integrated Safedraw®-P Stopcock Septum (Blue)	3.5cc	1 x Inline (Red) 1 x Stopcock (Blue)	5/20
688572	2		A: 233cm (143/90) V: 205cm (180/25)	A: Red/Clear V: Blue/Clear	Trifurcated Admin Set. Double Transducer Set with Integrated Safedraw®-P Septum (Red) and 2 x Sampling Stopcocks (Blue)	3cc	1 x Inline (Red)	5/20
688481	3		A: 213cm (183/30) V: 126cm (123/3) PA: 120cm	A: Red/Clear V: Blue/Clear PA: Yellow	Trifurcated Admin Set. Triple Transducer Set with Integrated Safedraw®-P Septum (Red) and Integrated Safedraw®-P Stopcock Septum (Blue)	3cc	1 x Inline (Red) 1 x Stopcock (Blue)	5/20
688396	1		193cm (163/30)	Red/Clear	Single Admin Set. Single 'Flexi' Transducer Set with non-shunting stopcock and Integrated Safedraw®-P Septum (Red)	2.5cc	1 x Inline	5/20
688389	1		193cm (163/30)	Red/Clear	Single Admin Set. Single 'Flexi' Transducer Set with Integrated Safedraw®-P Stopcock Septum (Red)	3cc	1 x Stopcock	5/20
688420/ 689311*	1		213cm (183/30)	Red/Clear	Single Admin Set. Single 'Flexi' Transducer Set with Integrated Safedraw®-P Septum (Red)	2.5cc	1 x Inline	5/25
688406	1		203cm (53/120/30)	Red/Clear	Single Admin Set. Single 'Flexi' Transducer Set with 2 x Integrated Safedraw®-P Septums (Red)	2.5cc	2 x Inline	5/20
688513	1		193cm (73/90/30)	Red/Clear	Single Admin Set. Single 'Flexi' Transducer Set with 2 x Integrated Safedraw®-P Stopcock Septums (Red)	3cc	2 x Stopcock	5/25
688408/ 689423*	1		143cm (128/15)	Clear	Single Paediatric Transducer Set with Integrated Safedraw®-P Septum (Red) <b>For use with Infusion Pump</b>	2.5cc	1 x Inline	5/25
688550/ 689312*	1		158cm (153/5)	Clear	Single Paediatric Transducer Set with Integrated Safedraw®-P Septum (Red) <b>For use with Infusion Pump</b>	2cc	1 x Inline	5/25
688600	1		158cm (153/5)	Clear	Single Paediatric Transducer Set with Integrated Safedraw®-P Septum (Red) <b>For use with Infusion Pump</b>	2cc	1 x Inline	5/25

\* Safedraw®-P DEHP-free sets

REUSABLE TRANSDUCER SETS



## RT2000™ Reusable Transducer Systems Overview

RT2000™ Reusable Transducer and Disposable Domes and Dome Systems are engineered to improve cost efficiencies by reducing disposal costs and waste while maintaining pressure waveform accuracy and reliability.

With the RT2000™ system, only the sterile disposable dome components are discarded after each patient use, while the reusable transducer is retained ready for the next patient’s use. This enables the cost-per-patient set-up to be lower than that of a disposable transducer system, and minimises the environmental impact by reducing the volume of hazardous waste.

The one-handed ‘snap-on’ dry coupling system for engagement of the disposable sterile dome set is simple and quick to perform, thereby minimising initial set-up time. The straight and visible fluid pathway is designed for the disposable dome to fill, bubble-free with minimal effort. The integrated flush device and venting EasyVent zeroing stopcock provide quick, effective fluid control.

RT2000™ systems are available in Flexi format and with the Safedraw® and Safedraw®-P safe arterial blood sampling systems. These systems are designed to reduce unnecessary blood discard and blood spillage. Whole blood samples are collected in the ‘closed-loop’ Safedraw® and Safedraw®-P sets either via the Safedraw® self-sealing silicone septum or via direct luer-slip access into the unique Safedraw®-P split-septum blood sampling port. The ‘closed-loop’ system design facilitates easier blood sampling.

Available in Standard and ‘Flexi’ Disposable Dome configurations, with Adult Safedraw® and Safedraw®-P options.

**See product table on next page.**

### RT2000™ REUSABLE TRANSDUCER SPECIFICATIONS

DETAIL	RT2000™ TRANSDUCER SPECIFICATIONS
Operating pressure:	-30 to 300 mmHg
Overpressure withstand	-400 to 5.170mmHg
Nominal sensitivity (gain):	5 μV/V/mmHg
Linearity:	± 1% over 0-50 mmHg
Hysteresis:	±3% over 50-300 mmHg
Unbalance (zero balance):	+40 mmHg
Thermal coefficient of sensitivity (typical):	+0.15%/ °C
Thermal coefficient of zero (typical):	± 1.0% mmHg / °C
Leakage:	<5μA at 220V 50Hz
Defibrillator withstand:	Withstand five repeated discharges in 5 minutes of 400J sinusoidal energy delivered into 50 ohms when all transducer wires are shorted and saline is in fluid column
Maximum Excitation Voltage:	10V DC or AC



## RT2000™ Reusable Transducer and Disposable Domes and Dome Systems

Sterile, single use

PRODUCT CODE	NO. OF DISPOSABLE DOMES	FLUSH/ FLOW RATE (mL/HR)	TOTAL PRESSURE LINE LENGTH (CM)	PRESSURE LINE(S) COLOUR	ADDITIONAL INFORMATION	SAFEDRAW SYRINGE VOLUME (mL)	MODEL	BOX/CASE CONTENTS (UNITS)
682266	N/A	N/A	N/A	N/A	RT2000™ reusable pressure transducer	N/A	RT2000™	10/100
682263	1	3	N/A	N/A	RT2000™ Disposable dome with integral flow/flush device (3cc/hr nominal flow rate) and Zeroing Stopcock	N/A	DD2003	10/100
682265	1	N/A	N/A	N/A	RT2000™ Disposable dome without integral flow/flush device or Zeroing Stopcock	N/A	DD2000	10/100
686958	1	3	175cm (150/25)	Clear	Single Admin Set. Single RT2000™ 'Universal' Dome Set	N/A	N/A	5/25
687886	1	3	175cm (150/25)	Red	Single Admin Set. Single 'Flexi' RT2000™ Arterial Dome Set with non-shunting stopcock	N/A	N/A	5/20
686961	1	3	175cm (150/25)	A: Red V: Blue	Single Admin Set. Single 'Flexi' RT2000™ Arterial Dome Set with non-shunting stopcock (with Separately Packaged 150cm Blue pressure line)	N/A	N/A	5/20
687887	1	3	175cm (160/15)	Blue	No Admin Set. Single 'Flexi' RT2000™ Add-On Venous Dome Set	N/A	N/A	5/25
687888	1	3	175cm (175)	Yellow	No Admin Set. Single 'Flexi' RT2000™ Add-On PA Dome Set	N/A	N/A	5/25
688397	1	3	193cm (163/30)	Red/Clear	Single Admin Set. Single 'Flexi' RT2000™ Arterial Dome Set with non-shunting stopcock and Integrated Safedraw®-P Septum (Red)	2.5cc	Safedraw®-P 1 x Inline	5/20
688546	1	3	A: 180cm (150/30) V: 150	A: Red/Clear V: Blue	Single Admin Set. Single RT2000™ Arterial Dome Set with non-shunting stopcock (with 150cm Blue pressure line) and Integrated Safedraw®-P Septum (Red)	2.5cc	Safedraw®-P 1 x Inline	5/20

### RT2000™ SAFEDRAW® AND SAFEDRAW®-P ACCESSORIES

PRODUCT CODE	PRODUCT DESCRIPTION	BOX/CASE CONTENT (UNITS)
688446	Safedraw®-P Luer-lock Adapter System	50/500
682158	Safedraw® Vacutainer Transfer Device (TA-STV)	25/125
682159	Safedraw® Luer-Lok Safe-Needle (TA-BPN)	50/500

## DISPOSABLE PRESSURE TRANSDUCER SETS



## Meritans DTXPlus® Paediatric and Neonatal

### Disposable Transducer Sets



Sterile, single use

Meritans DTXPlus® Paediatric and Neonatal Transducer Systems are a range of 'Standard' and 'Closed loop' Safedraw® and Safedraw®-P safe arterial blood sampling invasive pressure monitoring sets specifically designed with the Paediatric and Neonatal patient in mind. Whole blood samples are collected in the 'closed-loop' Safedraw® and Safedraw®-P sets either via the Safedraw® self-sealing silicone septum or via direct luer-slip access into the unique Safedraw®-P split-septum blood sampling port. Safedraw® and Safedraw®-P safe arterial blood sampling systems are designed to reduce unnecessary blood discard and facilitate the safe return of diluted blood back to the infant. These systems are configured for use with an Infusion Pump.

PRODUCT CODE	NO. OF DISPOSABLE TRANSDUCERS	FLUSH/ FLOW RATE (mL/HR)	TOTAL PRESSURE LINE LENGTH (CM)	PRESSURE LINE(S) COLOUR	ADDITIONAL INFORMATION	SAFEDRAW® SYRINGE VOLUME (mL)	SEPTUM MODEL	BOX/CASE CONTENTS (UNITS)
682022	1	30 (Paed Set)	160cm (150/10)	Clear	Single Paediatric Transducer Set <b>For use with Infusion Pump</b>	N/A	N/A	5/50
682025	1	30 (Paed Set)	35cm (25/10)	Clear	Single Paediatric Transducer Set <b>For use with Infusion Pump</b>	N/A	N/A	5/50
682052	1	30 (Paed Set)	135cm (125/10)	Clear	Single Paediatric Transducer Set with Integrated Safedraw® Septum (Red) <b>For use with Infusion Pump</b>	1cc	Safedraw® 1 x Inline	5/25
688077	1	3 (Paed Set)	135cm (125/10)	Clear	Single Paediatric Transducer Set with Integrated Safedraw® Septum (Red) <b>For use with Infusion Pump</b>	2.5cc	Safedraw® 1 x Inline	5/25
688849	1	3 (Paed Set)	29cm (22/7)	Clear	Single Paediatric Set with Integrated Safedraw® Septum (Red) <b>For use with Infusion Pump</b>	0.5cc	Safedraw® 1 x Inline	5/25
688408	1	3 (Paed Set)	143cm (128/15)	Clear	Single Paediatric Transducer Set with Integrated Safedraw®-P Septum (Red) <b>For use with Infusion Pump</b>	2.5cc	Safedraw®-P 1 x Inline	5/25
688550/689312*	1	3 (Paed Set)	158cm (153/5)	Clear	Single Paediatric Transducer Set with Integrated Safedraw®-P Septum (Red) <b>For use with Infusion Pump</b>	2cc	Safedraw®-P 1 x Inline	5/25
688600	1	30 (Paed Set)	158cm (153/5)	Clear	Single Paediatric Transducer Set with Integrated Safedraw®-P Septum (Red) <b>For use with Infusion Pump</b>	2cc	Safedraw®-P 1 x Inline	5/25

### SAFEDRAW® AND SAFEDRAW®-P ACCESSORIES

PRODUCT CODE	PRODUCT DESCRIPTION	BOX/CASE CONTENT (UNITS)
688446	Safedraw®-P Luer-lock Adapter System	50/500
682158	Safedraw® Vacutainer Transfer Device (TA-STV)	25/125
682159	Safedraw® Luer-Lok Safe-Needle (TA-BPN)	50/500



# Meritans DTXPlus® Cardiac Catheterisation Laboratory Transducer Set Configurations

Sterile, single use

With ease of priming and de-bubbling, with 'standard' and 106cm (42") cable option and a choice of flush and no-flush versions, the Meritans DTXPlus® transducer sets facilitate accurate pressure monitoring within the Cath Lab environment.

High pressure manifolds are supplied with a choice of 'Open' or 'Finger grip' (Block) design, 2, 3 or 4 ports, and having either a 'Left' or 'Right' rotating adapter to best meet specific procedural needs and physician preferences.

Available in clear, semi-rigid PVC and Flexible Braided Polyurethane with or without Rotating Adapter, High Pressure Injection Lines are designed to provide optimal flexibility, easy debubbling and secure connections and are rated up to 1200 psi.

Control Syringes are available with choice of 'Thumb Ring' or 'Palm Pad' design, 'Finger Ring' vs 'Finger Grip', and 10mL vs 12mL volume options to best suit most clinical needs and hand sizes.

Choice of syringe Palm Pad designs, including 'soft touch' options, allows for full force palm injections and thumb ring design facilitates aspirations and injections whilst choice of finger grip or ring style to meet most clinical needs. Contrast Saver Administration set options available to reduce contrast wastage and reduce risk of air entering the system.

PRODUCT CODE	DISPOSABLE TRANSDUCER MODEL	FLUSH/ FLOW RATE (mL/HR)	MANIFOLD	PRESSURE LINE(S) COLOUR	ADDITIONAL INFORMATION	CHECK VALVES	CONTRAST SYRINGE	BOX/CASE CONTENTS (UNITS)
688759	Meritans DTXPlus® TNF-R with 106cm transducer cable	No Flush	3-way, 360° 'off-directed', right orientated, block body, Rotating Adapter	25cm High Pressure braided	Cath Lab set with Meritans DTXPlus® TNF-R transducer, 3-way 'off-directed' right orientated manifold with 25cm High pressure braided tubing, administration set (15 micron filter), Dual Check Valve and 12ml Contrast Syringe. (For use with Contrast Saver set P/N: 688758)	1 x Dual Check Valve	12ml with thumb/ finger rings and rotating adapter	5/20
688758	N/A	N/A	N/A	N/A	Contrast Saver Kit with 30ml burette and white stopcock. (For use with Meritans DTXPlus® DT Manifold Cath Lab Kits)	N/A	N/A	5/25
686331	Meritans DTXPlus® DT-XX	3cc/hr	2-way 'on-directed', right orientated, block body, Rotating Adapter	75cm Clear	Cath Lab set with Meritans DTXPlus® DT-XX pressure transducer, 75cm Clear pressure line, 2-way 'on-directed' right orientated manifold with Contrast administration set (200 micron filter), 12ml Contrast Syringe and administration set (15 micron filter).	N/A	12ml with thumb/ finger rings and rotating adapter	5/20
686549	Meritans DTXPlus® DT-XX	3cc/hr	N/A	150cm Red	Cath Lab set with Meritans DTXPlus® DT-XX pressure transducer with administration set (15 micron filter), 150cm Red pressure line, stopcock and High Pressure Rotating Adapter. (FOR PAEDIATRIC CATHETERISATION PROCEDURES)	N/A	N/A	5/25
686132	Meritans DTXPlus® TNF-R	No Flush	2-way 360° 'off-directed', right orientated, block body, Rotating Adapter	150cm Clear	Cath Lab set with Dual Check Valve, 2-way 'off-directed' right orientated manifold with 2 x Contrast administration sets (200 micron filter), 150cm Clear pressure tubing and 10ml Contrast Syringe. Separately packaged Meritans DTXPlus® TNF-R transducer with tap and ROSE (Resonance OverShoot Device).	1 x Dual Check Valve	10ml with thumb ring/finger grips and 0.5ml 'stop' ring	5/20
686237	Meritans DTXPlus® DT-XX	3cc/hr	N/A	25cm, 30cm and 75cm Clear	Cath Lab set with Meritans DTXPlus® DT-XX pressure transducer with administration set (15 micron filter), 75cm Clear pressure line, High Pressure stopcock with 30cm Clear pressure line and 25cm Clear pressure line and High Pressure Rotating Adapter	N/A	N/A	5/25
688744	Meritans DTXPlus® DT-XX	3cc/hr	2-way 'on-directed', right orientated, block body, Rotating Adapter	150cm Clear	Cath Lab set with Meritans DTXPlus® DT-XX pressure transducer with administration set (15 micron filter), 150cm Clear pressure line and 2-way 'on-directed' right orientated manifold.	N/A	N/A	5/25
688709	Meritans DTXPlus® TNF-R	No Flush	3-way 'on-directed', right orientated, block body, Rotating Adapter	130cm Clear, 120cm High Pressure	Cath Lab set with Meritans DTXPlus® TNF-R transducer with tap, 130cm Clear pressure tubing, 3-way 'on-directed' right orientated manifold with administration set (15 micron filter) with tap, 120cm High Pressure tubing and 12ml Contrast Syringe. (For use with Contrast Saver set P/N: 688758)	N/A	12ml with thumb/ finger rings and rotating adapter	5/20
688143	Meritans DTXPlus® TNF-R	No Flush	N/A	130cm Clear	Cath Lab set with Meritans DTXPlus® TNF-R pressure transducer with 3 way stopcock, 130cm Clear pressure line, 2 x 3 way stopcocks and administration set (15 micron filter).	N/A	N/A	5/25
687699	Meritans DTXPlus® DT-XX	3cc/hr	N/A	100cm Clear	Cath Lab set with 10ml Contrast Syringe. Separately packaged Meritans DTXPlus® DT-XX pressure transducer with administration set (15 micron filter), 100cm Clear pressure line, 3 way stopcock with High Pressure Rotator Adapter and Dual Check Valve with contrast administration set (15 micron filter).	1 x Dual Check Valve	10ml with thumb ring/finger grips	5/20

ACCESSORIES



## Meritans DTXPlus® Accessories

Sterile, single use

PRODUCT CODE	MALE/FEMALE OR MALE/MALE LINE CONNECTION	TUBING LENGTH	LINE COLOUR CODING	BOX/CASE CONTENT (UNITS)	PRIMING VOLUME (CM <sup>3</sup> /CM)	TUBE DIAMETER INNER (MM)	OUTER (MM)
<b>PRESSURE LINES (PVC)</b>							
682226	M/F	25cm	Clear	25/125	0.019	1.575 (± 0.076)	3.175 (± 0.051)
682229	M/F	100cm	Clear	25/125	0.019	1.575 (± 0.076)	3.175 (± 0.051)
682230	M/F	150cm	Clear	25/125	0.019	1.575 (± 0.076)	3.175 (± 0.051)
682232	M/F	200cm	Clear	25/125	0.019	1.575 (± 0.076)	3.175 (± 0.051)
682234	M/M	100cm	Clear	25/125	0.019	1.575 (± 0.076)	3.175 (± 0.051)
682079	M/M	150cm	Clear	25/125	0.019	1.575 (± 0.076)	3.175 (± 0.051)
682236	M/M	200cm	Clear	25/125	0.019	1.575 (± 0.076)	3.175 (± 0.051)
682237	M/F	100cm	Red	25/125	0.019	1.575 (± 0.076)	3.175 (± 0.051)
685185	M/F	150cm	Red	25/125	0.019	1.575 (± 0.076)	3.175 (± 0.051)
685190	M/M	125cm	Blue	25/125	0.019	1.575 (± 0.076)	3.175 (± 0.051)
685171	M/M	150cm	Blue	25/250	0.019	1.575 (± 0.076)	3.175 (± 0.051)
686619	M/F	150cm	Yellow	25/125	0.019	1.575 (± 0.076)	3.175 (± 0.051)
<b>CONNECTA PE PRESSURE MONITORING LINES (POLYETHYLENE)</b>							
682194	M/F	15cm	Clear (Red/Blue Stickers)	25/250	0.011	1.20 (± 0.05)	2.50 (+ 0.1/- 0.05)
682185	M/F	50cm	Clear (Red/Blue Stickers)	25/250	0.011	1.20 (± 0.05)	2.50 (+ 0.1/- 0.05)
682189	M/F	100cm	Clear (Red/Blue Stickers)	25/250	0.011	1.20 (± 0.05)	2.50 (+ 0.1/- 0.05)
682191	M/F	150cm	Clear (Red/Blue Stickers)	25/250	0.011	1.20 (± 0.05)	2.50 (+ 0.1/- 0.05)
682192	M/F	200cm	Clear (Red/Blue Stickers)	25/250	0.011	1.20 (± 0.05)	2.50 (+ 0.1/- 0.05)
682186	M/M	50cm	Clear (Red/Blue Stickers)	25/250	0.011	1.20 (± 0.05)	2.50 (+ 0.1/- 0.05)
682190	M/M	100cm	Clear (Red/Blue Stickers)	25/250	0.011	1.20 (± 0.05)	2.50 (+ 0.1/- 0.05)
682193	M/M	200cm	Clear (Red/Blue Stickers)	25/250	0.011	1.20 (± 0.05)	2.50 (+ 0.1/- 0.05)

**PRODUCT CONTENTS**

- Pressure Lines (PVC) for accurate transmission of the pressure waveform
- PVC Pressure Lines maximum pressure tolerance (fluid filled) 200 psi
- Connecta PE (PolyEthylene) lines as a non PVC alternative
- Colour coded stickers supplied with Connecta PE to facilitate line identification
- Connecta PE Pressure Lines maximum pressure tolerance (fluid filled) 43 psi

**PRESSURE ACCESSORIES (NON-STERILE)**

PRODUCT CODE	MODEL	DESCRIPTION	BOX/CASE CONTENT (UNITS)
682093	TBG	Transducer BackGround for pole mounting (up to 4 transducers)	10/100
682274	TMM	Snap-On TMM Cable Organiser	5/50
682275	TMM/TBG	TBG with Snap-On TMM Cable Organiser (P/Ns 682093, 682274)	5/50
10026379	'O' Ring	Replacement 'O' Ring for Interface Cable	50
684003	Cable Checker	Transducer tester for monitor calibration	1
682104	UMM	Universal Manifold Mount	5/50





## Meritrans DTXPlus® Transducer Interface Cables


Merit supplies compatible interface cables for the majority of currently available blood pressure monitors. Below are listed the most common interface cables for use with Merit's Meritrans DTXPlus® transducers. If your requirements are not listed, please contact your local Merit office for assistance.

When purchasing cables, please specify:

1. Monitor type.
2. Model number and/or serial number.

PRODUCT CODE	MONITOR	MODEL
684078 / 684266	Datascope / Mindray	All Models (Passport, Accutorr, CS100/CS300 IABP)
684081 / 684267	Mindray	BeneView T5, T8, T1; IPM (8, 10, 12); BeneHeart D6
684081 / 684267	Philips	78341A, 78342A, etc. / 78354C, 78205A-D, 12-pin / 78352, CMS (Merlin line), M1006A/B, M1026A / M1166A / Viridia / HP Scop Eurelie 24C, CMS2001, V24, V26, M2, M3, M4, Intelivue
684082 / 684268	Dräger / Siemens	400 Series, 960, 961, Sirecust 4041, 700, 1260, 342, 732, SC7000, SC8000, SC9000, Infinity Gamma, Infinity Delta, Infinity Kappa, Infinity Vista all with 10 pin connector
684085 / 684269	Tektronics	400 Series
684085 / 684269	Vitalert	2000
684085 / 684269	Vittek	Series 400, 500
684085 / 684269	Welch Allyn	1500
684085 / 684269	Zoll	Zoll Defibrillators
684085 / 684269	Spacelabs	Series 400, 500, 90308/11/31, MOD-CD
684085 / 684269	Protocol Propaq	Protocol Propaq series; Propaq 106EL
684085 / 684269	Nellcor	Escort 300 MDE
684085 / 684269	Criticon / GE Marquette	Dynamap
684085 / 684269	Mindray	PM8000, PM7000, PM9000, Spectrum; MEC Series; Passport/Spectrum
684090 / 684277	Nihon Kohden	MFV1100, 1200, 1230A; Lifescope 8, 9, 11, 12, 14, BSM-8300J/K, 3000, BSM4100
684096	Marquette (GE)	7000 Series, 7010 Mod 7200
684102 / 684270	GE / Marquette	GE Monitor with Marquette 'front end': 7009 / Tramscope Series(600SL) - Eagle-Solar 8000, Solar 9000, Eagle3000/1000, Mab Lab 6Fx3F, Tram 450SL, 6200, Tram600A, EK36, Dash 3000, Dash 4000, Dash 5000
684104 / 684271	GE / Datex Ohmeda	GE Monitor with Datex-Ohmeda 'front end': SYS 800, MOD 8641, SIVA 3000, Cardiocap, Cardiocap II, AS3, S5, Aestiva, ADU(S5) and Aisys
684109 / 684272	No Monitor Connection	(Flying Leads) Shield: Dependent on Monitor. Black: + Signal; White: - Signal. Green: + Excitation; Red: - Excitation
684128	Dräger	PM8040, PM8060, Dräger Cato, Cicero & Julian Anes. Machine, Vitara
684128	E for M	RM300 (E for M/Honeywell)
684129 / 684278	Kontron, Arrow	Minimon 7131 etc., 7132 etc., K2000, 7272 / Supermon 7210 / Cart 2 Plus Balloon Pump, Oxford ACAT1
684138	Jostra	
684178 / 684273	Dräger / Siemens	SC6000P, 6002, 5000, 6004, 7000, Infinity Gamma, Infinity Delta, Infinity Kappa, Infinity Vista all with 7 pin connector
684181	Prucka Eng. Inc.	Cardiolab II / 1.1, AMP II
684188	Stöckert (Sorin)	S3 / S5
684227	Spacelabs (Vittek) / BCI	400, 500, 90485, Alpha-PC, 90603A, BCI 9100, 90303B-Mod90470-05, CD-CV, 90369, 90385, 90470, 90496, Ultraview 1600
684230	Physiocontrol/Medtronic	Lifepack 12 (6 pin)
684236	Dräger / Siemens	Hemo Pod Adapter (One adapter for each Hemo2 Pod and two for each Hemo4/HemoMed Pod)
684243	Stöckert (Sorin)	Heart Lung Machine
JP-900P	Nihon Kohden	NK9000

ACCESSORIES



**684109 / 684272 – Pin N/A - LF**

**No Monitor Connection** ('flying lead')

Wiring Guide:

(+) Excitation Voltage

(-) Excitation Voltage

(+) Signal

(-) Signal

Dependant on monitor specification



**684085 / 684269 – 6 Pin M - LF**

**Criticon/GE Marquette** - Dynamap

**Mindray** - PM/MEC Series, Spectrum

**Nellcor** - Escort 300 MDE

**Physiocontrol (Medtronic)** - Lifepack12

**Protocol Propaq** - Protocol Propaq series; Propaq 106EL;

**Spacelabs** - Series 400, 500, 90308/11/31,

MOD-CD, CD-CV Tektronics - 400 Series

**Vitalert** - 2000

**Vitatek** - Series 400, 500

**Welch Allyn** - 1500

**Zoll** - Zoll Defibrillators



**684078 / 684266 - 6 Pin F - LF**

**Datascope/Mindray** - All Models

Passport, Accutorr, CS100/CS300 IABP



**684128 – 10 Pin M - LF**

**E for M/Honeywell**

RM300

**Dräger**

PM8040, PM8060, Dräger Vitara, Cato,

Cicero & Julian Anes. machine

**Philips**

CM140



**684082 / 684268 – 10 Pin M - LF**

**Dräger/Siemens**

400 Series, 960, 961, Sirecust 4041,700,1260, 342, 732, Sirecust 342, 4XX, 7XX, 9XX, 12XX


SC7000, SC8000, SC9000, Infinity Gamma, Delta, Kappa, Vista all with 10 pin connector



**684178 / 684273 – 7 Pin M - LF**

**Dräger/Siemens**

SC6000P, 6002, 5000, 6004, 7000, Infinity Gamma, Delta, Kappa, Vista all with 7 pin connector

**684236 – 15 Pin M****Dräger/Siemens**

Hemo Pod Adapter (One adapter for each Hemo2 Pod & two for each Hemo4/HemoMed Pod)

**684104 / 684271 – 10 Pin F - LF****GE / Datex-Ohmeda**

GE Monitor with Datex-Ohmeda 'front end':  
SYS 800, MOD 8641, SIVA 3000, Cardiocap,  
Cardiocap II, AS3, S5, Aestiva, ADU(S5) and Aisys

**684102 / 684270 - 11 Pin M - LF****GE / Marquette**

GE Monitor with Marquette 'front end': 7009 / Tramscope Series(600SL), Eagle-Solar 8000, 9000, 3000/1000, Mab Lab 6Fx3F, Tram 450SL, 6200, Tram600A, EK36, Dash 3000, Dash 4000, Dash 5000, Careescape B650

**684138 – 5 Pin M****Jostra**

(Heart pump)

**684129 / 684278 – 12 Pin M - LF****Kontron, Arrow**

Minimon 7131, 7132 etc, K2000, 7272 / Supermon 7210 / Cart 2 Plus, Balloon Pump, Oxford ACAT1

**684096 – 8 Pin M - LF****GE / Marquette**

7000 series, 7010, Mod 7200

**684090 / 684277 – 5 Pin M - LF****Nihon Kohden**

MFV1100, 1200, 1230A; Lifescope 8, 9, 11, 12, 14, BSM-8300J/K, 3000, BSM4100

**JP - 900P (OTCN KB559)****Nihon Kohden**

NK9000

ACCESSORIES



**684081 / 684267 – 12 Pin M - LF**

**Philips**

*CMS (Merlin line), M1006A/B, M1026A / M1166A / Viridia / HP Scop Eurelie 24C, CMS2001, V24, V26, M2, M3, M4, Intelivue*

**Mindray**

*BeneView T5, T8, T1; IPM (8, 10, 12); BeneHeart D6*



**684181 – 8 Pin M**

**Prucka Eng Inc.**

*CARDIOLAB II / 1.1, AMP II*



**684227 – 6 Pin M - LF  
(identical to P/N 684085)**

**Spacelabs (Vitatek)**

*BCI - 400, 500, 90485, Alpha-PC, 90603A, BCI 9100, 90303B-Mod90470-05, CD-CV, 90369, 90385, 90470, 90496, Ultraview 1600*



**684188 - 4 Pin M**

**Stöckert (Sorin)**

*S3, S5*



**684243 – 4 Pin M - LF**

**Stöckert (Sorin)**

*Heart Lung Machine*

**684056 Model: SP4042****In-Line Injectate Sensor Cable**

For use with Abbott Oximetrix, Arrow, Baxter, Edwards 9520A or Com1, Mennen Horizon 2000 series (Cable 800-030-200), Spacelebs Alpha PC or PCMS (bifurcated cable model 306655-001), Datex AS3 and Critikit Injectate Set or SP5045 Sensor Housing

**684062 Model: SP4045****In-Line Injectate Sensor Cable**

For use with Nihon Kohden, Spacelabs (bifurcated cable model 306655-002), Fukuda Denshi IB5006 (DSCP), Dynascope DS3300 and Critikit Injectate Set or SP5045 Sensor Housing

**684057 Model: SP4043****In-Line Injectate Sensor Cable**

For use with Siemens bifurcated cable, Siemens Pod system (model 8419145) and Critikit Injectate Set or SP5045 Sensor Housing

**684063 Model: SP4046****In-Line Injectate Sensor Cable**

For use with Philips Cardiac Output module (trifurcated cable M1641A or M1642A) and Critikit Injectate Set or SP5045 Sensor Housing


**DISPOSABLE PRESSURE TRANSDUCER SETS**

# Meritans® Pressure Transducer Sets with Cardiac Catheterisation Kit Configuration

## MERITANS® PRESSURE TRANSDUCERS

With reliable precision, customised configurations, and a clear fluid path, Meritans® pressure transducers can be customised to fit the needs of your physicians and labs. Our transducers are available in both hand-held and Paley setups; Meritans® can also be configured with a broad range of Merit accessories.

**MERITANS® BLOOD PRESSURE TRANSDUCERS** with winged housing. Packaged 25 units per box, 4 boxes per case.



12" (30 cm)  
Pigtail Cable


48" (122 cm)  
Pigtail Cable

**CATALOGUE NUMBER** MER100

MER148



**MERITANS® BLOOD PRESSURE TRANSDUCER** with non-winged housing. Packaged 25 units per box, 4 boxes per case.



12" (30 cm) Pigtail Cable


48" (122 cm) Pigtail Cable

**CATALOGUE NUMBER**

MER212

MER200

**MERITANS® MER100** with bonded standard pressure 1-way and 3-way Marquis® stopcocks. Packaged 25 kits per box, 4 boxes per case.




12" (30 cm) Pigtail Cable

**CATALOGUE NUMBER**

K11-00096

**MERITANS® MER100** with two bonded standard pressure 1-way Marquis stopcocks. Packaged 15 kits per box, 4 boxes per case.




12" (30 cm) Pigtail Cable

**CATALOGUE NUMBER**

K11-00098

## Meritrans® Pressure Transducers


**MERITRANS® MER100** with two bonded standard pressure 3-way Marquis® stopcocks. Packaged 15 kits per box, 4 boxes per case.



12" (30 cm) Pigtail Cable

CATALOGUE NUMBER K11-00099


**MERITRANS® MER100** with bonded standard pressure 3-way and 1-way Marquis stopcocks. Packaged 15 kits per box, 4 boxes per case.



12" (30 cm) Pigtail Cable

CATALOGUE NUMBER K11-00100


**MERITRANS® MER200** with bonded standard pressure 1-way Marquis stopcock. Packaged 25 kits per box, 4 boxes per case.



48" (122 cm) Pigtail Cable

CATALOGUE NUMBER K11-00137

**MERITRANS® MER200** with a 3-way Marquis stopcock bonded to both ends. Packaged 15 kits per box, 4 boxes per case.




48" (122 cm) Pigtail Cable

CATALOGUE NUMBER K11-00144

DISPOSABLE PRESSURE TRANSDUCERS

# Meritrans® Pressure Transducers


**MERITRANS® MER100** with 3-way Marquis® stopcocks bonded to the female luer, flush device and 72" (183 cm) small bore fluid administration line with vented/filtered drip chamber. Packaged 15 kits per box, 4 boxes per case.



12" (30 cm) Pigtail Cable

CATALOGUE NUMBER K11-00280


**MERITRANS® MER200** with a 3-way Marquis Stopcock (medium pressure, rotating male adapter) and a 1-way Marquis stopcock. Packaged 15 kits per box, 4 boxes per case.



48" (122 cm) Pigtail Cable

CATALOGUE NUMBER K11-00392


**MERITRANS® MER200** with a 3-way Marquis stopcock and a 1-way Marquis stopcock. Packaged 15 kits per box, 4 boxes per case.



48" (122 cm) Pigtail Cable

CATALOGUE NUMBER K11-00202

**MERITRANS® MER100** with a 1-way Marquis stopcock. Packaged 15 kits per box, 4 boxes per case.




48" (122 cm) Pigtail Cable

CATALOGUE NUMBER K11-01321



# Meritans® Pressure Transducers

MERITANS® MER200 with a 3-way Marquis® Stopcock and a 1-way Marquis stopcock. Packaged 15 kits per box, 4 boxes per case.



48" (122 cm) Pigtail Cable


CATALOGUE NUMBER K11-00419



DISPOSABLE PRESSURE TRANSDUCER SETS

# TRAM® Manifolds with Integrated Transducer

The TRAM® Manifold with Integrated Transducer gives a single operator the convenience to de-bubble, zero and rezero within the sterile field. The TRAM® has optimised clinician comfort by reducing hand fatigue with low torque, easy-to-turn wide-spaced handles, even when utilising higher injection pressures. The Devos™ TRAM® manifold has a check relief valve. Packaged 10 units per box, 4 boxes per case.



1 Port

2 Port

3 Port

4 Port

**200 PSI (14 BAR)**


ON Handles Wide				
Standard	201LTPHN-R		203LTPHN-R	
OFF Handles Wide		202LTPHWF-R	203LTPHWF-R	204LTPHWF-R
	Devos			2V4LTPHWF-R
	Devos with 360° handle			2V4LTPHWF-RT
	Standard		203LTPHF-R	

**500 PSI (34 BAR)**

ON Handles Wide		502LTPHWN-R	503LTPHWN-R	504LTPHWN-R
Standard	501LTPHN-R		503LTPHN-R	
OFF Handles Wide		502LTPHWF-R	503LTPHWF-R	504LTPHWF-R
	Standard		503LTPHF-R	

## Tram-P™ Manifolds with Integrated Transducer

The TRAM-P™ features a pressure transducer on the first port. By positioning the compensation line on an IV pole rather than on a table, the TRAM-P™ maintains a constant calibrated reading. With this configuration, the manifold can be moved to any position without affecting the accuracy of the pressure reading. Comfortable wide-spaced finger grip option. The Devos™ TRAM-P™ manifold has a check relief valve. Packaged 10 units per box, 4 boxes per case.




### 200 PSI (14 BAR)

ON Handles Wide		202LTCHWN-R	203LTCHWN-R	204LTCHWN-R
Standard	201LTCHN-R		203LTCHN-R	
OFF Handles Wide		202LTCHWF-R	203LTCHWF-R	204LTCHWF-R
Devos				22V4LTCHWF-R
Devos with 360° handle				2V4LTCHWF-RT
Standard			203LTCHF-R	

### 500 PSI (34 BAR)

ON Handles Wide		502LTCHWN-R	503LTCHWN-R	504LTCHWN-R
Standard	501LTCHN-R		503LTCHN-R	
OFF Handles Wide		502LTCHWF-R	503LTCHWF-R	504LTCHWF-R
Standard	501LTCHF-R		503LTCHF-R	




TRAM-P™  
Transducer and Manifold

PRESSURE LINES & ACCESSORIES


# Meritrans® Accessories

## MERITRANS® ORGANIZER

Plate with coloured labels. Packaged 5 units per box.



2 Port



4 Port


CATALOGUE NUMBER

MER110


MER111

## CLIPPER POLE MOUNT

Packaged 1 unit per box.



Rotates to Vertical or Horizontal Position



Fixed Position

CATALOGUE NUMBER


MER112

MER113

NOTE: Autoclavable Clipper® (MER113A) is available upon special request. Please allow six weeks lead time for delivery.

## MERITRANS® REUSABLE CABLE SERIES


Includes 19 standard cables for use with most monitor configurations. Compatible cables for your monitoring system available through your Merit representative.



## Meritrans® Cables

### MERITRANS® REUSABLE CABLE SERIES

Includes 19 standard cables for use with most monitor configurations. Compatible cables for your monitoring system available through your Merit representative.



## PRESSURE LINES &amp; ACCESSORIES

# The Mentor™ Pressure Transducer Simulator/Tester

**Accuracy:**

±1% of Reading or ±1 mmHg, whichever is greater

**Sensor Overpressure:**

1000mmHg

**Mentor Working Pressure Range:**

+300 mmHg, displays "OPr" when outside this range

**Transducer Simulator:**

- Sensitivity: 5  $\mu\text{V/V/mmHg}$ , +/- 1%
- Excitation: 2 to 10 vdc, or vac rms to 5 khz
- Input Impedance: 2300  $\Omega$ , 330  $\Omega$  with Reusable Cable
- Output Impedance: 300  $\Omega$
- Asymmetry: Less than 1%

Electronic, hand-held simulator to help you troubleshoot your haemodynamic system

**MER500**

With 0.2 Micron Microbiological Filter

**MER555**





# THERMODILUTION & PULMONARY ARTERY CATHETERS


---

**32** Exacta® Percutaneous Sheath Introducer Sets

**33** Criticath® PA Catheters

**34** Critikit™ Cardiac Output Injectate Kits and Accessories

## PERCUTANEOUS SHEATH INTRODUCER SETS



# Exacta® Percutaneous Sheath Introducer Sets

Sterile, single use

### PERCUTANEOUS SHEATH INTRODUCERS FOR THE INTRODUCTION OF PULMONARY ARTERY CATHETERS

PRODUCT CODE	SIZE*	INTRODUCER LENGTH	CANNULA INTRODUCER		NEEDLE INTRODUCER		GUIDEWIRE OD & LENGTH	CONTAMINATION SHIELD LENGTH	HIGH-FLOW STOPCOCK	SYRINGE	SCALPEL	OBTURATOR	BOX/CASE CONTENT (UNITS)
			GAUGE	LENGTH	GAUGE	LENGTH							
680128**	6Fr	100mm	16G	70mm	17G	70mm	0.9 x 450mm	800mm	Yes	Yes	Yes	Yes	5/50
680140**	7Fr	100mm	16G	70mm	17G	70mm	0.9 x 450mm	800mm	Yes	Yes	Yes	Yes	5/50
680126**	8Fr	100mm	16G	70mm	17G	70mm	0.9 x 450mm	800mm	Yes	Yes	Yes	Yes	5/50
680132	Link Kit 8.5Fr	100mm	16G	70mm	17G	70mm	0.9 x 450mm	-	Yes	Yes	Yes	Yes	5/50
680124**	8.5Fr	100mm	16G	70mm	17G	70mm	0.9 x 450mm	800mm	Yes	Yes	Yes	Yes	5/50

\* For introduction of 1Fr smaller size catheters

\*\* Exacta® – latex free


#### PRODUCT DESCRIPTION

- Radiopaque FEP sheath with integral sidearm
- Tuohy Borst adapter and haemostasis valve
- Matched dilator with luer lock collar
- Exacta® sets include a contamination shield with Tuohy Borst lock
- Two separately packaged sterile obturators supplied per box of 5 Exacta® sheath systems

#### How Merit percutaneous introducers can benefit you

- Radiopaque FEP sheath is designed to reduce the risk of thrombus
- Flexible hub minimises potential for kinking of introducers
- Suture point facilitates skin fixation
- Profiled sheath tip and matched dilator is designed to facilitate a smooth atraumatic insertion
- Dilator luer locking collar secures dilator to sheath hub during insertion
- Latex free, self-sealing haemostasis valve is designed to minimise the risk of blood loss
- Exacta® sets include the Venaguide™ Seldinger wire advancing device, which allows single-handed wire insertion





## Criticath® Pulmonary Artery Catheters

Sterile, single use

### PULMONARY ARTERY CATHETERS FOR MONITORING RIGHT HEART PRESSURES AND CARDIAC OUTPUT MEASUREMENT

PRODUCT CODE	MODEL	SIZE	LENGTH	NO. OF LUMENS	PA WEDGE PRESSURE	CV PRESSURE	CARDIAC OUTPUT	CVP MEDIC. PORT	ADDITIONAL INFORMATION	BOX/CASE CONTENT (UNITS)
680349	SP5105	5Fr	80cm	4	Yes	Yes	Yes		PVC Catheter Material	5/20
680103	SP5107U	7Fr	110cm	4	Yes	Yes	Yes			5/20
680381	SP5107U-14	7Fr	110cm	4	Yes	Yes	Yes		Pre-mounted Luer-Lok™ AC shield	5/20
680382	SP5107U+SP5045	7Fr	110cm	4	Yes	Yes	Yes		Pre-mounted In-Line Inject, sensor	5/20
680378	SP5507U	7.5Fr	110cm	5	Yes	Yes	Yes	Yes		5/20
680379	SP5507U-14	7.5Fr	110cm	5	Yes	Yes	Yes	Yes	Pre-mounted Luer-Lok™ AC shield	5/20
680380	SP5507U+SP5045	7.5Fr	110cm	5	Yes	Yes	Yes	Yes	Pre-mounted In-Line Inject, sensor	5/20

### MONITORING CATHETER

PRODUCT CODE	MODEL	SIZE	LENGTH	NO. OF LUMENS	PA WEDGE PRESSURE	CV PRESSURE	CARDIAC OUTPUT	CVP MEDIC. PORT	ADDITIONAL INFORMATION	BOX/CASE CONTENT (UNITS)
680372	SP5325	5Fr	80cm	2	Yes	No	No	No	Monitoring Catheter	5/20
680373	SP5325L	5Fr	110cm	2	Yes	No	No	No	Monitoring Catheter	5/20
680374	SP5327	7Fr	110cm	2	Yes	No	No	No	Monitoring Catheter	5/20
680375	SP5337	7Fr	110cm	3	Yes	Yes	No	No	Monitoring Catheter	5/20
680376	SP5337-13	7Fr	110cm	3	Yes	Yes	No	No	Monitoring Catheter with Pre-mounted push fit A/C shield	5/20
680377	SP5337-14	7Fr	110cm	3	Yes	Yes	No	No	Monitoring Catheter with Pre-mounted Luer Lok A/C shield	5/20

### PRODUCT DESCRIPTION

- A range of flow directed, balloon tipped catheters for both adult and (7F/7.5F) and paediatric use (5F)
- 7Fr and 7.5Fr catheter sizes with polyurethane body catheter material
- Use for monitoring right heart pressures and obtaining venous blood samples
- Thermodilution technique used to measure Cardiac Output
- Individually packed, sterile and non-pyrogenic

### How Criticath® can benefit you

- Rounded catheter tip design allows for smooth insertion and is designed to reduce the risk of vessel trauma
- Strong circular lumen feature is designed to reduce kinking and enhance pressure signal transfer
- Catheter depth markings facilitate accurate placement and catheter manipulation
- Smooth balloon/catheter junction is designed to reduce the risk of balloon rupture and thrombus formation
- Some models have a pre-mounted in-line injectate sensor housing or an anticontamination (AC) sheath on the catheter for user convenience



## Critikit™ Ice/Room Temperature Injectate Systems and Accessories

Sterile, single use

### CLOSED LOOP IN-LINE INJECTATE SYSTEMS FOR IMPROVED ACCURACY OF CARDIAC OUTPUT MEASUREMENT

PRODUCT CODE	MODEL	DESCRIPTION	STERILE	BOX/CASE CONTENT (UNITS)
680000*	SP4500**	Ice/Room temperature Injectate Critikit (with in-line sensor, control syringe & cooling coil)	Yes	4/8
680178*		Ice/Room temperature Injectate Critikit (with in-line sensor and control syringe, without cooling coil)	Yes	5/20
680006*	SP5045	In-line injectate sensor housing	Yes	25/250
682252*	SY10-RP	10 mL control syringe	Yes	10/50

#### PRODUCT DESCRIPTION

- Available in both iced and room temperature injectate options
- Closed circuit design offers an aseptic approach to cardiac output measurement
- Set comprises an injectate sensor housing with three-way stopcock, one-way check valve, cooling coil, delivery tubing and IV spike

\* Latex free

\*\* DEHP-free Critikit™

#### How Critikit™ can benefit you

- Closed loop system is designed to reduce the risk of line contamination
- Insulated volume-restricted syringe facilitates precise and reproducible injectate
- In-line injectate sensor optimises accuracy of temperature readings
- Syringe design allows smooth injectate action improving data validity

#### NON-STERILE ACCESSORIES

PRODUCT CODE	MODEL	DESCRIPTION	STERILE	BOX/CASE CONTENT (UNITS)
684056	SP4042	In-line injectate sensor cable. For use with Abbott Oximetrix, Arrow, Baxter Edwards 9520A or Com1, Mennen Horizon 2000 Series (cable 800-030-200), Spacelabs Alpha PC or PCMS (bifurcated cable model 306655-001), Datex AS3 and Critikit™ Injectate set or SP5045 sensor housing.	No	1
684057	SP4043	In-line injectate sensor cable. For use with Siemens bifurcated cable, Siemens Pod system (Model 8419145) and Critikit™ Injectate™ set or SP5045 sensor housing.	No	1
684062	SP4045	In-line injectate sensor cable. For use with HemoPro (SP1445), Nikon Kohden, Spacelabs (bifurcated cable model 306655-002), Fukuda Denshi IB5006 (DSCP), Dynascope DS3300 and Critikit™ Injectate set or SP5045 sensor housing.	No	1
684063	SP4046	In-line injectate sensor cable. For use with Hewlett Packard Cardiac Output Module (HP trifurcated cable M1641A or M1642A) and Critikit™ Injectate set or SP5045 sensor housing.	No	1



# CENTRAL VENOUS CATHETERS

---

- 36** Floswitch®
- 37** Secalon-T™ Emergency Catheters
- 38** Careflow™ Central Venous Catheters and Catheter Sets
- 43** Hydrocath Assure™ Antimicrobial Central Venous Catheters and Catheter Sets



**Floswitch®**  
Sterile, single use

- Floswitch® acts as an ON/OFF device on infusion lines, stopping the flow in order to carry out line procedures in closed circuit.
- Designed to reduce risk of patient emboli and contamination of hospital staff to blood.
- The Floswitch® is designed to increase safety.
- Floswitch® can be connected to any catheter fitted with a Luer connector.
- Winged to facilitate attachment.

**Floswitch® closed**

External compression of the extension tubing by means of a steel ball, provides risk-free unambiguous opening and closing of the infusion line.



**Floswitch® open**




**ORDERING INFORMATION**

PRODUCT CODE	DESCRIPTION	BOX/CASE CONTENT (UNITS)
681033	Floswitch®	25/250

Latex free  
CE marked

Class IIa according to the European Medical Devices Directive 93/42/EEC as amended by Directive 2007/47/EC Shelf life: 5 years  
Sterilisation method: Ethylene Oxide



## Secalon-T™

Sterile, single use

- Secalon-T™ emergency catheter designed for high flow rate infusions.
- PTFE (Polytetrafluoroethylene) material construction and incorporating a Floswitch®, Secalon-T™ is ideal for rapid short-term insertion into central veins or large peripheral vessels.
- Complete kit contains pre-puncture stylus and stabilising wing.

### ORDERING INFORMATION

PRODUCT CODE	GAUGE	LENGTH	CATHETER EXTERNAL DIAMETER	COLOUR	BOX/CASE CONTENT (UNITS)	FLOW RATE (ML/MIN)
681000	18G	90mm	1.2mm	Green	25/250	55
681002	16G	130mm	1.7mm	Grey	25/250	139
681004	16G	160mm	1.7mm	Grey	25/250	128
681006	14G	160mm	2.0mm	Orange	25/250	206



## Careflow™ Central Venous Catheters

### Single Lumen without Extension

Sterile, single use

- Careflow™ is a range of paediatric and adult polyurethane central venous catheters designed for catheterisation via the Seldinger technique.
- Range of sets contains a Nitinol guidewire and the Venaguide™ Seldinger wire-advancing device facilitating single-handed guidewire insertion.
- Tapered atraumatic catheter tip designed to reduce the risk of venous damage while improving catheter insertion.
- The streamlined junction hub enables ease of catheter suture and comfortable positioning on the patient's skin.
- Careflow™ catheters offer 1cm and 5cm interval depth markings to facilitate catheter placement and clear lumen size and hub colour coding indication.
- Floswitch® and Floswitch® Introducer Catheter options available throughout range.
- Pressure rated with a maximum pressure tolerance of 45 psi.

#### CAREFLOW™ SINGLE LUMEN SETS WITHOUT EXTENSION

CODE WITHOUT FLOSWITCH®	CODE WITH FLOSWITCH®	EXTERNAL DIAMETER	LENGTH	CANNULA INTRODUCER		FLOSWITCH™ INTRODUCER		NEEDLE INTRODUCER		GUIDEWIRE OD & LENGTH	SECONDARY FIXATION (SFD)	BOX/CASE CONTENT (UNITS)	FLOW RATE (ML/MIN)
				GAUGE	LENGTH	GAUGE	LENGTH	GAUGE	LENGTH				
681614	681623	2.5Fr (22G)	100mm	N/A	N/A	20G	45mm	22G	40mm	0.5 x 350mm		10/100	10
681639	681640 <sup>≠</sup>	3Fr (20G)	150mm	N/A	N/A	20G	45mm	22G	40mm	0.5 x 350mm		10/100	16
681643	681644	3Fr (20G)	200mm	N/A	N/A	20G	45mm	22G	40mm	0.5 x 450mm		10/100	15
681649	681650	4Fr (18G)	200mm	N/A	N/A	18G	60mm	17G	70mm	0.6 x 450mm		10/100	36
681669	681670	5Fr (16G)	150mm	N/A	N/A	15G	70mm	17G	70mm	0.9 x 450mm	Yes	10/100	71
681672	681673	5Fr (16G)	200mm	N/A	N/A	15G	70mm	17G	70mm	0.9 x 450mm	Yes	10/100	63
681706	681707	7Fr (14G)	150mm	N/A	N/A	14G	70mm	17G	70mm	0.9 x 450mm	Yes	10/100	190
681709	681710	7Fr (14G)	200mm	N/A	N/A	14G	70mm	17G	70mm	0.9 x 450mm	Yes	10/100	136

#### PRODUCT CONTENTS


##### Single lumen Careflow™ with Floswitch® set contents:

Careflow™ Polyurethane catheter, steel needle introducer, Floswitch® overneedle introducer, vessel dilator, Venaguide™ Seldinger wire-advancing device with Nitinol guidewire, Floswitch® revolving Luer-Lok™ system, luer lock plug and secondary fixation device (where shown).

##### Single lumen Careflow™ without Floswitch® set contents:

Careflow™ Polyurethane catheter, steel needle introducer, vessel dilator, Venaguide™ Seldinger wire-advancing device with Nitinol guidewire and secondary fixation device (where shown).

≠ Note: Part number 681640 contains 0.5 x 450mm guidewire.



## Careflow™ Central Venous Catheter Sets

Single Lumen  
Sterile, single use

CODE WITHOUT FLOSWITCH*	CODE WITH FLOSWITCH*	EXTERNAL DIAMETER	LENGTH	CANNULA INTRODUCER		FLOSWITCH* INTRODUCER		NEEDLE INTRODUCER		GUIDEWIRE OD & LENGTH	SECONDARY FIXATION (SFD)	BOX/CASE CONTENT (UNITS)	FLOW RATE (ML/MIN)
				GAUGE	LENGTH	GAUGE	LENGTH	GAUGE	LENGTH				
<b>1. CAREFLOW™ SINGLE LUMEN SETS WITH EXTENSION</b>													
681612		2.5Fr (22G)	60mm	22G	38mm	N/A	N/A	22G	40mm	0.5 x 350mm		5/50	12
681622		2.5Fr (22G)	100mm	22G	38mm	N/A	N/A	22G	40mm	0.5 x 350mm		5/50	10
681634		2.5Fr (22G)	200mm	22G	38mm	N/A	N/A	22G	40mm	0.5 x 650mm		5/50	7
681647		3Fr (20G)	150mm	20G	45mm	N/A	N/A	22G	40mm	0.5 x 450mm		5/50	18
681638		3Fr (20G)	200mm	20G	45mm	N/A	N/A	22G	40mm	0.5 x 650mm		5/50	15
681624		4Fr (18G)	60mm	20G	45mm	N/A	N/A	22G	40mm	0.5 x 350mm		5/50	47
681658		4Fr (18G)	200mm	20G	45mm	N/A	N/A	22G	40mm	0.5 x 650mm		5/50	35
681662		5Fr (16G)	150mm	16G	70mm	N/A	N/A	17G	70mm	0.9 x 450mm	Yes	5/50	65
681666		5Fr (16G)	200mm	16G	70mm	N/A	N/A	17G	70mm	0.9 x 700mm	Yes	5/50	60
681679		5Fr (16G)	300mm	16G	70mm	N/A	N/A	17G	70mm	0.9 x 700mm	Yes	5/50	46
681698		7Fr (14G)	150mm	16G	70mm	N/A	N/A	17G	70mm	0.9 x 450mm	Yes	5/50	159
681699		7Fr (14G)	200mm	16G	70mm	N/A	N/A	17G	70mm	0.9 x 700mm	Yes	5/50	149
681701		7Fr (14G)	300mm	16G	70mm	N/A	N/A	17G	70mm	0.9 x 700mm	Yes	5/50	136
<b>2. CAREFLOW™ SINGLE LUMEN SET WITH FLOSWITCH*</b>													
	681633	2.5Fr (22G)	200mm	N/A	N/A	20G	45mm	22G	40mm	0.5 x 450mm		10/100	7
<b>3. CAREFLOW™ SINGLE LUMEN WITH INTEGRAL FLOSWITCH*</b>													
	681660	5Fr (16G)	150mm	N/A	N/A	15G	70mm	17G	70mm	0.9 x 450mm	Yes	10/100	65
	681664	5Fr (16G)	200mm	N/A	N/A	15G	70mm	17G	70mm	0.9 x 700mm	Yes	10/100	60
	681678	5Fr (16G)	300mm	N/A	N/A	15G	70mm	17G	70mm	0.9 x 700mm	Yes	10/100	46
<b>4. CAREFLOW™ SINGLE LUMEN WITH EXTENSION AND FLOSWITCH*</b>													
	681615	2.5Fr (22G)	60mm	N/A	N/A	20G	45mm	22G	40mm	0.5 x 350mm		5/50	12
	681625	4Fr (18G)	60mm	N/A	N/A	20G	45mm	22G	40mm	0.5 x 350mm		5/50	47

### PRODUCT CONTENTS

#### 1. Single lumen Careflow™ set with extension contents:

Careflow™ Polyurethane catheter with clear extension, extension slide clamp, cannula introducer, steel needle introducer, vessel dilator, Venaguide™ Seldinger wire-advancing device with Nitinol guidewire, Luer-Lok™ plug, secondary fixation device (where shown) and 5ml syringe.

#### 2. Single lumen Careflow™ set with Floswitch® contents:


Careflow™ Polyurethane catheter with clear extension, steel needle introducer, Floswitch® overneedle introducer, PTFE coated soft-tip straight guidewire and Floswitch® revolving Luer-Lok system, luer lock plug.

#### 3. Single lumen Careflow™ set with integral Floswitch® contents:

Careflow™ Polyurethane catheter with clear extension and integral Floswitch®, Floswitch® overneedle introducer, steel needle introducer, vessel dilator, Venaguide™ Seldinger wire-advancing device with Nitinol guidewire, Luer-Lok plug and secondary fixation device.

#### 4. Single lumen Careflow set with extension and Floswitch® contents:

Careflow™ Polyurethane catheter with clear extension, steel needle introducer, Floswitch® overneedle introducer, vessel dilator, Venaguide™ Seldinger wire-advancing device with Nitinol guidewire, Floswitch® revolving Luer-Lok system and luer lock plug.



## Careflow™ Central Venous Catheter Sets

Double & Triple Lumen

Sterile, single use

### CAREFLOW™ DOUBLE LUMEN SETS

CODE WITHOUT FLOW SWITCH*	CODE WITH FLOW SWITCH*	EXTERNAL DIAMETER	LENGTH	CANNULA INTRODUCER		FLOSWITCH® INTRODUCER		NEEDLE INTRODUCER		GUIDEWIRE OD & LENGTH	SECONDARY FIXATION (SFD)	BOX/CASE CONTENT (UNITS)	FLOW RATE (ML/MIN)		
				GAUGE	LENGTH	GAUGE	LENGTH	GAUGE	LENGTH				D	M	P
681616	681617	4Fr (20G, 22G)	60mm	20G	45mm	20G	45mm	22G	40mm	0.5 x 350mm		5/50	20		8
681626	681627	4Fr (20G, 22G)	100mm	20G	45mm	20G	45mm	22G	40mm	0.5 x 350mm		5/50	19		6
681653	681652	4Fr (20G, 22G)	150mm	20G	45mm	20G	45mm	22G	40mm	0.5 x 450mm		5/50	16		4
681656#	681655	4Fr (20G, 22G)	200mm	20G	45mm	20G	45mm	22G	40mm	0.5 x 450mm		5/50	17		4
681619	681628	5Fr (18G, 20G)	60mm	18G	60mm	18G	60mm	17G	70mm	0.6 x 450mm	Yes	5/50	40		10
681683	681682	5Fr (18G, 20G)	150mm	18G	60mm	18G	60mm	17G	70mm	0.6 x 450mm	Yes	5/50	29		6
681687	681686	5Fr (18G, 20G)	200mm	18G	60mm	18G	60mm	17G	70mm	0.6 x 650mm	Yes	5/50	30		6
681703	681702	7Fr (16G, 18G)	150mm	16G	70mm	15G	70mm	17G	70mm	0.9 x 450mm	Yes	5/50	74		24
681714	681713	7Fr (16G, 18G)	200mm	16G	70mm	15G	70mm	17G	70mm	0.9 x 700mm	Yes	5/50	73		21
681718	681717	7Fr (16G, 18G)	300mm	16G	70mm	15G	70mm	17G	70mm	0.9 x 700mm	Yes	5/50	61		19

### PRODUCT CONTENTS

#### Double lumen Careflow™ with Floswitch® set contents:

Careflow™ double lumen Polyurethane catheter with clear extensions, steel needle introducer, Floswitch® overneedle introducer, vessel dilator, Venaguide™ Seldinger wire-advancing device with Nitinol guidewire, Floswitch® revolving Luer-Lok system, luer lock plugs and secondary fixation device (where shown).

#### Double lumen Careflow™ without Floswitch® set contents:

Careflow™ double lumen Polyurethane catheter with clear extensions, cannula introducer, steel needle introducer, vessel dilator, Venaguide™ Seldinger wire-advancing device with Nitinol guidewire, extension slide clamps, 5ml syringe, luer lock plugs and secondary fixation device (where shown). ≠ Note: Part number 681656 contains 0.5 x 650mm guidewire.

### CAREFLOW™ TRIPLE LUMEN SETS

CODE WITHOUT FLOW SWITCH*	CODE WITH FLOW SWITCH*	EXTERNAL DIAMETER	LENGTH	CANNULA INTRODUCER		FLOSWITCH® INTRODUCER		NEEDLE INTRODUCER		GUIDEWIRE OD & LENGTH	SECONDARY FIXATION (SFD)	BOX/CASE CONTENT (UNITS)	FLOW RATE (ML/MIN)		
				GAUGE	LENGTH	GAUGE	LENGTH	GAUGE	LENGTH				D	M	P
681620	681629	5Fr (18G, 20G, 20G)	60mm	18G	60mm	18G	60mm	17G	70mm	0.6 x 450mm	Yes	5/50	40	10	12
681630	681618	5Fr (18G, 20G, 20G)	100mm	18G	60mm	18G	60mm	17G	70mm	0.6 x 450mm	Yes	5/50	32	7	8
681691	681690	5Fr (18G, 20G, 20G)	150mm	18G	60mm	18G	60mm	17G	70mm	0.6 x 450mm	Yes	5/50	29	5	6
681695	681694	5Fr (18G, 20G, 20G)	200mm	18G	60mm	18G	60mm	17G	70mm	0.6 x 650mm	Yes	5/50	26	5	5
681721	681720	7Fr (16G, 18G, 18G)	150mm	16G	70mm	15G	70mm	17G	70mm	0.9 x 450mm	Yes	5/50	66	22	23
681725	681724	7Fr (16G, 18G, 18G)	200mm	16G	70mm	15G	70mm	17G	70mm	0.9 x 700mm	Yes	5/50	59	16	17
681729	681728	7Fr (16G, 18G, 18G)	300mm	16G	70mm	15G	70mm	17G	70mm	0.9 x 700mm	Yes	5/50	53	11	9

### PRODUCT CONTENTS


#### Triple lumen Careflow™ with Floswitch set contents:

Careflow™ triple lumen Polyurethane catheter with clear extensions, steel needle introducer, Floswitch® overneedle introducer, vessel dilator, Venaguide™ Seldinger wire-advancing device with Nitinol guidewire, Floswitch® revolving Luer-Lok system, luer lock plugs and secondary fixation device.

#### Triple lumen Careflow™ without Floswitch set contents:

Careflow™ triple lumen Polyurethane catheter with clear extensions, cannula introducer, steel needle introducer, vessel dilator, Venaguide™ Seldinger wire-advancing device with Nitinol guidewire, extension slide clamps, 5ml syringe, luer lock plugs and secondary fixation device.





## Careflow™ Central Venous Catheter Sets

### Quad & Quint Lumen

Sterile, single use

#### CAREFLOW™ QUAD LUMEN SETS

CODE WITHOUT FLOSWITCH®	CODE WITH FLOSWITCH®	EXTERNAL DIAMETER	LENGTH	CANNULA INTRODUCER		FLOSWITCH® INTRODUCER		NEEDLE INTRODUCER		GUIDEWIRE OD & LENGTH	BOX/CASE CONTENT (UNITS)	FLOW RATE (ML/MIN)				
				GAUGE	LENGTH	GAUGE	LENGTH	GAUGE	LENGTH			D	M1	M2	P1	P2
681731	681732	8.5Fr (14G, 16G, 16G, 18G)	150mm	16G	70mm	15G	70mm	17G	70mm	0.9 x 450mm	5/50	64	82	57	20	
681733	681734	8.5Fr (14G, 16G, 16G, 18G)	200mm	16G	70mm	15G	70mm	17G	70mm	0.9 x 700mm	5/50	58	76	45	19	

#### PRODUCT CONTENTS

##### Quad lumen Careflow™ with Floswitch® set contents:

Careflow™ quad lumen Polyurethane catheter with clear extensions, steel needle introducer, Floswitch® overneedle introducer, vessel dilator, Venaguide™ Seldinger wire-advancing device with Nitinol guidewire, Floswitch® revolving Luer-Lok™ system, scalpel, luer lock plugs and secondary fixation device.

##### Quad lumen Careflow™ without Floswitch® set contents:

Careflow™ quad lumen Polyurethane catheter with clear extensions, cannula introducer, steel needle introducer, vessel dilator, Venaguide™ Seldinger wire-advancing device with Nitinol guidewire, extension slide clamps, 5ml syringe, scalpel, luer lock plugs and secondary fixation device.

#### CAREFLOW™ QUINT LUMEN SETS

CODE WITHOUT FLOSWITCH®	CODE WITH FLOSWITCH®	EXTERNAL DIAMETER	LENGTH	CANNULA INTRODUCER		FLOSWITCH® INTRODUCER		NEEDLE INTRODUCER		GUIDEWIRE OD & LENGTH	BOX/CASE CONTENT (UNITS)	FLOW RATE (ML/MIN)				
				GAUGE	LENGTH	GAUGE	LENGTH	GAUGE	LENGTH			D	M1	M2	P1	P2
681735	681736	9.5Fr (14G, 16G, 16G, 18G, 18G)	150mm	16G	70mm	15G	70mm	17G	70mm	0.9 x 450mm	5/50	73	61	65	19	21
681737	681738	9.5Fr (14G, 16G, 16G, 18G, 18G)	200mm	16G	70mm	15G	70mm	17G	70mm	0.9 x 700mm	5/50	68	58	60	13	12

#### PRODUCT CONTENTS

##### Quint lumen Careflow™ with Floswitch® set contents:

Careflow™ quint lumen Polyurethane catheter with clear extensions, steel needle introducer, Floswitch® overneedle introducer, vessel dilator, Venaguide™ Seldinger wire-advancing device with Nitinol guidewire, Floswitch® revolving Luer-Lok system, scalpel, luer lock plugs and secondary fixation device.

##### Quint lumen Careflow™ without Floswitch® set contents:

Careflow™ quint lumen Polyurethane catheter with clear extensions, cannula introducer, steel needle introducer, vessel dilator, Venaguide™ Seldinger wire-advancing device with Nitinol guidewire, extension slide clamps, 5ml syringe, scalpel, luer lock plugs and secondary fixation device.

## CENTRAL VENOUS CATHETER SETS



# Careflow™ Central Venous Catheter Sets with Guiding Syringe

Single & Multi Lumen

Sterile, single use

PRODUCT CODE	EXTERNAL DIAMETER	LENGTH	CANNULA INTRODUCER		NEEDLE INTRODUCER		GUIDEWIRE OD & LENGTH	SECONDARY FIXATION DEVICE (SFD)	BOX/CASE CONTENT (UNITS)	FLOW RATE (ML/MIN)					WITH GUIDING SYRINGE
			GAUGE	LENGTH	GAUGE	LENGTH				D	M1	M2	P1	P2	
<b>CAREFLOW™ SINGLE LUMEN SETS WITH GUIDING SYRINGE</b>															
681788	5Fr (16G)	150mm	16G	70mm	17G	70mm	0.9 x 450mm	Yes	5/50	65					Yes
681789	5Fr (16G)	200mm	16G	70mm	17G	70mm	0.9 x 600mm	Yes	5/50	60					Yes
681790	5Fr (16G)	300mm	16G	70mm	17G	70mm	0.9 x 700mm	Yes	5/50	46					Yes
681791	7Fr (14G)	150mm	16G	70mm	17G	70mm	0.9 x 500mm	Yes	5/50	159					Yes
681792	7Fr (14G)	200mm	16G	70mm	17G	70mm	0.9 x 600mm	Yes	5/50	149					Yes
681793	7Fr (14G)	300mm	16G	70mm	17G	70mm	0.9 x 700mm	Yes	5/50	136					Yes
<b>CAREFLOW™ DOUBLE LUMEN SETS WITH GUIDING SYRINGE</b>															
681781	5Fr (18G, 20G)	60mm	18G	60mm	17G	70mm	0.6 x 450mm	Yes	5/50	40			10		Yes
681782	5Fr (18G, 20G)	150mm	18G	60mm	17G	70mm	0.6 x 450mm	Yes	5/50	29			6		Yes
681783	5Fr (18G, 20G)	200mm	18G	60mm	17G	70mm	0.6 x 650mm	Yes	5/50	30			6		Yes
681771	7Fr (16G, 18G)	150mm	16G	70mm	17G	70mm	0.9 x 500mm	Yes	5/50	74			24		Yes
681772	7Fr (16G, 18G)	200mm	16G	70mm	17G	70mm	0.9 x 600mm	Yes	5/50	73			21		Yes
681773	7Fr (16G, 18G)	300mm	16G	70mm	17G	70mm	0.9 x 700mm	Yes	5/50	61			19		Yes
<b>CAREFLOW™ TRIPLE LUMEN SETS WITH GUIDING SYRINGE</b>															
681784	5Fr (18G, 20G, 20G)	60mm	18G	60mm	17G	70mm	0.6 x 450mm	Yes	5/50	40	10		12		Yes
681785	5Fr (18G, 20G, 20G)	100mm	18G	60mm	17G	70mm	0.6 x 450mm	Yes	5/50	32	7		8		Yes
681786	5Fr (18G, 20G, 20G)	150mm	18G	60mm	17G	70mm	0.6 x 450mm	Yes	5/50	29	5		6		Yes
681787	5Fr (18G, 20G, 20G)	200mm	18G	60mm	17G	70mm	0.6 x 650mm	Yes	5/50	26	5		5		Yes
681774	7Fr (16G, 18G, 18G)	150mm	16G	70mm	17G	70mm	0.9 x 500mm	Yes	5/50	66	22		23		Yes
681775	7Fr (16G, 18G, 18G)	200mm	16G	70mm	17G	70mm	0.9 x 600mm	Yes	5/50	59	16		17		Yes
681776	7Fr (16G, 18G, 18G)	300mm	16G	70mm	17G	70mm	0.9 x 700mm	Yes	5/50	53	11		9		Yes
<b>CAREFLOW™ QUAD LUMEN SETS WITH GUIDING SYRINGE</b>															
681777	8.5Fr (14G, 16G, 16G, 18G)	150mm	16G	70mm	17G	70mm	0.9 x 500mm	Yes	5/50	64	82	57	20		Yes
681778	8.5Fr (14G, 16G, 16G, 18G)	200mm	16G	70mm	17G	70mm	0.9 x 600mm	Yes	5/50	58	76	45	19		Yes
<b>CAREFLOW™ QUINT LUMEN SETS WITH GUIDING SYRINGE</b>															
681779	9.5Fr (14G, 16G, 16G, 18G, 18G)	150mm	16G	70mm	17G	70mm	0.9 x 500mm	Yes	5/50	73	61	65	19	21	Yes
681780	9.5Fr (14G, 16G, 16G, 18G, 18G)	200mm	16G	70mm	17G	70mm	0.9 x 600mm	Yes	5/50	68	58	60	13	12	Yes

### PRODUCT CONTENTS

#### Single and Multi Lumen Careflow™ with Guiding Syringe set contents:

Careflow™ Polyurethane catheter with clear extensions and pinch clamps, cannula introducer, steel needle introducer, vessel dilator, Venaguide™ Seldinger wire-advancing device with Nitinol guidewire, 5ml Guiding Syringe, luer lock plugs and secondary fixation device.



## Hydrocath Assure™ Central Venous Catheter Sets

Antimicrobial-Coated  
Sterile, single use

Hydrocath Assure™ is a range of single, double and triple lumen polyurethane central venous catheters, designed for catheterisation by the Seldinger technique.

Every catheter is coated with a hydrophilic substance activated by biological fluids, and is designed to facilitate catheter insertion.

In addition, the Hydromer hydrophilic matrix is impregnated with the antimicrobial agent Benzalkonium Chloride (BZC).

CODE WITHOUT FLOSWITCH®	CODE WITH FLOSWITCH®	EXTERNAL DIAMETER	LENGTH	CANNULA INTRODUCER		FLOSWITCH® INTRODUCER		NEEDLE INTRODUCER		GUIDEWIRE OD & LENGTH	SECONDARY FIXATION (SFD)	BOX/CASE CONTENT (UNITS)	FLOW RATE (ML/MIN)		
				GAUGE	LENGTH	GAUGE	LENGTH	GAUGE	LENGTH				D	M2	P2
<b>HYDROCATH ASSURE™ SINGLE LUMEN SETS</b>															
681161		5Fr (16G)	200mm	16G	70mm	N/A	N/A	17G	70mm	0.9 x 700mm		5/50	56		
<b>HYDROCATH ASSURE™ DOUBLE LUMEN SETS</b>															
681152	681150	7Fr (16G, 18G)	200mm	16G	70mm	15G	70mm	17G	70mm	0.9 x 700mm	Yes	5/50	69		22
681153	681151	7Fr (16G, 18G)	300mm	16G	70mm	15G	70mm	17G	70mm	0.9 x 700mm	Yes	5/50	61		17
<b>HYDROCATH ASSURE™ TRIPLE LUMEN SETS</b>															
681157	681154	7Fr (16G, 18G, 18G)	150mm	16G	70mm	15G	70mm	17G	70mm	0.9 x 450mm	Yes	5/50	68	22	24
681158	681155	7Fr (16G, 18G, 18G)	200mm	16G	70mm	15G	70mm	17G	70mm	0.9 x 700mm	Yes	5/50	58	18	19
681159	681156	7Fr (16G, 18G, 18G)	300mm	16G	70mm	15G	70mm	17G	70mm	0.9 x 700mm	Yes	5/50	50	13	13

### PRODUCT CONTENTS

**All of the above sets contain:**

All sets without Floswitch® contain a cannula introducer, needle introducer, sliding clamps and 5ml syringe.

All Floswitch® sets include an Floswitch® revolving Luer-Lok™ system, Floswitch® cannula introducer and needle introducer.

Moss MA, Tebbs SE, Faroqui MH, Herbst T, Isaac JL, Brown J, Elliott TSJ, "A central venous catheter coated with Benzalkonium Chloride for the prevention of catheter-related microbial colonisation". Eur J Anaesth 2000, 17, 680-687.

CENTRAL VENOUS CATHETERS



## Hydrocath Assure™ Central Venous Catheter Sets with Guiding Syringe

Antimicrobial-Coated

Sterile, single use

Hydrocath Assure™ central venous catheter sets are also available with a guiding syringe, designed for catheterisation by the Seldinger technique.

Every catheter is coated with a hydrophilic substance activated by biological fluids, and is designed to facilitate catheter insertion.

In addition, the Hydromer hydrophilic matrix is impregnated with the antimicrobial agent Benzalkonium Chloride (BZC).

PRODUCT CODE	EXTERNAL DIAMETER	LENGTH	CANNULA INTRODUCER		NEEDLE INTRODUCER		GUIDEWIRE OD & LENGTH	SECONDARY FIXATION (SFD)	BOX/CASE CONTENT (UNITS)	FLOW RATE (ML/MIN)			WITH GUIDING SYRINGE
			GAUGE	LENGTH	GAUGE	LENGTH				D	M2	P2	
<b>HYDROCATH ASSURE™ DOUBLE LUMEN SETS WITH GUIDING SYRINGE</b>													
681924	7Fr (16G, 18G)	150mm	16G	70mm	17G	70mm	0.9 x 500mm	Yes	5/50	71		23	Yes
681925	7Fr (16G, 18G)	200mm	16G	70mm	17G	70mm	0.9 x 600mm	Yes	5/50	69		22	Yes
681926	7Fr (16G, 18G)	300mm	16G	70mm	17G	70mm	0.9 x 700mm	Yes	5/50	61		17	Yes
<b>HYDROCATH ASSURE™ TRIPLE LUMEN SETS WITH GUIDING SYRINGE</b>													
681934	7Fr (16G, 18G, 18G)	150mm	16G	70mm	17G	70mm	0.9 x 500mm	Yes	5/50	68	22	24	Yes
681935	7Fr (16G, 18G, 18G)	200mm	16G	70mm	17G	70mm	0.9 x 600mm	Yes	5/50	58	18	19	Yes
681936	7Fr (16G, 18G, 18G)	300mm	16G	70mm	17G	70mm	0.9 x 700mm	Yes	5/50	50	13	13	Yes

### PRODUCT CONTENTS

**Multi Lumen Hydrocath Assure™ with Guiding Syringe set contents:**

Hydrocath Assure™ BZC impregnated and hydrophilic coated Polyurethane catheter with clear extensions and pinch clamps, cannula introducer, steel needle introducer, vessel dilator, Venaguide™ Seldinger wire-advancing device with Nitinol guidewire, 5ml Guiding Syringe, luer lock plugs and secondary fixation device.

Moss MA, Tebbs SE, Faroqui MH, Herbst T, Isaac JL, Brown J, Elliott TSJ, "A central venous catheter coated with Benzalkonium Chloride for the prevention of catheter- related microbial colonisation". Eur J Anaesth 2000, 17, 680-687.



# CATHETERS AND HAEMOSTASIS

---

46 Arterial Catheters

47 Haemostasis

## ARTERIAL CATHETERS



## Careflow™ Arterial Catheters

Sterile, single use

Merit Careflow™ Polyurethane catheters are designed for access to radial or femoral arteries, for the determination of blood parameters or monitoring blood pressure.

### CAREFLOW™ ARTERIAL CATHETER SETS

PRODUCT CODE	PRODUCT	CATHETER		NEEDLE INTRODUCER		GUIDEWIRE OD & LENGTH	BOX/CASE CONTENT (UNITS)	PRIMING VOLUME (ML)	FLOW RATE (ML/MIN)
		EXTERNAL DIAMETER	LENGTH	GAUGE	LENGTH				
681631	Radial artery set	2.5Fr (22G)	60mm	22G	40mm	0.5 x 250mm	10/100	0.08	12.4
681646	Radial artery set	3Fr (20G)	60mm	22G	40mm	0.5 x 250mm	10/100	0.08	24.0
681642	Femoral artery set	3Fr (20G)	200mm	22G	40mm	0.5 x 450mm	10/100	0.14	15.4
681648	Femoral artery set	4Fr (18G)	200mm	17G	70mm	0.6 x 450mm	10/100	0.24	35.8

Merit Careflow™ Polyurethane catheter set contents:

Merit Careflow™ radiopaque polyurethane catheter, introducer needle, metal guidewire and Luer Lok™ plug.



## PreludeSYNC™ Radial Compression Device

The Next Generation of Radial Compression  
Sterile, single use

### ***Not just a product, an experience.***

The PreludeSYNC™ is a next generation Radial Compression Device that is personal, effective, comfortable, and secure. The PreludeSYNC™ was designed to address the market need for improved patient comfort without compromising safety. This innovative radial compression device comes packaged with creative designs printed directly on the bands to personalise the patient experience.

COMPRESSION DEVICE WITH SYRINGE	PRELUDESYNC™ COMPRESSION DEVICE SIZE/DESCRIPTION	UNITS/BOX
SRB24AC	24 cm 20 mL (Specialised Connection)	5
SRB29AC	29 cm 20 mL (Specialised Connection)	5
STAND ALONE SYRINGES	SYRINGE PRELUDESYNC™	UNITS/BOX
SRB20AC	20 mL (Specialised Connection)	10

## Clo-sur P.A.D.™

The Clo-Sur P.A.D.™ is a non-woven topical pad made with a proprietary polyprolate, a positively-charged material that reacts with negatively-charged red blood cells to accelerate clotting. The Clo-Sur P.A.D.™ promotes haemostasis and the Clo-Sur PLUS P.A.D.™ promotes haemostasis and is indicated as an antimicrobial barrier.

### Catalogue

Number	Description	Quantity
MVP44	Clo-sur P.A.D.™	10 per box
SCV22P-AB	Clo-sur Plus P.A.D.™	10 per box
HPAD22	Clo-sur P.A.D.™ w/Tegaderm	10 per box



## SafeGuard® Pressure Assisted Device

The SafeGuard® pressure assisted device delivers passive compression and enables immediate pressure adjustment. It also maintains consistent pressure on the site when obtaining haemostasis with manual compression, and during patient recovery, positioning and transport. The SafeGuard® pressure assisted device reduces demands on the clinical staff, maximises valuable resources, and enhances patient comfort.



82050

82000

**Catalogue**

Number	Size	Description
82000	24 cm	SafeGuard® Pressure Assisted Device
82050	12 cm	SafeGuard® Pressure Assisted Device

Packaged 10 units per box.

## SafeGuard® Radial™ Compression Device

The SafeGuard® Radial™ compression device is a 26 cm long, self-adhesive wristband that provides focused compression and assists patent haemostasis at the conclusion of radial procedures.

**Catalogue**

Number	Description	Quantity
SGR26	SAFE GUARD® Radial™	10 per box

Band length = 26 cm







# FLUID MANAGEMENT


---

- 50** Manifolds
- 55** Stopcocks
- 56** Adapters
- 57** Fluid Administration Sets
- 60** Contrast Management Systems
- 64** Flow Control Switch

# Merit® Manifolds

Merit® Manifolds have been designed to meet numerous clinician preferences. The clear polycarbonate material allows easy visibility for de-bubbling and the handles have been designed for ease of torque. In addition, options include one to five ports, 200 psi or 500 psi, 6" (15 cm) extension tubing and rotating luer, with block body or half body styles. Packaged 25 units per box, 4 boxes per case.

## MANIFOLDS WITH RIGHT PORT ORIENTATION AND ROTATING LUER



### 200 PSI (14 BAR)


ON Handles	202HN-R	203HN-R	204HN-R	204BN-R	203BN-R	202BN-R
Wide Port	202LHWN-R	203LHWN-R	204LHWN-R			
OFF Handles	202HF-R	203HF-R	204HF-R	204BF-R	203BF-R	202BF-R
Wide Port	202LHWF-R	203LHWF-R	204LHWF-R			
360/OFF Handles				204BF-RT	203BF-RT	
360/ON Handles					203BN-RT	

### 500 PSI (34 BAR)

ON Handles	502HN-R	503HN-R	504HN-R	504BN-R	503BN-R	502BN-R	501BN-R
Wide Port	502LHWN-R	503LHWN-R	504LHWN-R				
360/ON Handles		503HN-RT					
OFF Handles	502HF-R	503HF-R	504HF-R	504BF-R	503BF-R	502BF-R	501BF-R
Wide Port	502LHWF-R	503LHWF-R	504LHWF-R				
360/OFF Handles	502HF-RT				503BF-RT		

## MANIFOLDS WITH LEFT PORT ORIENTATION AND ROTATING LUER

Packaged 25 units per box, 4 boxes per case.



### 200 PSI (14 BAR)


ON Handles			203BN-L
OFF Handles	<b>202HF-L</b>	<b>203HF-L</b>	203BF-L

### 500 PSI (34 BAR)

OFF Handles	<b>502HF-L</b>	<b>503HF-L</b>
-------------	----------------	----------------

**MANIFOLDS WITH ALL FEMALE PORTS**

Packaged 25 units per box, 4 boxes per case.




2 Port Half Body    3 Port Half Body    3 Port Block Body    4 Port Block Body    5 Port Block Body

**200 PSI (14 BAR)**

ON Handles			203BN-F		
OFF Handles	202HF-F	203HF-F			205BF-F
360/OFF Handles			203BF-FT	204BF-FT	205BF-FT

**500 PSI (34 BAR)**


OFF Handles			503BF-F		
-------------	--	--	---------	--	--



205BF-F

**MANIFOLDS WITH RIGHT PORT ORIENTATION, 6" (15 CM) TUBING EXTENSION AND ROTATING LUER**

Packaged 25 units per box, 4 boxes per case.



2 Port Half Body    3 Port Half Body    4 Port Half Body    4 Port Block Body    3 Port Block Body    2 Port Block Body

**200 PSI (14 BAR)**


ON Handles	202HN-RX	203HN-RX			
OFF Handles	202HF-RX	203HF-RX	204HF-RX		203BF-RX
360/OFF Handles					203BF-RXT

**500 PSI (34 BAR)**

ON Handles		503HN-RX			502BN-RX
OFF Handles	502LHWF-R	503HF-RX	504HF-RX	503BF-RX	502BF-RX
360/OFF Handles				503BF-RXT	

# Devos™ Manifolds

Devos™ Manifolds with check relief valve, right port orientation and rotating luer. TRAM® and TRAM-P™ manifold ports have the added benefit of an integrated transducer on the port closest to the rotating adapter with comfortable wide-spaced finger grips. Packaged 25 units per box, 4 boxes per case. TRAM® and TRAM-P™ manifolds packaged 10 units per box, 4 boxes per case.



**4 Port  
Half Body**

**3 Port  
Half Body**

**2 Port  
Half Body**

**4 Port  
TRAM**

**4 Port  
TRAM-P\***

**200 PSI (14 BAR)**

OFF Handles	2V4LHWF-R	2V3HF-R	2V2HF-R	2V4LTPHWF-R	2V4LTCHWF-R
Tubing Extension	2V4LHWF-RXT				
360/OFF Handles	2V4LHWF-RT	2V3HF-RT	2V2HF-RT	2V4LTPHWF-RT	2V4LTCHWF-RT


\*TRAM-P™ includes 48" (122 cm) Compensation Line and Clip



**2V3HF-R**

## TRAM<sup>®</sup> Manifolds with Integrated Transducer

The TRAM<sup>®</sup> Manifold with Integrated Transducer gives a single operator the convenience to de-bubble, zero and rezero within the sterile field. The TRAM<sup>®</sup> has optimised clinician comfort by reducing hand fatigue with low torque, easy-to-turn wide-spaced handles, even when utilising higher injection pressures. The Devos<sup>™</sup> TRAM<sup>®</sup> manifold has a check relief valve. Packaged 10 units per box, 4 boxes per case.



1 Port

2 Port

3 Port

4 Port

### 200 PSI (14 BAR)


ON Handles Wide				
Standard	201LTPHN-R		203LTPHN-R	
OFF Handles Wide		202LTPHWF-R	203LTPHWF-R	204LTPHWF-R
Devos <sup>™</sup>				2V4LTPHWF-R
Devos <sup>™</sup> with 360° handle				2V4LTPHWF-RT
Standard			203LTPHF-R	

### 500 PSI (34 BAR)

ON Handles Wide		502LTPHWN-R	503LTPHWN-R	504LTPHWN-R
Standard	501LTPHN-R		503LTPHN-R	
OFF Handles Wide		502LTPHWF-R	503LTPHWF-R	504LTPHWF-R
Standard			503LTPHF-R	

## Tram-P™ Manifolds with Integrated Transducer

The TRAM-P™ features a pressure transducer on the first port. By positioning the compensation line on an IV pole rather than on a table, the TRAM-P™ maintains a constant calibrated reading. With this configuration, the manifold can be moved to any position without affecting the accuracy of the pressure reading. Comfortable wide-spaced finger grip option. The Devos™ TRAM-P™ manifold has a check relief valve. Packaged 10 units per box, 4 boxes per case.




### 200 PSI (14 BAR)

ON Handles Wide		202LTCHWN-R	203LTCHWN-R	204LTCHWN-R
Standard	201LTCHN-R		203LTCHN-R	
OFF Handles Wide		202LTCHWF-R	203LTCHWF-R	204LTCHWF-R
Devos				22V4LTCHWF-R
Devos with 360° handle				2V4LTCHWF-RT
Standard			203LTCHF-R	

### 500 PSI (34 BAR)

ON Handles Wide		502LTCHWN-R	503LTCHWN-R	504LTCHWN-R
Standard	501LTCHN-R		503LTCHN-R	
OFF Handles Wide		502LTCHWF-R	503LTCHWF-R	504LTCHWF-R
Standard	501LTCHF-R		503LTCHF-R	




# Marquis® Series Stopcocks

Merit Medical's large, easy-grip handle makes a Marquis Stopcock effortless to grip and turn. Available in pressure ratings of 200, 500 and 1050 psi with easy-turning handles. Luer options include rotating, swivel nut, fixed and slip. Packaged 50 units per box, 4 boxes per case.



H1RC



**200 PSI (14 BAR)**


1-way (white handle)	S1SLC	S1FLC	S1SNC
----------------------	-------	-------	-------

**500 PSI (34 BAR)**

1-way (light blue handle)	M1SLC	M1FLC	M1SNC
---------------------------	-------	-------	-------

**1050 PSI (72 BAR)**

1-way (dark blue handle)		H1FLC	H1SNC	H1RC
--------------------------	--	-------	-------	------



**50 PSI (3.5 BAR)**

3-way Large Bore		L3SNC	
4-way Large Bore			L4SNC

**200 PSI (14 BAR)**

3-way (white handle)	S3FC	S3SLC	S3FLC	S3SNC
4-way (white handle)		S4SLC	S4FLC	S4SNC

**500 PSI (34 BAR)**

3-way (light blue handle)	M3SLC	M3FLC	M3SNC	M3RRC	M3RLC
3-way (red handle)				M3RRC-R	
4-way (light blue handle)		M4FLC	M4SNC	M4RRC	

**1050 PSI (72 BAR)**

3-way (dark blue handle)		H3FLC	H3SNC	H3RRC	H3RLC
4-way (dark blue handle)		H4FLC	H4SNC	H4RRC	

**1200 PSI (83 bar)**

3-way (red handle)				U3RRC
--------------------	--	--	--	-------

## Airless Rotating Adapters

Used as connectors between two devices. Luer rotates a full 360° preventing torque transfer between the two attaching devices. Enhances patient safety by minimising the possibility of trapped air. Packaged 25 units per box, 4 boxes per case.



Male-Female  
Rotating Adapter

Male-Male  
Rotating Adapter

---


<b>Catalogue Number</b>	MFA100	MMA100
-------------------------	--------	--------



## Fluid Administration Sets


Merit offers a large variety of products for fluid delivery, including simple, single-line spike systems, complex multiple-line drip chambers with bifurcated or trifurcated lines, and large and small bore tubing with roller clamps or slide clamps. Products are available with micro- or macro-drip chambers, and with or without a micro filter. Custom fluid administration sets available.

**FLUID ADMINISTRATION SET** with user-selectable vented macro-drip chamber and small bore tubing. Packaged 25 units per box, 4 boxes per case.



	48" Length (122 cm)	60" Length (152 cm)	72" Length (183 cm)	108" Length (274 cm)
<b>Catalogue Number for without Micro Filter</b>	FAS1048	FAS1060	FAS1072	FAS1108
<b>Catalogue Number for with Micro Filter</b>		FAS1160	FAS1172	

**FLUID ADMINISTRATION SET** with non-vented macro-drip chamber, micro filter and small bore tubing. Packaged 25 units per box, 4 boxes per case.




	60" Length (152 cm)	72" Length (183 cm)
<b>Catalogue Number</b>	FAS1260	FAS1272

**FLUID ADMINISTRATION SET** with user-selectable vented spike. Packaged 25 units per box, 4 boxes per case.



	48" Length (122 cm)	60" Length (152 cm)	72" Length (183 cm)	108" Length (274 cm)
<b>Catalogue Number for Small Bore Tubing</b>	FAS4048	FAS4060	FAS4072	FAS4108
<b>Catalogue Number for Large Bore Tubing</b>	FAS5048	FAS5060	FAS5072	FAS5108

**FLUID ADMINISTRATION SET** with spike and large bore tubing. Packaged 25 units per box, 4 boxes per case.



20" Length (51 cm)  
With single disk check valve  
and user-selectable vented spike




48" Length (122 cm)  
with pinch clamp  
and non-vented spike


<b>Catalogue Number</b>	K08-02181	FAS5248
-------------------------	-----------	---------

## FLUID ADMINISTRATION SETS

**FLUID ADMINISTRATION SPIKE** with user-selectable vent, 1" (2.5 cm) small bore tubing and port protectors. Packaged 25 units per box, 4 boxes per case.



Female luer



3-way Marquis® Stopcock

**Catalogue Number**

CAS100

CAS200

**FLUID ADMINISTRATION SET** with vented macro-drip chamber with ball, 72" large bore tubing. Single-patient use for fluid or contrast management. Packaged 25 units per box, 4 boxes per case.




72" Length (183 cm)

**Catalogue Number**

K08-02385


**FLUID MANAGEMENT TUBE** 20" (50 cm) small bore tubing with small bore male slip luer. Packaged 15 units per box, 4 boxes per case.



**Catalogue Number**

K08-01859

**FLUID ADMINISTRATION SET** with user-selectable vented macro-drip chamber, micro filter, roller clamps and small bore tubing. Packaged 25 units per box, 4 boxes per case.



Bifurcated  
72" Length (183 cm)


Trifurcated  
72" Length (183 cm)

**Catalogue Number**

FAS2572


FAS3072

**FLUID ADMINISTRATION SET** with user-selectable vented spike, roller clamps and small bore tubing. Packaged 25 units per box, 4 boxes per case.



Bifurcated  
72" Length (183 cm)


Catalogue Number FAS6072




FAS6072

## Pressurised Contrast Administration Systems

Features one-way piercing Sosa spike and large bore tubing. Packaged 15 units per box, 4 boxes per case.



72" Length (183 cm)



6.5" Length (17 cm)

Catalogue Number K08-00006

K08-00685A

# Contrast Management Systems


Spike with large bore tubing, bonded Marquis® Stopcock and vented sideport protectors. Packaged 25 units per box, 4 boxes per case.

**Catalogue Number**

---

K08-02438

Vented spike, large bore with user-selectable vent. 1" large bore tubing, 1-way stopcock with distal female luer.




**Catalogue Number**

---

K08-02384

Vented spike, large bore with user-selectable vent. 2" large bore tubing, 3-way stopcock with distal female luer.




**Catalogue Number**

---

K08-00417B

Vented spike, large bore with user-selectable vent. 6.5" large bore tubing, 3-way stopcock with distal female luer.




**Catalogue Number**

---

K08-00077B

Vented spike, large bore with user-selectable vent. 2" large bore tubing, 1-way stopcock with distal male luer.




**Catalogue Number**

---

K08-00078A

Vented spike, large bore with user-selectable vent. 2" large bore tubing, 3-way stopcock with distal male luer.




# Contrast Management Systems

**Catalogue Number**

---

K08-00577

Macro Drip Chamber with filter, large bore, with 6.5" large bore tubing, 3-way stopcock with distal female luer.



**Catalogue Number**

---

K08-01878B

Merit CT Contrast Management spike with 20" (50 cm) large bore tubing, large bore vented spike, dual check valve, and fixed male luer. Packaged 20 units per box, 4 boxes per case.




**Catalogue Number**

---

K08-00719

Vented spike, large bore with vent. 1" large bore tubing, 3-way stopcock with distal female luer.




## CT Transfer Sets

**Catalogue Number**

---

K08-01961B

Vented, valved spike with 20" (50 cm) small bore tubing and nexus Luer. Packaged 25/box, 4 boxes/case.




**Catalogue Number**

---

K08-02195A

Vented, valved spike with 20" (50 cm) small bore blue stripe tubing and nexus Luer. Packaged 25/box, 4 boxes/case.



# Miser® Contrast Management Systems

The Miser® Contrast Management System’s exclusive orange ball and burette design prevents air from entering fluid lines during angiography, ensuring patient safety. The dual one-way check valve’s microbial barrier blocks cross contamination. Automatic refill during each contrast withdrawal allows for quick and easy preparation.

## MISER CONTRAST MANAGEMENT SYSTEM

Includes 72" (183 cm) distal tubing, 2 one-way check valves, 30 mL burette, air-inlet valve. Configurations available with and without a distal stopcock and with large or small bore distal tubing. Packaged 10 units per box, 4 boxes per case.

**CATALOGUE NUMBER**

**DESCRIPTION**

MISER1	Large Bore Tubing without Stopcock
MISER2	Large Bore Tubing with Stopcock
MISER3	Small Bore Tubing without Stopcock
MISER4	Small Bore Tubing with Stopcock



## MISER® CONTRAST MANAGEMENT BURETTE SPIKE

Includes large bore piercing spike, 30 mL burette, air-inlet valve, large bore distal 3-way stopcock with female luers. MUST BE USED IN CONJUNCTION WITH K08-00422D (shown below). Packaged 10 units per box, 4 boxes per case.

**CATALOGUE NUMBER MISER5**



## MERIT IN-LINE CONTRAST MANAGEMENT SYSTEM

Includes 54" (137 cm) large bore tubing, 2 one-way check valves, standard pressure 3-way stopcock and fixed male luers. Packaged 10 units per box, 4 boxes per case.

**CATALOGUE NUMBER K08-00422D**



# Flow Control Switch

The Merit Medical Flow Control Switch gives you the confidence you need when stopping and starting flow during critical procedures. This feature-rich, integrated, one-handed, single-channel device is designed with clinician and patient safety in mind. Merit’s High Pressure Flow Control Switch can be used with any application that involves liquid flow and requires inline shut off capability.

- Audible switch ensures ease of use and confirmation of open/close
- High pressure capable of pressures up to 1200 psi in dynamic (ON) applications
- On/Off slider switch one-handed actuation
- Customisation options available


## CATALOGUE

NUMBER	DESCRIPTION	SLIDER COLOUR	BOX QUANTITY
U1SHC	Flow Control Switch	Blue	24
U1SHC-W	Flow Control Switch	White	24
U1SHC-R	Flow Control Switch	Red	24
U1SHC-Y	Flow Control Switch	Yellow	24
U1SHC-G	Flow Control Switch	Green	24

Pressure Rating: Open 1200 psi, Closed 150 psi

Lumen I.D.: 0.065 inches

Connectors: Male to Female







# TUBING

---

**66** High Pressure Tubing

**67** Pressure Monitoring Tubing


## TUBING

## High Pressure Tubing

Merit offers three specific material configurations ranging from highly flexible polyurethane to clear PVC to meet your specific pressure needs. Merit's high pressure tubing configurations offer both a fixed male connector and a rotating adapter connector.

### CONTRAST INJECTION TUBING, POLYURETHANE, FLEXIBLE BRAIDED


Rated at 1200 psi (83 bar). Packaged coiled, 25 units per box, 4 boxes per case.



Male to Female	10" Length (25 cm)	20" Length (51 cm)	30" Length (76 cm)	48" Length (122 cm)	60" Length (152 cm)	72" Length (183 cm)
Airless Rotating Adapter	HPF100E	HPF200E	HPF300E	HPF480E		HPF720E
With KEEP™ Tubing Clip		HPF200K	HPF300K	HPF480K	HPF600K	HPF720K
Fixed Male Connector	HPF101E	HPF201E	HPF301E	HPF481E		HPF721E
With KEEP Tubing Clip		HPF201K	HPF301K	HPF481K	HPF601K	HPF721K

### CONTRAST INJECTION TUBING, PVC, CLEAR

Rated at 1200, 900 and 500 psi (83, 62 and 34 bar). Packaged straight, 25 units per box, 4 boxes per case.




Male to Female	10" Length (25 cm)	20" Length (51 cm)	30" Length (76 cm)	48" Length (122 cm)	59" Length (150 cm)	72" Length (183 cm)
<b>Airless Rotating Adapter</b>						
1200 psi (83 bar)	HPT100E	HPT200E	HPT300E	HPT480E	HPT590E	
900 psi (62 bar)	HP9100E	HP9200E	HP9300E	HP9480E		HP9720E
500 psi (34 bar)	HP5100E	HP5200E	HP5300E	HP5480E		

#### Fixed Male Connector

1200 psi (83 bar)	HPT101E	HPT201E	HPT301E	HPT481E	HPT591E	
900 psi (62 bar)	HP9101E	HP9201E	HP9301E	HP9481E		HP9721E
500 psi (34 bar)	HP5101E	HP5201E	HP5301E	HP5481E		

### CONTRAST INJECTION TUBING, FLEXIBLE, NYLON / POLYURETHANE, EXCITE™

Rated at 1200 psi (83 bar). Packaged straight, 25 units per box, 4 boxes per case.




Male to Female	10" Length (25 cm)	20" Length (51 cm)	30" Length (76cm)	39" Length (100 cm)	48" Length (122 cm)	59" Length (150 cm)	72" Length (183 cm)
<b>Airless Rotating Adapter</b>	HPC100E	HPC200E	HPC300E	HPC390E	HPC480E	HPC590E	HPC720E
<b>Fixed Male Connector</b>	HPC101E	HPC201E	HPC301E	HPC391E	HPC481E	HPC591E	HPC721E


# Pressure Monitoring Tubing

PRESSURE MONITORING TUBING PVC, RIGID NON-EXPANDING


Packaged coiled, 25 units per box, 4 boxes per case.



<b>Male to Male</b>	6" Length (15 cm)	12" Length (30 cm)	24" Length (61 cm)	36" Length (91 cm)	48" Length (122 cm)	60" Length (152 cm)	72" Length (183 cm)
<b>Catalogue Number</b>	PM6006	PM6012	PM6024	PM6036	PM6048	PM6060	PM6072



<b>Male to Female</b>	6" Length (15 cm)	12" Length (30 cm)	24" Length (61 cm)	36" Length (91 cm)	48" Length (122 cm)
<b>Catalogue Number</b>	PM6106	PM6112	PM6124	PM6136	PM6148



<b>Male to Female</b>	60" Length (152 cm)	72" Length (183 cm)	84" Length (213 cm)	96" Length (244 cm)
<b>Catalogue Number</b>	PM6160	PM6172	PM6184	PM6196



# SAFETY & INFECTION CONTROL

---

- 69** DualCap® Disinfection & Protection System
- 70** Syringes
  - 70** Coronary Control Syringes
  - 72** Vaclok®
  - 73** Medallion®
- 76** Safety and Waste Management
  - 76** Futura® Safety Scalpel
  - 77** BackStop® & BackStop Plus™
  - 77** MiniStop® & MiniStop Plus™
  - 77** ShortStop® & ShortStop Plus® & Grandstand®
  - 79** DriPAD™
- 80** Pen & Labels



## DualCap®

### Disinfection & Protection System

Sterile, single use

The DualCap® System is a uniquely designed disinfecting cap system for use on both the male luer connector at the end of the IV tubing and the needle-free valve. The DualCap® System consists of two caps designed to help prevent intraluminal contamination, as well as device cross-contamination, especially important for IV catheter use over extended dwell times.

- Disinfects rapidly and automatically at both IV catheter points
- The only male luer connector cap that prevents IPA from entering the fluid path
- Locks in ongoing protection against intravascular infections
- Easy to grasp
- Durable IV Pole Strips™ designed to withstand handling and provide convenient bedside delivery
- Supports infection control goals




CATALOGUE NUMBER	DESCRIPTION	CAPS PER BOX	BOX PER CASE	TOTAL CAPS PER CASE
450-LB	DualCap® Solo Cap with foil lid for needle-free valve, packed in dispenser box.	350	12	4200
550-LB	DualCap® IV Pole Strips For needle-free valves. 35 strips/10 caps per strip.	350	12	4200
660-DB	DualCap® IV Pole Strips For male luer connectors. 35 strips/5 caps per strip.	175	6	1050

# CCS™ Coronary Control Syringes

The CCS™ Coronary Control Syringe has a clear barrel for clarity and smoothness and a solid plunger body to maintain stability and durability under pressure. The rotating and fixed luers allow flexibility and confidence in the luer connection. The specially designed safety space limits the potential for introducing air bubbles into the catheter. The CCS™ locking ring permits vacuum or negative pressure. The Inject8™ and Inject10™ offer clinicians more pressure with less applied force. Packaged 25 units per box, 4 boxes per case.

SAFETY & INFECTION CONTROL



- |                                 |                                   |                                     |                           |                                  |                           |                           |
|---------------------------------|-----------------------------------|-------------------------------------|---------------------------|----------------------------------|---------------------------|---------------------------|
| 6 mL<br>Locking Ring<br>Plunger | 8 mL<br>Solid Plunger<br>Inject8™ | 10 mL<br>Solid Plunger<br>Inject10™ | 10 mL<br>Solid<br>Plunger | 12 mL<br>Locking Ring<br>Plunger | 12 mL<br>Solid<br>Plunger | 20 mL<br>Solid<br>Plunger |
|---------------------------------|-----------------------------------|-------------------------------------|---------------------------|----------------------------------|---------------------------|---------------------------|

**Rotating Adapter**

CCS460	CCS880	CCX010	CCS800	CCS100	CCS600	CCS320
--------	--------	--------	--------	--------	--------	--------


**Fixed Male Connector**

CCS461	CCS881	CCX011	CCS801	CCS101	CCS601	CCS321
--------	--------	--------	--------	--------	--------	--------



ST800

**SMART TIP™ CONTROL SYRINGE** with blue silicone tip, ring plunger and ring grip. Packaged 25 units per box, 4 boxes per case.



- |                        |                        |                        |
|------------------------|------------------------|------------------------|
| 10 mL<br>Solid Plunger | 12 mL<br>Solid Plunger | 20 mL<br>Solid Plunger |
|------------------------|------------------------|------------------------|


**Rotating Adapter**

ST800	ST600	ST320
-------	-------	-------

**Fixed Male Connector**

ST801	ST601	ST321
-------	-------	-------


**CORONARY CONTROL SYRINGE** with palm pad plunger and ring grip. Packaged 25 units per box, 4 boxes per case.



10 mL Solid Plunger    10 mL Solid Plunger    12 mL Solid Plunger    12 mL Solid Plunger

**Rotating Adapter**    CCSB900    CCSB700

**Fixed Male Connector**    CCSB901    CCSB701



8 mL Solid Plunger Inject8™    8 mL Solid Plunger Inject8™    10 mL Solid Plunger Inject10™    10 mL Solid Plunger Inject10™


**Rotating Adapter**    CCSB880    CCXB010

**Fixed Male Connector**    CCSB881    CCXB011



CCSB900


**CORONARY CONTROL SYRINGE** with palm pad plunger and wing grip. Packaged 25 units per box, 4 boxes per case.



10 mL Solid Plunger    10 mL Solid Plunger    12 mL Solid Plunger    12 mL Solid Plunger

**Rotating Adapter**    CCSB910    CCSB710

**Fixed Male Connector**    CCSB911    CCSB711



8 mL Solid Plunger Inject8™    8 mL Solid Plunger Inject8™    12 mL Locking Ring Plunger    12 mL Locking Ring Plunger

**Rotating Adapter**    CCSW880    CCS200

**Fixed Male Connector**    CCSW881    CCS201

SAFETY & INFECTION CONTROL

# VacLok® Vacuum Pressure Syringes

The VacLok® vacuum pressure syringe comes with a clear polycarbonate barrel, white plunger, four- or six-position locking mechanism and fixed male luer. Packaged 20 or 25 units per box, 4 boxes per case.



**10 mL**  
four position  
25/box

**20 mL**  
four position  
20/box

**30 mL**  
six position  
20/box

**60 mL**  
six position  
25/box

**Catalogue Number**

VAC110

VAC120


VAC130

VAC160




# Medallion® Syringes

Merit Medical's Medallion® Syringes come in seven colour options with custom printing or a frosted field for writing information directly on the syringe barrel. The Medallion® Syringe features a clear polycarbonate barrel, coloured plunger, and fixed male Luer. Packaged 25 units per box, 4 boxes per case.



MSS011-LG



Plunger Colours

1 mL      3 mL      6 mL      10 mL      20 mL      30 mL      60 mL

White	MSS011	MSS031	MSS061	MSS111	MSS121	MSS131	MSS160
Light Blue	MSS011-LB	MSS031-LB	MSS061-LB	MSS111-LB	MSS121-LB	MSS131-LB	
Yellow	MSS011-Y	MSS031-Y	MSS061-Y	MSS111-Y	MSS121-Y	MSS131-Y	
Red	MSS011-R	MSS031-R	MSS061-R	MSS111-R	MSS121-R	MSS131-R	
Light Green	MSS011-LG	MSS031-LG	MSS061-LG	MSS111-LG	MSS121-LG	MSS131-LG	
Dark Green	MSS011-DG	MSS031-DG	MSS061-DG	MSS111-DG	MSS121-DG	MSS131-DG	
Purple	MSS011-PR	MSS031-PR	MSS061-PR	MSS111-PR	MSS121-PR	MSS131-PR	

**MEDALLION® SYRINGES** with clear polycarbonate barrel, coloured plunger, black medication printing, and fixed male Luer. Packaged 25 units per box, 4 boxes per case.



Plunger Colours

NITRO (10 mL)      LIDOCAINE (10 mL)      HEPARIN (10 mL)      CONTRAST (10 mL)      SALINE (10 mL)

Light Blue	K01-60085	K01-60081	K01-60082	K01-60086	K01-60083
Yellow	K01-60101	K01-60097	K01-60098	K01-60102	K01-60099
Red	K01-60109	K01-60105	K01-60106	K01-60110	
White	K01-60077			K01-60078	K01-60075
Light Green	K01-60093	K01-60089	K01-60090	K01-60094	
Dark Green				K01-00786	

Plunger Colours

FENTANYL (3 mL)      XYLOCAINE (10 mL)      VERSED (3 mL)      FLUSH (10 mL)


Light Blue		K01-00917	K01-05185	K01-05093
Red		K01-00953		
Light Green	K01-03243			

SAFETY & INFECTION CONTROL

## SYRINGES

# Medallion® Sword handle Syringes

Medallion® Sword Handle Syringes feature sword handles, clear polycarbonate barrels, coloured plungers, and fixed male Luers. Packaged 25 units per box, 4 boxes per case.



Plunger Colours	10 mL	20 mL
White	MSSW11	MSSW21
Light Blue	MSSW11-LB	MSSW21-LB
Yellow	MSSW11-Y	MSSW21-Y
Red	MSSW11-R	MSSW21-R
Light Green	MSSW11-LG	MSSW21-LG
Dark Green	MSSW11-DG	MSSW21-DG
Purple	MSSW11-PR	MSSW21-PR




MSSW11-Y

# Custom Medallion® Syringes


Merit Medical's signature Medallion® Syringes can be custom ordered to your specifications. Medallion® Syringes come with a clear polycarbonate barrel, coloured plunger, customised medication printing, and Luer options including fixed male or slip Luer styles.

Select one option from each category and identify your custom medication to be printed on the barrel. Your Merit Representative will provide you with the corresponding catalogue number.

Custom Medallion® Syringes are packaged 25 units per box, 4 boxes per case.



## SIZE OPTIONS



## PLUNGER COLOUR OPTIONS:

- White
- Light Blue
- Light Green
- Dark Green
- Yellow
- Red
- Purple

## PRINT COLOUR OPTIONS:

- Black
- White

## LUER OPTIONS:


- Fixed Male Luer
- Slip Luer

SAFETY & INFECTION CONTROL

# Futura® Safety Scalpels

These retractable safety scalpels are designed to protect clinicians and patients from accidental cuts during passing and disposal. The one-handed operation combines convenience with rapid response for any procedure. Packaged 50 units per box, 2 boxes per case.

SAFETY & INFECTION CONTROL



Blade Style #11, Blue  
SMS211



Blade Style #15, Yellow  
SMS215



Blade Style #10, Green  
SMS210

Catalogue Number	Description	Colour
SMS210	#10 Blade Style	green
SMS211	#11 Blade Style	blue
SMS215	#15 Blade Style	yellow

## BACKSTOP® AND MINISTOP® DISPOSAL SYSTEMS

For quick and safe fluid waste elimination. Includes absorbent material and particulate-free, die-cut sponge and seal-tight lid.

Catalogue Number	Description
MBS100	BackStop 500 mL fluid capacity, packaged 10 per box
MMS100	MiniStop 100 mL fluid capacity, packaged 20 per box
MMSF200	MiniStop filter, packaged 20 per box



MBS100

## BACKSTOP+™ AND MINISTOP+™ DISPOSAL SYSTEMS

For quick and safe fluid waste elimination. Has a seal-tight lid to prevent fluid leakage after a procedure and includes a temporary sharps holder.

Catalogue Number	Description
MBSP100	BackStop+ 500 mL fluid capacity, packaged 10 per box
MMSP100	MiniStop+ 100 mL fluid capacity, packaged 20 per box



MBSP100

## DUGOUT® CLOSED WASTE BASINS

With 1000 mL fluid capacity and absorbent material. Includes particulate-free, die-cut sponge, secondary compartment for storage of ancillary products or saline, and seal-tight lid. Packaged 6 units per box, 4 boxes per case.

Catalogue Number	Description
DUG100	DugOut 1000 mL fluid capacity
DUG1SS	DugOut 1000 mL fluid capacity with ShortStop® Temporary Sharps Holder



DUG100

## SHORTSTOP® TEMPORARY SHARPS HOLDER

Reduce the possibility of accidental sharps injuries. Rigid, puncture-resistant material reduces potential of needle sticks and increases safety during disposal. Coring-resistant foam provides easy removal and minimises particulate.

Catalogue Number SS10

Packaged 25 per box, 4 boxes per case.



SS10

## SHORTSTOP ADVANTAGE® TEMPORARY SHARPS HOLDER

The ShortStop Advantage® temporary sharps holder features a new side access to rest a fully loaded syringe attached to a needle.

**Catalogue Number** SS20

Packaged 25 per box, 4 boxes per case.



## SHORTSTOP GRANDSTAND® TEMPORARY SHARPS HOLDER

The ShortStop Grandstand® holds sharp devices upright for convenient access and safety. The compact configuration is a space-saving solution for the safe storage of needles, syringes and scalpels.

**Catalogue Number** SS30

Packaged 10 per box, 4 boxes per case.



# DriPAD™ Super Absorbent Pad

## SUPER ABSORBENT TECHNOLOGY THAT HELPS YOU MANAGE THE RISK OF BIOHAZARDS

The DriPAD™ can be used underneath patients to protect linens and beds. It can also be conveniently placed on the floor to absorb spills, trapping fluids inside.

### Super Absorbent Technology

- DriPADs™ absorb up to 35x their weight in fluids,\* making them more efficient than the leading absorbent pad on the market. That's up to 10x more absorption power than alternative pads.
- A specially treated absorbent core draws fluid away from the surface and traps it in the interior, so that the surface remains dry.

### Containment & Odor Control


- Each DriPAD™ features a robust polypropylene backing to provide a barrier that prevents leakage.

## ORDERING INFORMATION


Part Number	Description	Quantity
<b>Sterile</b>		
<b>PAD0909E</b>	9" x 9" Soaker Pad; absorbs and solidifies .47 liters (16oz.) of fluid	10 units
<b>PAD2320E</b>	23" x 20" Utility Pad; absorbs and solidifies 1.89 liters (64oz.) of fluid	5 units
<b>PAD2332E</b>	23" x 32" Large Pad; absorbs and solidifies 2.84 liters (96oz.) of fluid	5 units
<b>PAD2360E</b>	23" x 60" Stretcher-Size Gurney Pad; absorbs and solidifies 5.68 liters (192oz.) of fluid	5 units

Part Number	Description	Quantity
<b>Non-sterile</b>		
<b>PAD0909NS</b>	9" x 9" Soaker Pad; absorbs and solidifies .47 liters (16oz.) of fluid	10 units
<b>PAD2320NS</b>	23" x 20" Utility Pad; absorbs and solidifies 1.89 liters (64oz.) of fluid	5 units
<b>PAD2332NS</b>	23" x 32" Large Pad; absorbs and solidifies 2.84 liters (96oz.) of fluid	5 units
<b>PAD2360NS</b>	23" x 60" Stretcher-Size Gurney Pad; absorbs and solidifies 5.68 liters (192oz.) of fluid	5 units

### Available in 4 sizes



### Absorption Capacities



\*Data on file



## PAL™ - PEN AND LABELS

Strategically designed medication labeling system. The labels are made of moisture resistant polypropylene material. In addition, each PAL includes a sterile, smear resistant, permanent marker. Packaged 25 units per box, 4 boxes per case.


### Catalogue

Number	Description
PAL100	Contains pen and label set of 35
PAL200	Contains pen and label set of 24
PAL225	Contains pen, label set of 24 and guide wire bowl label set of 5
PAL250	Contains pen, label set of 24 and large bowl label set of 7
PAL300	Contains pen and label set of 12
PAL325	Contains pen, label set of 12 and guide wire bowl label set of 5
PAL350	Contains pen, label set of 12 and large bowl label set of 7
PAL400	Contains pen and blank label set of 12
PAL425	Contains pen, blank label set of 12 and guide wire bowl label set of 5
PAL450	Contains pen, blank label set of 12 and large bowl label set of 7


SAFETY & INFECTION CONTROL



### PAL 200 SERIES



### PAL 300 SERIES




### PAL 400 SERIES



# PAL™ - Pen & Labels

## PAL500


Contains pen and guide wire bowl label set of 5



Catalogue Number PAL500

## PAL600



Contains pen and large bowl label set of 7



Catalogue Number PAL600

## PAL700

Contains pen, label set of 35 and large bowl label set of 7



Catalogue Number PAL700

All label pictures are not to scale.

## PAL™ PEN

Sterile, smear resistant, permanent marker. Packaged 25 units per box, 4 boxes per case.



Catalogue Number K13-00476

SAFETY & INFECTION CONTROL

# Custom PAL™ Labels



Customers now have the option to have a custom label set created to meet their specific clinical needs.


Comply with the Joint Commission’s National Patient Safety Goals\* using Merit Medical’s strategically designed custom medication labeling system.

Quickly and easily order all your custom medication labels with Merit’s revolutionary iPad app. Through the app, you’ll be able to select critical drug information, including drug name, dosing, date and signatures while maintaining complete compliance with the Joint Commission.

Contact your Merit Medical Sales Representative today for all your medication labeling needs.

## STEP 1: CHOOSE A TEMPLATE

Label Set Template options for syringes, cups and small bowls:




Label Set of 12

Label Set of 24

Label Set of 35

Label Set Template options for large bowls:



Label Set of 5 (Rectangle Design)

Label Set of 7 (Curved Label Design)


## STEP 2: IDENTIFY THE MEDICATION NAMES YOU WOULD LIKE TO HAVE PRINTED ON EACH LABEL

## STEP 3: CHOOSE THE APPROPRIATE COLOURS FROM THE OPTIONS LISTED BELOW

### Colour Options:

Choose Label Color

- Red
- Dark Blue
- Blue
- Aqua
- Lt. Green
- Green
- Yellow
- Purple
- Pink
- Lt. Orange
- Orange
- White
- Transparent
- Gray
- Brown



Your Merit Representative will provide you with the corresponding Catalogue number.

Most custom label sets are packaged 25 units per box, 4 boxes per case.

To better meet your clinical needs, each custom label set kit can contain a sterile permanent marker. The sterile permanent marker is also available stand alone.

NOTE: Custom label sets are available upon special request. Please allow up to six weeks for delivery.

\* When label is completed at the time the medication or solution is prepared.



# COMPLETE PROCEDURAL SOLUTIONS

---

**84** Merit Value Packs™

**89** Radial Drapes

**90** Scatter Radiation Drapes

# Angiography Tray 1

Includes


Quantity	Description
2	Band Bags 36" x 48" (91 x 122 cm)
1	Set Up Cover 44" x 77" (112 x 196 cm)
1	Dome Band Bag 22" Deep (56 cm)
1	Needle 18 G x 1.5" (40 mm)
1	Needle 21 G x 1.5" (40 mm)
1	Needle 22 G x 1.5" (40 mm)
1	Needle 25 G x 1.5" (40 mm)
1	#11 Futura® Safety Scalpel
1	Medicine Cup 2oz, 60 mL
1	Bowl 16oz, 500 mL
1	Bowl 32oz, 1000 mL
1	Guide Wire Bowl 2000 mL
1	Tray 9" x 5" x 2" (23 x 12 x 5 cm)
1	Curved Forceps 5" (13 cm)
50	4" x 4" (10 x 10 cm) Gauze Pads (12 Ply)
6	Cloth Towels (Blue)
4	Towel Clamps
3	Sponge Sticks 8" (20 cm)
2	XL NRF Gowns
1	Back Table Cover 50" x 90" (127 x 229 cm)
1	Femoral Drape 86" x 135" (206 x 343 cm) with 2 Windows
2	Line Identification Labels (Red/Yellow/Blue/White)
1	Merit Disposal Depot™ Waste Bag (MDD200HP)
1	Angiographic Needle 18 G x 2 3/4" (7 cm) Thin Wall
3	Medallion® Syringes 20 mL L/L, White, Saline
2	Medallion® Syringes 10 mL L/L, Lt Blue, Lidocaine
2	Medallion® Syringes 10 mL L/L, Yellow, Contrast
1	ShortStop® Temporary Sharps Holder
1	Manifold Setup With OFF Handles*

**\*MANIFOLD SETUP CONTAINS:**

- 203HN-R - 3-Port Manifold
- 1st Port - PM6172 - Pressure Monitoring Tubing
- 2nd Port - FAS5072 - Tubing w/ Large Bore Spike
- 3rd Port - FAS1272 - Tubing w/ Non-Vented Macro Drip Chamber, Male Luer
- CCX010 - Coronary Control Syringe
- M3SNC - 3-Way Stopcock

Packaged 4 per case.

**Catalogue Number** K09T-11906



# Custom Angiography Kits

When your clinical needs are unique, Merit Medical can customise an angiography kit to meet your specific requirements. A Merit Sales Representative can assist you in combining Merit components into a single kit.

This drawing is property of Merit Medical Systems Inc. Reproduction or use without written authorization is prohibited.

REVISIONS				
REV	ECO	DESCRIPTION	DATE	APPROVED
A	KO449	DELETED MEDALLIONS, 500PSI	02-14-94	


TOLERANCE	SIGNATURES	DATE	MERIT MEDICAL SYSTEM INC. 79 W. 4500 SOUTH SALT LAKE CITY, UTAH	
FRAC. ±	DRAWN BY: K. DENNING		HOSPITAL:	
X.XX ±	SALES REP.: B. WOELPPER		KIT NAME: CUSTOM CATH LAB KIT	
X.XXX ±	CUSTOMER:		SIZE A	DRAWING NUMBER K09-X X X X
ANGLE ±	MFG. ENG.:		CUSTOMER REV A	
DO NOT SCALE BREAK ALL EDGES	QA APPROVAL T. MYRUP		SCALE NONE	SHEET 1 OF 1
				DWG REV A

NOTE: All custom kits require a signed drawing (example above) and contract.

# Pericardiocentesis Kits

## 8.3F PERICARDIOCENTESIS PROCEDURE KIT

Designed to assist the clinician in draining fluid from the pericardial sac.



Includes

Quantity	Description	Catalogue Number	Description
1	8.3F Pigtail or straight catheter, 40 cm with one marker band	PC801	Pigtail Shape
1	8.5F Dilator 22 cm length	PC802	Straight Shape
1	Guide wire: 0.035" (0.89 mm) x 80 cm, 3 mm J, double ended standard PTFE	PC803EU	Pigtail Shape; kit does not include ChloraPrep®
1		PC804EU	Straight Shape; kit does not include ChloraPrep®
2	3-way Stopcocks large bore		
1	18G x 1 1/2" (3.8 cm) Medication Needle		
1	25G x 1 1/2" (3.8 cm) Medication Needle		
1	18G x 3" (9 cm) Access needle		
1	18G x 5.9" (15 cm) Chiba Access needle		
1	10 mL Luer lock syringe		
1	20 mL Luer lock syringe		
1	60 mL Luer lock syringe		
1	1000 mL Merit Disposal Depot™ waste bag with drainage spout		
10	4" x 4" (10 cm x 10 cm) Gauze sponges		
1	61 cm x 61 cm Fenestrated drape		
1	46 cm x 66 cm Absorbent barrier drape		
1	Double-ended alligator clip cable		
1	CSR Drape		
1	#11 Futura® Safety Scalpel		
1	3.0 silk suture with curved needle		
1	Medication label set		

Pericardiocentesis Kits PC801 and PC802 do not have the CE Marking. ChloraPrep® is a registered trademark of CareFusion 2200, Inc.

# Pericardiocentesis Kits

## 6F PERICARDIOCENTESIS PROCEDURE KIT

Designed to assist the clinician in draining fluid from the pericardial sac.



Includes

Quantity	Description
1	6F High-flow pigtail or straight catheter 60 cm
1	6F Dilator 22 cm length
1	3-way stopcock
1	Guide wire 0.035" (0.89 mm) x 80 cm
1	Double ended Heavy Duty PTFE (J & straight)
1	22G x 2 3/8" (6 cm) Medication Needle
1	19G x 1 1/2" (3.8 cm) Medication Needle
1	18G x 3" (9 cm) Access Needle
1	18G x 5.9" (15 cm) Chiba Access Needle
1	10 mL Luer lock syringe
1	60 mL Luer lock syringe
1	1400 Merit Disposal Depot™ waste bag with connecting tube
10	4" x 4" (10 cm x 10 cm) Gauze sponges
1	40 cm x 40 cm Fenestrated drape
1	46 cm x 66 cm Absorbent barrier drape
1	Double-ended alligator clip cable
1	CSR wrap
1	#11 Futura® Safety Scalpel

Packaged 4 kits per box.

Catalogue Number	Description
PC101	6F High-flow pigtail catheter
PC201	6F High-flow straight catheter


COMPLETE PROCEDURAL SOLUTIONS

## Rapid Response Right Heart Kits

Right Heart Kits are pre-assembled and include necessary pressure monitoring and fluid management components.

### RAPID RESPONSE RIGHT HEART ON KIT


2-port manifold, ON handle with transducer, 60" (152 cm) pressure monitoring tubing, and 60" (152 cm) non-vented small bore tubing flush line. Packaged 10 units per box, 4 boxes per case.



**Catalogue Number** K09-10471

### RAPID RESPONSE RIGHT HEART OFF KIT

2-port manifold, OFF handle with transducer, 60" (152 cm) pressure monitoring tubing, and 60" (152 cm) non-vented small bore tubing flush line. Packaged 10 units per box, 4 boxes per case.



**Catalogue Number** K09-10472


Rapid Response Kits do not have CE Marking.




# Radial Drapes

The 4-Hole Drapes are standalone convenience products for radial procedures that Merit Medical has made available for clinical settings where a fenestrated drape outside of the procedure pack is necessary or desired. Packaged 10 per case.


Catalogue Number	Description
K12T-04696	Femoral Radial Angiography Drape, Angled 4-hole
K12T-05047	Standard 4-hole drape
K12T-07173	Femoral Radial Angiography Drape



K12T-04696



K12T-07173



K12T-05047


COMPLETE PROCEDURAL SOLUTIONS

## Drape Armour® - Scatter Radiation Drapes


These positional protection shields are designed to protect at the procedure site where scatter radiation is most prevalent. They can help protect surgeons and other members of the clinical staff, and offer substantial attenuation as demonstrated by direct beam studies. Packaged 10 per case.

**Catalogue**


Number	Description
K12T- 09374	41x 25cm
K12T- 09375	41x 41 cm
K12T- 09376	41x 35cm w/slit and circ hole



K12T- 09374



K12T- 09375



K12T- 09376

Drape Armour® is a registered trademark of Microtek Medical, Inc.



# ACCESSORIES

---

**92** Pressure Infusor Bags™


**93** Flow Control Switch

**94** Floswitch®

PRESSURE INFUSOR BAGS

# PIB - Pressure Infusor Bags™

Nylon bags for pressurisation of fluids. Dual-setting pressure relief valve easily switches between standard and overpressure mode with the click of a button (PIB500 and PIB1000 only).



Catalogue Numbers	Description
PIB500	500 mL 12 per box, 48 per case
PIB1000	1000 mL 12 per box, 48 per case
MTM330E	3000 mL 2 per box, 10 per case



PIB500

PIB1000

ACCESSORIES

# Flow Control Switch

The Merit Medical Flow Control Switch gives you the confidence you need when stopping and starting flow during critical procedures. This feature-rich, integrated, one-handed, single channel device is designed with clinician and patient safety in mind. Merit’s High Pressure Flow Control Switch can be used with any application that involves liquid flow and requires inline shut off capability.

- Audible switch ensures ease of use and confirmation of open/close
- High pressure capable of pressures up to 1200 psi in dynamic (ON) applications
- On/Off slider switch one-handed actuation
- Customisation options available


## CATALOGUE

NUMBER	DESCRIPTION	SLIDER COLOUR	BOX QUANTITY
U1SHC	Flow Control Switch	Blue	24
U1SHC-W	Flow Control Switch	White	24
U1SHC-R	Flow Control Switch	Red	24
U1SHC-Y	Flow Control Switch	Yellow	24
U1SHC-G	Flow Control Switch	Green	24

Pressure Rating: Open 1200 psi, Closed 150 psi

Lumen I.D.: 0.065 inches

Connectors: Male to Female





**Floswitch®**  
Sterile, single use


- Floswitch® acts as an ON/OFF device on infusion lines, stopping the flow in order to carry out line procedures in closed circuit.
- Designed to reduce risk of patient emboli and contamination of hospital staff to blood.
- The Floswitch® is designed to increase safety.
- Floswitch® can be connected to any catheter fitted with a Luer connector.
- Winged to facilitate attachment.

**Floswitch® closed**

External compression of the extension tubing by means of a steel ball, provides risk-free unambiguous opening and closing of the infusion line.



**Floswitch® open**



**ORDERING INFORMATION**

PRODUCT CODE	DESCRIPTION	BOX/CASE CONTENT (UNITS)
681033	Floswitch	25/250

**PRODUCT SPECIFICATIONS**

Latex free

CE marked






















Class IIa according to the European Medical Devices Directive 93/42/ EEC as amended by Directive 2007/47/EC

Sterilisation method: Ethylene Oxide

## Quality Management with Merit

Merit products are designed and manufactured to meet applicable international Quality Management System requirements in accordance with the European Medical Devices Directive 93/42/EEC (MDD) and FDA QSR where applicable. Hence, all Merit manufacturing sites hold a valid certification against ISO 9001:2000 and/or ISO 13485:2003. Likewise, all Merit manufacturing sites are subject of periodic external as well as internal audits to ensure compliance with aforementioned Quality Management System requirements.


## Symbol Guidance

	The CE mark, signifying compliance with European Medical Device Directive.	
	Catalogue number/ reference number/ re-order number	 Serial number
	Lot number/batch code/batch number	 Fragile
	Authorised representative in the European Community	 Keep dry
	Sterilised using Ethylene Oxide	 Do not re-use/single use/use only once
	Use by	 Keep away from sunlight
	Date of manufacture	 Do not use if package is damaged
	Manufacturer	 Green Dot - Recycle in Europe
	This way up	 Caution/attention/see instructions for use
	Does not contain natural rubber latex	 Contains or presence of natural rubber latex
	Contains or presence of phthalates: Bis(2-ethylhexyl) phthalate (DEHP)	 Temperature limitation
<b>Rx ONLY</b>	Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician	 Caution, consult accompanying documents


# Cannulae/Catheter Conversions

Gauge	Outer Ø mm	Inch	Charriere (French) 1Ch = 0.33mm
	0.5	0.019	2
24	0.6	0.024	
23	0.7	0.027	
22	0.8	0.031	
21	0.9	0.035	3
20	1.0	0.039	
19	1.1	0.043	
18	1.2	0.047	4
	1.3	0.051	
17	1.4	0.055	
	1.5	0.059	
16	1.6	0.063	5
	1.7	0.067	
15	1.8	0.071	
14	1.9	0.075	6
	2.0	0.079	
13	2.1	0.083	
	2.2	0.086	
13	2.3	0.091	7
	2.4	0.094	
	2.5	0.098	

The definition of gauge measurement comes from the AAMI (Association for the Advancement of Medical Instrumentation - USA). This states how many cannulae fit into a tube with an inner diameter of 1" (1 inch = 25.4mm).




The AAMI uses only even numbers (18, 20, 22, etc.) as measuring units. A similar measurement definition exists in Great Britain and this is known as the SWG (Standard Wire Gauge). The SWG uses consecutive numbers, ie 13 to 24 and is the more common method of measurement for catheters in Europe. The Charriere, commonly known as French, relates directly to the catheter size, i.e. 1 Ch = 0.33mm. In the following tables the measurements of the AAMI - SWG - French are shown in relation to the outer Ø in mm.



Merit's multi lumen catheter (Careflow™) hub colour coding differentiates the lumens and not the catheter gauge, i.e.


white	○	proximal
brown	●	tip
blue	●	middle





# UNDERSTAND. INNOVATE. DELIVER.™

Merit Medical Systems is a worldwide leader in the development, manufacturing, and distribution of quality medical devices used in interventional and diagnostic procedures. With a global workforce and cutting-edge manufacturing facilities equipped with the most advanced technology, Merit Medical consistently delivers quality products that are backed by excellent service and professional expertise.



Understand. Innovate. Deliver.™

**Merit Medical Systems, Inc.**  
 1600 West Merit Parkway  
 South Jordan, Utah 84095  
 1.801.253.1600  
 1.800.35.MERIT

**Merit Medical Europe, Middle East, & Africa (EMEA)**  
 Amerikalaan 42, 6199 AE Maastricht-Airport  
 The Netherlands  
 +31 43 358 82 22

**Merit Medical Ireland Ltd.**  
 Parkmore Business Park West  
 Galway, Ireland  
 +353 (0) 91 703 733

[merit.com](http://merit.com)

**Australia**  
 +61 (0) 3 8373 4943

**Austria**  
 0800 295 374

**Belgium**  
 0800 72 906 (Dutch)  
 0800 73 172 (Français)

**Brazil**  
 +55 11 4561 - 0062

**Canada**  
 1.800.364.4370

**China**  
 +86-10 8561 0788

**Denmark**  
 80 88 00 24

**Finland**  
 0800 770 586

**France**  
 0800 91 60 30

**Germany**  
 0800 182 0871

**Hong Kong**  
 +852-2207 6888

**India**  
 +91-80-41223376

**Ireland (Republic)**  
 1800 553 163

**Italy**  
 800 897 005

**Luxembourg**  
 8002 25 22

**Netherlands**  
 0800 022 81 84

**Norway**  
 800 11629

**Portugal**  
 308 801 034

**Russia**  
 +7 495 221 89 02

**Singapore**  
 +65 67536311

**South Korea**  
 +8210 2205 0707

**Spain**  
 +34 911238406

**Sweden**  
 020 792 445

**Switzerland**  
 (Deutsch)  
 +41 22 518 02 30

(Français)  
 +41 22 518 02 52

(Italiano)  
 +41 22 518 00 35

**UK**  
 0800 973 115