

# PHILIPS

Image guided therapy

Zenition 50 Series

Mobile C-arms



# Philips Zenition 50 Series

## Specifications

Point and shoot simplicity and reliable uptime make the Zenition 50 Series, the ideal image-intensifier-based systems for intensive use in interventions and surgeries. The Zenition harmonized range of systems opens new opportunities for you to reduce operational costs, simplify use and streamline fleet management. All systems offer the same intuitive control and handling for every user. Every Zenition system shares a standard Windows® platform and is ready to embrace new clinical applications and service and support technologies as they evolve. The Technology Maximizer program helps you manage upgrade costs for years to come.







# Contents

1	Image intensifier	4
2	Geometry	5
3	Imaging	6
4	X-ray generation	10
5	Workflow	14
6	Clinical extensions	15
7	System overview	16
8	Dimensions	19
9	Options	20
10	Service	22

# 1 Image Intensifier



Zenition 50 Series with 9" image intensifier



Zenition 50 Series with 12" image intensifier

Specifications	9" image intensifier	12" image intensifier
<b>Image intensifier</b>	The Zenition 50 Series comes with a 9" image intensifier and can go wherever you need it - surgery, intensive care and the emergency room.	The Zenition 50 Series comes with a 12" image intensifier and can go wherever you need it - surgery, intensive care and the emergency room.
<b>Image intensifier type</b>	Triple-mode HRC	Triple-mode
<b>Nominal II formats</b>	23, 17, and 13 cm (9", 7", and 5")	31, 23, and 17 cm (12", 9", and 7")
<b>Grid type</b>	Focused Grid, Circular, Carbon fibre cover; 60 lines/cm Ratio = 1:10 SID = 100 cm	Focused Grid, Circular, Carbon fibre cover; 60 lines/cm Ratio = 1:10 SID = 100 cm
<b>TV camera type</b>	CCD; high-resolution 1 k2	CCD; high-resolution 1 k2
<b>Image rotation</b>	Digital, live and on LIH	Digital, live and on LIH
<b>Image reversal</b>	Yes	Yes
<b>Mirror up/down</b>	Digital, live and on LIH	Digital, live and on LIH
<b>Mirror left/right</b>	Digital, live and on LIH	Digital, live and on LIH
<b>Automatic anatomical measuring field</b>	Yes with 'BodySmart' and 'MetalSmart'	Yes with 'BodySmart' and 'MetalSmart'
<b>Add central resolution in all magnifications</b>	Typical central visual resolution (lp/mm) 23 cm: 4.2 17 cm: 5.0 13 cm: 5.7	Typical central visual resolution (lp/mm) 31 cm: 4.0 23 cm: 4.8 13 cm: 6.0



# 2 Geometry

## 2.1 C-arm stand

<b>Angulation</b>	115° rotation (+90°/-25°)
<b>Extended Angulation (option)</b>	135° rotation (+90°/-45°) for increased projection flexibility
<b>Motorized height movement Speed</b>	490 mm/19.3" typical: up ~1.8 cm/s, down ~2.4 cm/s
<b>Longitudinal movement</b>	200 mm/7.9"
<b>Panning movement (swivel)</b>	± 10°
<b>Rotation</b>	± 180°, with safety stop at ± 135°
<b>Lowest lateral position</b>	1034 mm (C-arc under the table) 740 mm (C-arc over the table)
<b>Source to image distance (SID)</b>	98.3 cm
<b>Free space in C-arm</b>	76.6 cm
<b>C-arm depth</b>	610 mm
<b>Parallel movement</b>	Dedicated parallel movement with rear-wheel steering, for easy positioning along operating table
<b>C-arm stand length</b>	2060 mm (81.1")
<b>Weight</b>	310 kg for 9" / 312 kg for 12"
<b>C-arm stand width</b>	815 mm
<b>C-arm stand height</b>	1680 mm for 9" / 1770 mm for 12"
<b>Brakes on all movements</b>	Yes, manual
<b>Steering</b>	Rear wheels
<b>Cable deflectors</b>	Yes

## 2.2 Mobile View Station specifications

<b>Depth</b>	701 mm
<b>Width</b>	702 mm (monitors folded) 943 mm (monitors unfolded)
<b>Height</b>	1850 mm
<b>Weight (including options)</b>	140 kg
<b>Monitor rotation</b>	180°
<b>Monitor height movement</b>	230 mm



# 3 Imaging

---

## Specifications

---

### SmartVision

- Unique BodySmart software allows free positioning of the anatomy, even at the edge of the image detector. It automatically detects anatomy and adjusts parameters to produce high quality images
  - Contrast and brightness can be adjusted automatically in real time or manually
  - Our fully automatic MetalSmart feature excludes metal artifacts caused by metal implants to provide higher image quality and efficient dose control for orthopedic and other procedures, compared to systems without metal exclusion
  - Philips premium imaging technologies correct for patient or accidental table motion, automatically and in real-time on live images; reduce noise and artifacts, also on moving structures and objects; and enhance images and sharpen edges
- 





---

## Specifications

---

### DoseWise

- Imaging modes:
    - Fluoro modes ranging from Low X-ray dose to High Level, enabling dose savings when desired or enhanced image quality when necessary
    - Three different pulse rates for fluoro modes; the lower pulse rate can help to manage X-ray dose
  - Collimation:
    - Graphical shutter and iris and image orientation on Last Image Hold image on C-arm stand tablet-like UI without applying radiation
    - Real lead asymmetrical shutters
    - Independent shutter positioning
    - Automatic Electronic Blanking following the lead shutters and iris to enhance image quality
  - Automatic Shutter Positioning (ASP) sets shutters to the anatomy of interest for excellent image quality at the touch of a button
  - The II side laser (optional) and the tube side laser (optional) lets staff position the C-arm without using radiation
  - Several features contribute to increased dose awareness, including dose reporting, dose display, and an alert when exceeding a pre-defined procedure dose level
- 

### Acquisition settings

Preset acquisition settings apply dedicated fluoroscopy settings to obtain superb image quality for the anatomy of interest without applying more X-ray dose than necessary. Within each program there are different X-ray modes available (depending on anatomy of interest):

- Low Dose Fluoroscopy
- Normal Dose Fluoroscopy
- Medium Dose Fluoroscopy
- High Dose Fluoroscopy
- Exposure run to produce high quality images
- Single shot exposure, for extra-sharp, single snapshot images
- Hiding choices on the UI
- Auto Contrast Brightness (ACB) on/off settings
- Subtraction
- Trace

Blur reduction and noise reduction buttons to further adjust the level of temporal noise reduction level to the amount of movement in the region of interest

---

Removable grid



---

## Specifications

---

### Real time processing functions

---

<b>Feature</b>	<ul style="list-style-type: none"><li>• Feed-forward gain control</li><li>• White compression</li><li>• Adaptive temporal recursive noise reduction</li><li>• Spatial noise reduction</li><li>• Adaptive multiresolution brightness / contrast / edge enhancement</li><li>• Blanking</li><li>• Video invert</li><li>• Digital image rotation</li><li>• Mirroring</li><li>• Flipping</li><li>• View trace</li><li>• Subtract</li><li>• Trace</li><li>• Manual/auto contrast/brightness</li></ul>
----------------	---

---

### Post processing functions

---

<b>Feature</b>	<ul style="list-style-type: none"><li>• 360° digital rotation mirror left/right and up/down without radiation</li><li>• Contrast and Brightness and edging</li><li>• Annotation (for a single image or all images in an examination)</li><li>• Video invert (negative)</li><li>• Zoom and roam (factor 2x real-time magnification, freely movable to any section of an image)</li><li>• Measurement (to precisely quantify lengths and angles in images)</li><li>• Electronic shutters (to block out over exposed areas outside the region of interest)</li><li>• View Trace</li><li>• Landmarking (mixing anatomy with vessel tree)</li></ul>
----------------	--

---

### Mobile view station monitors

---

<b>Resolution</b>	1280 x 1024 pixels
<b>Maximum light output</b>	650 cd/m <sup>2</sup> (High Brightness, optional)
<b>Maximum light output</b>	330 cd/m <sup>2</sup> (Standard Brightness)
<b>Contrast ratio</b>	>500:1 (standard brightness)
<b>Viewing angle</b>	170° in horizontal and vertical direction
<b>Tablet-like UI</b>	Offers easy access to post-processing of acquired images, patient demographics as well as patient administration and data export (live monitor) (optional)
<b>Monitors</b>	Two 19" Standard Brightness (high brightness optional) Color LCD monitors for diagnostic image quality display. Monitor LUT: <ul style="list-style-type: none"><li>• DICOM GSDF compliant</li></ul>

---



# 4 X-ray generation

<b>Specifications</b>	
<b>X-ray generator</b>	The Zenition 50 Series uses a monoblock architecture with the high tension transformer in the X-ray tank. With the monoblock there is no need to transmit pulses over high voltage cables, which can result in a ramping up and ramping down effect, due to the electrical impedance of the cables. Because the monoblock generator operates at high frequencies (78.125 kHz and 15.625 kHz), it produces sharp pulses, which results in fewer motion artifacts in the image. This also allows less soft radiation to be used and produces less heat.
<b>X-ray tube</b>	Zenition 50 series have a rotating anode and high power generator with excellent heat management to perform the most demanding interventional procedures.
<b>Tube type</b>	Rotating anode X-ray tube
<b>Nominal x-ray tube voltage</b>	120 kV
<b>Nominal focal spot values</b>	0.3 and 0.6 mm
<b>Maximum anode heat content</b>	315 kHU = 225kJ
<b>Maximum anode heat dissipation</b>	54 kJ/min = 75.6 kHU/min = 900 W
<b>Max. uninterrupted fluoro time</b>	10 minutes legal block
<b>Maximum continuous anode heat dissipation</b>	550 W
<b>Anode rotating speed</b>	3000 rpm
<b>Anode target angle</b>	10°
<b>Anode material</b>	RT-TZM (Rhenium-Tungsten-Titanium-Zirconium-Molybdenum)
<b>Maximum housing heat content</b>	1350 kJ = 1890 kHU
<b>Maximum housing heat dissipation</b>	192 W = 11.5kJ/min = 16.1kHU/min
<b>Maximum continuous housing heat dissipation</b>	120 W = 7.2kJ/min = 10.1 kHU/min
<b>Cooling method</b>	Active oil-circulation cooling
<b>Inherent filter</b>	0.75 mm Al eq at 75 kV
<b>Integrated beam filter</b>	1.0 mm Al. + 0.1 mm Cu
<b>Total beam filtration</b>	4.73 mm Al eq.
<b>Maximum Generator Output</b>	15 kW

<b>Specifications</b>	
<b>Operating values with Exposure Runs</b>	
kV range	40 - 120 kV
mA peak range	2.5 - 60 mA
Pulse width	7.407 - 31.11 ms
Pulse rate	4, 7.5, 15 pps (optional 30 pps)
<b>Operating values with Continuous Fluoroscopy</b>	
kV range	40-120 kV
mA range	0.10 - 20 mA
<b>Operating values with Pulsed Fluoroscopy</b>	
kV range	40-120 kV
mA peak range	0.5 - 60 mA
Pulse width	8 - 22.22 ms
Pulse rate	6.5, 12.5 pps
<b>Operating values with Single shot exposure (snapshot)</b>	
kV range	40-120 kV
mA peak range	2.5 - 100 mA (normal) 5.2 - 125 mA (high power)
Time range	100, 167, 300 ms (normal) 100 ms (high power) (Optional)
<b>X-ray collimation</b>	The Zenition 50 Series makes collimation easy. Its full lead (assymmetric) shutters can be rotated and moved independently, and the unique Philips Automatic Shutter Positioning (ASP) feature automatically positions shutters for high image quality at the touch of a button. You can position shutters or adjust the iris on the last X-ray image (Last Image Hold on Stand UI), eliminating additional X-ray dose during collimation.



## Specifications

<b>Shutters</b>	Two independent lead shutters with steel wedge: shutters can be coupled for rotation and translation, or moved individually for asymmetric collimation
<b>Automatic shutter positioning</b>	Automatic shutter placement based on image content
<b>Shutter material</b>	3 mm Pb
<b>Rotation of shutter</b>	360°
<b>Wedge material</b>	0.2 to 2.5 mm stainless steel
<b>Adjustment of shutter and iris</b>	Stepless
<b>Iris material</b>	Lead with 5% antimony: Pb (Sb5%)
<b>Iris diameter (at detector)</b>	12" <50 mm - 295 mm 9" <50 mm - 240 mm
<b>Position indication</b>	On screen and also on last image hold without radiation on Stand UI



Zenith 50 Series with 12 inch Image Intensifier

# 5 Workflow

---

## Specifications

---

### Unify workflow\*

Unify workflow brings intuitive control and handling to your system to reduce training and enhance teamwork.

• Tablet-like UI

Physicians and operators experience a whole new level of simplicity with our tablet-like user interface on the C-arm and intelligent workflow. Now you can just touch the screen with a finger to drag the shutters and iris into position on Last Image Hold. At each step you only see the features you need, making it easy to find the right selections.

• ClearGuide and color coding

Our unique ClearGuide in combination with color coding on the C-arm speeds up positioning. ClearGuide provides a uniform reference for the operator and physician to use during positioning. A set of numbers (3, 6, 9, 12) on the Detector corresponds to the same numbers displayed on the clinical image. The physician can say 'Rotate orange towards 3' and the operator knows exactly what to do. The numbers always match up, even if the image is rotated, flipped, or mirrored.

---

### Mobile View Station

The compact Mobile View Station fits perfectly in the surgical workflow. Its intelligent design provides the user with easy system setup, enhanced viewing capabilities and easy transportation. Its unique design also makes it easy to clean. All system controls are at your fingertips on the live monitor of the Mobile View Station. With this tablet-like UI, you can intuitively set up an exam, post-process images, or export a case to PACS. The operator can easily see the clinical image on this 15" stand UI screen without obstructing the view of the physician. Because we use advanced infra-red technology, the touchscreen monitor delivers the same image quality as non-touchscreen monitors.

---

### Connectivity

- The Wireless Data Transfer option allows users to connect to the RIS/HIS to send and retrieve images or other relevant data wirelessly and reduce the amount of cable clutter in the OR
  - Optional Digital Video out to display live and reference images on additional monitors (e.g. ceiling mounted) without a loss of resolution
  - Optional Video in allows you to conveniently display external video signals like endoscopy or ultrasound on the right C-arm monitor, so all the information needed is in one view
  - Optional USB storage provides a convenient way to store and images for use in reports or presentations
- 

\*Position memory not available for Zenition 50 Sereis



ClearGuide



## Specifications

<b>DICOM</b>	<p>DICOM is seamlessly integrated into the system for digital image to DICOM translation. A highly intuitive user interface simplifies use.</p> <ul style="list-style-type: none"> <li>· DICOM print (optional)</li> <li>· DICOM store (optional)</li> <li>· Modality Worklist Management (MWL) (optional)</li> <li>· Modality Performed Procedure Step (MPPS) (optional)</li> <li>· Storage Commit (optional)</li> <li>· DICOM Storage to DVD or USB memory (optional)</li> <li>· DICOM query/retrieve (optional)</li> <li>· DICOM Radiation Dose Structured Reports</li> </ul> <p>DICOM image formats:</p> <ul style="list-style-type: none"> <li>· Secondary Capture (SC) with/without text</li> <li>· Angiography (XA - multi frame)</li> <li>· Patient dose report</li> </ul>
<b>Integrated Healthcare Enterprise (IHE)</b>	<p>The Zenition 50 Series is compliant with the IHE Scheduled Workflow Integration Profile as an Acquisition Modality Actor.</p> <p>For an optimal viewing angle, the LCD monitors can be rotated 180° and adjusted in height (230 mm/9" )</p>
<b>Digital video out (optional)</b>	2 DVI connectors live and reference monitor
<b>Video in (optional)</b>	S-Video, DVI (digital and analog), SDI
<b>USB storage (optional)</b>	PNG, MP4, BMP
<b>IP addressing (optional)</b>	Static IP, DHCP
<b>Wireless standards supported</b>	IEEE 802.11a / b / g / n / ac (2.4 GHz and 5 GHz band) FIPS 140-2:compliant
<b>Number of antennas</b>	2 (embedded within the system, not visible)
<b>User configurable SSID support</b>	Up to 16 SSIDs, each with a unique MAC address and configurable SSID Broadcast
<b>Authentication protocols</b>	PSK, IEEE 802.1x EAP-TLS and PEAP
<b>Security</b>	AES, TKIP and WEP encryption FIPS 140-2 compliant whitelisting to prevent malware
<b>External room X-ray indication</b>	Yes (optional)



Color coding

<b>PC Hardware details</b>	
# of USB ports	2 USB 2.0 and 1 High speed USB 3.0 port
Storage	up to 140,000 images
DICOM STORE (DVD/USB/ PACS) and Retrieve (USB/DVD/PACS) both are optional	Yes
Embedded MMV	Image Viewer (optional)
Service Tools (PSC, Remote, LOTS)	On system service tool (Philips Support Connect) Remote Service Remote Assistance (Look Over The Shoulder, LOTS)
Boot up time	< 77 seconds
Operating System	Windows® 7 Embedded
• Processor speed	Intel Core™ i7-4790S (4 GHz)
• RAM	8 GB: 2x DDR3 1600 MHz 4 GB SO-DIMM
• Storage type	2x 1TB HDD
Image processing bits	II subsystem: 12 bits, system image processing: 14 bits
Storage capacity in GB	2x 1TB HDD, of which ± 800 GB or 140,000 images for image storage
Storage image matrix	1 k x 1 k
Storage image bits	14 bit image data + 1 bit measuring field




# 6 Clinical extensions

<b>Specifications</b>	
<b>Outlining</b>	The outlining tool allows users to draw an outline (in 5 different colors, 1 color at a time) digitally on an image on the touchscreen monitor using a stylus pen, for example, to mark a bifurcation or side branches on live fluoroscopy images. Drawings can be switched/toggled on and off. The undo button corrects mistakes and the delete function gets rid of the drawing with the touch of a button.
<b>Pain management extension</b>	The Pain management extension offers digital subtraction functionality to enable enhanced visualization of contrast medium injections. Pain interventionalists can work with confidence to avoid potentially dangerous intravascular injections when using digital subtraction to enhance visualization of contrast distribution.
<b>Cardiac extension</b>	This combination of dedicated cardiac exam types, high pulse rates and image memory of 140,000 images is the ideal package for cardiac interventions. The cardiac extension includes dedicated anatomical programmed fluoroscopy parameters for electrophysiology procedures, pacemaker placements, and cardiac surgery. With motion-stopping pulses of up to 30 pulses per second and 60 mA, the Zenition 50 Series captures sharp images of fast moving anatomy in the region of interest. Standard memory (140,000 images) provides capacity to record long cases at high frame rates.
<b>Vascular extension</b>	The vascular extension offers you the full support for vascular cases, providing an extensive range of vascular imaging tools. Most vascular functions can be controlled by handheld remote or at the user interface on the Mobile View Station. Image memory of 140,000 images available.
<b>Vascular Processing</b>	<ul style="list-style-type: none"> <li>• Subtracted fluoroscopy mode displays images in subtracted mode</li> <li>• Trace mode shows in real time the maximum opacification of the vessels (Iodine + CO<sub>2</sub>)</li> <li>• Roadmap images support catheter guidance</li> <li>• Remask lets you reselect the most suitable image in your run as a mask image for contrast runs</li> <li>• Smart Mask manages dose and contrast medium usage by re-using previously acquired images for roadmapping</li> <li>• Landmarking provides a non-subtracted background image for anatomical reference</li> <li>• Real time pixel shift compensates for movement artifacts</li> <li>• Subtraction on/off simplifies the orientation for subtracted images during roadmap procedures (controlled by remote control or User Interface on the Mobile View Station)</li> <li>• View Trace creates a trace image in post-processing (Iodine + CO<sub>2</sub>)</li> <li>• CO<sub>2</sub> mode for subtraction, trace white and roadmap with Smart Mask and view trace</li> </ul>
<b>Cardiovascular extension</b>	The cardiovascular extension provides high quality fluoroscopy, subtraction runs, and roadmap guidance to support challenging cardiac and vascular procedures. Vascular processing features include subtracted fluoroscopy mode, trace mode to show the maximum opacification of vessels in real time, and real-time pixel shift to compensate for movement artifacts. Extra support for cardiac examinations is provided with dedicated anatomically programmed fluoroscopy parameters for electrophysiology procedures, pacemaker placements, and cardiac surgery. With motion-stopping pulses of up to 30 pulses per second and 60 mA, the Zenition 50 Series captures high quality images of fast moving anatomy in the region of interest.



# 7 System overview



- 1 Tablet-like UI (Unify Workflow)
- 2 Application specific protocols and customizable presets
- 3 Adjust shutter and image orientation while on LIH
- 4 Color coding (Unify Workflow)
- 5 Fully Counterbalanced C-arc
- 6 2 image intensifier sizes (9" and 12")
- 7 ClearGuide (Unify Workflow)
- 8 Removable X-ray grid
- 9 Unique asymmetric shutters
- 10 A smart new way of filtering
- 11 BodySmart
- 12 MetalSmart
- 13 15kW system
- 14 Compact X-ray tank
- 15 Monoblock design
- 16 Technology Maximizer
- 17 Remote control
- 18 Point and Shoot design



- 19 Outlining
- 20 High Speed wireless data transfer / Wireless data handling
- 21 Image Viewer/MMV
- 22 Lightweight MVS
- 23 Advanced Image Processing algorithms/radiation management features
- 24 80 seconds boot-up
- 25 Remote Support + Philips Support Connect tool
- 26 Windows® Platform
- 27 Advances Connectivity and interoperability tool
- 28 External Video-In

# 8 Dimensions


## Specifications

### C-arm stand


- 9" IITV: height x width x depth: 45.7 cm x 28.8 cm x 36.3 cm
- 12" IITV: height x width x depth: 54.7 cm x 38.5 cm x 41.2 cm
- Tube Tank: height x width x depth: 39.4 cm x 18.7 cm x 46.7 cm
- SID (focus to II surface distance): 9" IITV: 982 mm, 12" IITV: 983 mm
- Focus to end of collimator assembly distance: 15.3 cm (source to diaphragm distance); 20.3 cm (source to IEC 20 cm spacer cover); 30.3 cm (source to HHS 30 cm spacer cover)

### Mobile View Station

- Monitors: Height x width x depth: 37.3 cm x 42.4 cm x 10 cm
- Find drawings on next page







### Mobile View Station







### Tube tank, new

### Tube tank, old



\* FOR EXTENDED ROTATION : 45°  
 \* FOR STANDARD ROTATION : 25°

Zenition 50 Series with 12" Image Intensifier



\* FOR EXTENDED ROTATION : 45°  
 \* FOR STANDARD ROTATION : 25°

Zenition 50 Series with 9 inch Image Intensifier



# 9 Options

## Specifications

<b>Tank laser aiming device</b>	Optional laser projects a crosshair from the X-ray tank towards the Image Intensifier, indicating the center of the X-ray beam and enabling alignment of the C-arm without X-ray.
<b>II laser aiming device</b>	Positioning device for use at the image intensifier side
<b>Video paper printer</b>	Thermal printer to print video images from the live (left) monitor on to paper during or after examinations. Print 1, 2, 4, or 6 images on one page in landscape or portrait format. Sony UP-97xAD printer.
<b>Video paper/transparency printer</b>	Thermal printer to print video images from the live (left) monitor on to paper or transparencies during or after examinations. Print 1, 2, 4, or 6 images on one page in landscape or portrait format. Sony UP-99xAD printer.
<b>DICOM and IHE</b>	Zenition 50 can be equipped with the Philips Integrated DICOM solution, which transfers images from the Zenition 50 onto the hospital network in a Secondary Capture DICOM SC or a DICOM XA format. The Basic DICOM package supports DICOM Print and DICOM Store. The advanced DICOM/IHE package (optional) supports: <ul style="list-style-type: none"><li>• Modality Worklist Management</li><li>• Modality Performed Procedure Step</li><li>• Storage Commit</li><li>• Full compliance to the IHE Scheduled Workflow integration profile as an Acquisition Modality Actor</li><li>• Supports DICOM Structured Dose Reports (standard feature)</li><li>• Query/Retrieve (Image Viewer option)</li></ul>



Handheld remote control



Wireless footswitch

---

## Specifications

---

<b>Image Viewer</b>	Offers an intuitive multi-purpose platform for retrieving and handling DICOM images from different modalities. It lets you compare pre-operative images side-by-side with the live fluoroscopy images. 500 Gbyte hard disk. <ul style="list-style-type: none"><li>• MIP/MPR - maximum intensity projection singles out high intensity areas for optimized 2D projection of a 3D volume</li></ul>
<b>Wireless data transfer</b>	The Wireless Data Transfer option allows users to connect to the RIS/HIS to send and retrieve images or other relevant data wirelessly and reduce the amount of cable clutter in the OR.
<b>Handheld remote control</b>	The remote control unit is a handheld infrared keypad used to control the main image handling functions. For sterile operation, it can be used in a transparent sterile plastic cover. The functions include: <ul style="list-style-type: none"><li>• Run loop</li><li>• Overview run/exam</li><li>• Retrieve previous image/run</li><li>• Retrieve next image/run</li><li>• Park image on Reference monitor (Smart mask)</li><li>• Retrieve image from Reference monitor (Smart mask)</li><li>• Protect image/release image</li><li>• Mode selection</li><li>• Detector-format selection</li><li>• Subtraction on/off</li><li>• Image grab</li></ul>
<b>Footswitch</b>	wired footswitch cable length: 3.5 m and Wireless

---



Wired footswitch

# 10 Service

<b>Specifications</b>	
<b>RightFit Customer Service Agreement</b>	Different RightFit Customer Service Agreements are available that allows you to leverage the capabilities of your in-house service teams to maintain and service your Zenition 50 series to reduce delays and improve uptime.
<b>Technology Maximizer</b>	A program that runs in tandem with your RightFit Customer Service Agreement. <sup>3</sup> When you opt into the program, you receive the latest available software and hardware technology releases for a predictable fee.
<b>Planning</b>	<ul style="list-style-type: none"> <li>• System installation project management</li> <li>• Room design services</li> </ul>
<b>Start-up</b>	<ul style="list-style-type: none"> <li>• Clinical application education</li> <li>• In-house service training</li> <li>• Online learning center</li> </ul>
<b>Peak usage</b>	<ul style="list-style-type: none"> <li>• Comprehensive, partner and flexible service agreements</li> <li>• Remote services</li> <li>• Service information portal</li> </ul>
<b>Renewal</b>	<ul style="list-style-type: none"> <li>• Installed base programs</li> <li>• System relocation services</li> <li>• Refurbished systems</li> </ul>
<b>Remote proactive support</b>	One valuable feature in our Service Agreements is remote support. It helps you get the most from your imaging system and maintain its peak performance every day. Philips Remote Services is an advanced, secure network that links your Zenition 50 Series to our global remote services customer care centers. Configuration, customization, log file analysis and other services that formerly required on-site visits are now available by connecting to our remote experts. Remote system log analysis allows our experts to detect anomalies and plan an on-site visit to increase efficiency. Remote service helps to maintain equipment performance – protecting you against lengthy downtime and unexpected costs. If a deteriorating situation is detected by our remote service, corrective action can be carried out quickly and decisively, often with no interruption to your busy schedule. Our application specialist can also use the remote viewing feature to resolve application-related issues more quickly, without an on-site visit. A global platform for system communication certifies that all service data is handled via best-in-class encrypted transmission technology.

<sup>3</sup> Eligible RightFit Service Agreements are available with Technology Maximizer.





PHILIPS

Arne







© 2019 Koninklijke Philips N.V. All rights reserved. Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

4522 991 40571 \* FEB 2019

**How to reach us**  
Please visit [www.philips.com](http://www.philips.com)  
[healthcare@philips.com](mailto:healthcare@philips.com)