

Bioline™ TSUTSUGAMUSHI

SCRUB TYPHUS ANTIBODY TEST

The Bioline™ Tsutsugamushi test is a solid phase, immunochromatographic assay for the rapid, qualitative detection of IgG, IgM or IgA antibodies to *Orientia tsutsugamushi* in human serum, plasma or whole blood.

- Disease: Scrub Typhus
- Specimen : Serum, Plasma (10 µl) / Whole blood (20 µl)
- Test result : 10-15 minutes
- Shelf life and storage temperature: 18 months from the date of manufacturing at 1-30 °C
- Performance: Sensitivity 99 %, Specificity 96 %




MATERIALS PROVIDED

- Test device
- Assay diluent

SIMPLE PROCEDURE


1 Add Specimen

Dispense 10 µl of serum, plasma, or 20 µl of whole blood into the specimen well "S".



2 Add Assay Diluent


Dispense 3-4 drops of the assay diluent.



Wait 10-15 mins.



RESULTS INTERPRETATION



The presence of any test line, no matter how faint, the result is considered positive.

ORDERING INFORMATION

PRODUCT	CAT. NO.	TYPE	SPECIMEN	PACK SIZE
Tsutsugamushi	18FK10	Device	Serum/Plasma/Whole blood	30T/Kit
Tsutsugamushi	18FK11	Multi-Device	Serum/Plasma/Whole blood	10Tx10/Kit