

Performa Diagnostic Cardiology Catheters

Introducing Merit Medical's Performa line of Diagnostic Cardiology Catheters, designed with all the torque you need, improved shaft strength for better pushability and large inner lumens for increased flow rates.

INNER-LUMEN SIZES

- 4F 0.042" (1.07 mm)
- 5F 0.046" (1.17 mm)
- 5F 0.052" (1.32 mm) high flow
- 6F 0.054" (1.37 mm)
- 6F 0.059" (1.49 mm) high flow
- Yields higher flow rates at lower pressures
- Superior guide wire trackability

WINGED POLYCARBONATE HUB

- Enhanced handling and control

STRAIN RELIEF


- Increased strength at catheter/hub interface
- Color coded for easy identification

TIP

- Soft tip for smooth atraumatic positioning
- Excellent tip memory



WIRE BRAID DESIGN

- Greater torque for maximum control
- Improved longitudinal stability
- High degree of maneuverability
- Provides superior pushability



PRESSURE


- High pressure limits at given flow rates



Performa Diagnostic Cardiology Catheters

Catalog Number	French Size	Internal Diameter	Length	Length of Tip	Description	
JUDKINS LEFT						
7700-10	4	0.042"/1.07mm	100cm	2.3cm	JL3.5	bumper tip
7700-20	4	0.042"/1.07mm	100cm	2.3cm	JL4.0	bumper tip
7700-30	4	0.042"/1.07mm	100cm	2.3cm	JL5.0	bumper tip
7700-40	4	0.042"/1.07mm	100cm	2.3cm	JL6.0	bumper tip
7521-13	5	0.046"/1.17mm	100cm	2.5cm	JL3.5	
7521-21	5	0.046"/1.17mm	100cm	2.5cm	JL4.0	
7521-F1	5	0.046"/1.17mm	100cm	2.5cm	JL4.5	
7521-33	5	0.046"/1.17mm	100cm	2.5cm	JL5.0	
7521-43	5	0.046"/1.17mm	100cm	2.5cm	JL6.0	
7521-53	5	0.046"/1.17mm	125cm	2.5cm	JL3.5	
7521-63	5	0.046"/1.17mm	125cm	2.5cm	JL4.0	
7521-73	5	0.046"/1.17mm	125cm	2.5cm	L5.0	
7521-83	5	0.046"/1.17mm	125cm	2.5cm	JL6.0	
7601-21	5	0.052"/1.32mm	100cm	2.5cm	JL4.0-high flow	non-braided
7601-33	5	0.052"/1.32mm	100cm	2.5cm	JL5.0-high flow	non-braided
7601-43	5	0.052"/1.32mm	100cm	2.5cm	JL6.0-high flow	non-braided
7501-13	6	0.054"/1.37mm	100cm	1.9cm	JL3.5	
7501-21	6	0.054"/1.37mm	100cm	1.9cm	JL4.0	
7501-H1	6	0.054"/1.37mm	100cm	1.9cm	JL4.5	
7501-33	6	0.054"/1.37mm	100cm	1.9cm	JL5.0	
7501-43	6	0.054"/1.37mm	100cm	1.9cm	JL6.0	
7501-53	6	0.054"/1.37mm	125cm	1.9cm	JL3.5	
7501-63	6	0.054"/1.37mm	125cm	1.9cm	JL4.0	
7501-73	6	0.054"/1.37mm	125cm	1.9cm	JL5.0	
7501-83	6	0.054"/1.37mm	125cm	1.9cm	JL6.0	
7574-13	6	0.059"/1.49mm	100cm	2.5cm	JL3.5-high flow	non-braided
7574-21	6	0.059"/1.49mm	100cm	2.5cm	JL4.0-high flow	non-braided
7574-33	6	0.059"/1.49mm	100cm	2.5cm	JL5.0-high flow	non-braided
7574-43	6	0.059"/1.49mm	100cm	2.5cm	JL6.0-high flow	non-braided
JUDKINS LEFT SHORT TIP						
7501-A3	6	0.054"/1.37mm	100cm	1.9cm	JL3.5 ST	
7501-B1	6	0.054"/1.37mm	100cm	1.9cm	JL4.0 ST	
7501-C3	6	0.054"/1.37mm	100cm	1.9cm	JL5.0 ST	
7501-D3	6	0.054"/1.37mm	100cm	1.9cm	JL6.0 ST	

Recommended guide wire is 0.038" (0.97 mm)




Performa Diagnostic Cardiology Catheters

Catalog Number	French Size	Internal Diameter	Length	Length of Tip	Description	
JUDKINS RIGHT						
7701-10	4	0.042"/1.07mm	100cm	2.5cm	JR3.5	bumper tip
7701-20	4	0.042"/1.07mm	100cm	2.5cm	JR4.0	bumper tip
7701-50	4	0.042"/1.07mm	100cm	2.5cm	JR4.5	bumper tip
7701-30	4	0.042"/1.07mm	100cm	2.5cm	JR5.0	bumper tip
7701-40	4	0.042"/1.07mm	100cm	2.5cm	JR6.0	bumper tip
7523-13	5	0.046"/1.17mm	100cm	2.5cm	JR3.5	
7523-21	5	0.046"/1.17mm	100cm	2.5cm	JR4.0	
7523-93	5	0.046"/1.17mm	100cm	2.5cm	JR4.5	
7523-33	5	0.046"/1.17mm	100cm	2.5cm	JR5.0	
7523-43	5	0.046"/1.17mm	100cm	2.5cm	JR6.0	
7523-53	5	0.046"/1.17mm	125cm	2.5cm	JR3.5	
7523-63	5	0.046"/1.17mm	125cm	2.5cm	JR4.0	
7523-73	5	0.046"/1.17mm	125cm	2.5cm	JR5.0	
7523-83	5	0.046"/1.17mm	125cm	2.5cm	JR6.0	
7503-13	6	0.054"/1.37mm	100cm	1.9cm	JR3.5	
7503-21	6	0.054"/1.37mm	100cm	1.9cm	JR4.0	
7503-93	6	0.054"/1.37mm	100cm	1.9cm	JR4.5	
7503-33	6	0.054"/1.37mm	100cm	1.9cm	JR5.0	
7503-43	6	0.054"/1.37mm	100cm	1.9cm	JR6.0	
7503-53	6	0.054"/1.37mm	125cm	1.9cm	JR3.5	
7503-63	6	0.054"/1.37mm	125cm	1.9cm	JR4.0	
7503-73	6	0.054"/1.37mm	125cm	1.9cm	JR5.0	
7503-83	6	0.054"/1.37mm	125cm	1.9cm	JR6.0	
JUDKINS RIGHT SHORT TIP						
7523-E3	5	0.046"/1.17mm	100cm	2.5cm	JR3.5 ST	
7523-F1	5	0.046"/1.17mm	100cm	2.5cm	JR4.0 ST	
7523-G3	5	0.046"/1.17mm	100cm	2.5cm	JR5.0 ST	
7523-H3	5	0.046"/1.17mm	100cm	2.5cm	JR6.0 ST	
7503-A3	6	0.054"/1.37mm	100cm	1.9cm	JR3.5 ST	
7503-B1	6	0.054"/1.37mm	100cm	1.9cm	JR4.0 ST	
7503-C3	6	0.054"/1.37mm	100cm	1.9cm	JR5.0 ST	
7503-D3	6	0.054"/1.37mm	100cm	1.9cm	JR6.0 ST	
JUDKINS RIGHT MODIFIED						
7502-13	6	0.054"/1.37mm	100cm	1.9cm	JR3.5	MOD
7502-23	6	0.054"/1.37mm	100cm	1.9cm	JR4.0	MOD
7502-33	6	0.054"/1.37mm	100cm	1.9cm	JR5.0	MOD
7502-43	6	0.054"/1.37mm	100cm	1.9cm	JR6.0	MOD



Performa Diagnostic Cardiology Catheters

Catalog Number	French Size	Internal Diameter	Length	Length of Tip	Description	
AMPLATZ LEFT						
7703-10	4	0.042"/1.07 mm	100 cm	5.0 cm	AL1	bumper tip
7703-20	4	0.042"/1.07 mm	100 cm	2.2 cm	AL2	bumper tip
7703-30	4	0.042"/1.07 mm	100 cm	2.3 cm	AL3	bumper tip
7703-40	4	0.042"/1.07 mm	100 cm	2.5 cm	AL4	bumper tip
7527-13	5	0.046"/1.17 mm	100 cm	5.0 cm	AL1	
7527-23	5	0.046"/1.17 mm	100 cm	5.0 cm	AL2	
7527-33	5	0.046"/1.17 mm	100 cm	2.0 cm	AL3	
7527-43	5	0.046"/1.17 mm	100 cm	2.5 cm	AL4	
7507-13	6	0.054"/1.37 mm	100 cm	5.0 cm	AL1	
7507-23	6	0.054"/1.37 mm	100 cm	5.0 cm	AL2	
7507-33	6	0.054"/1.37 mm	100 cm	2.4 cm	AL3	
7507-43	6	0.054"/1.37 mm	100 cm	2.5 cm	AL4	
AMPLATZ RIGHT						
7702-40	4	0.042"/1.07 mm	100 cm	3.3 cm	AR MOD	bumper tip
7702-10	4	0.042"/1.07 mm	100 cm	3.3 cm	AR1	bumper tip
7702-20	4	0.042"/1.07 mm	100 cm	5.0 cm	AR2	bumper tip
7702-30	4	0.042"/1.07 mm	100 cm	5.0 cm	AR3	bumper tip
7526-11	5	0.046"/1.17 mm	100 cm	3.3 cm	AR MOD	
7528-13	5	0.046"/1.17 mm	100 cm	3.3 cm	AR1	
7528-23	5	0.046"/1.17 mm	100 cm	5.0 cm	AR2	
7528-33	5	0.046"/1.17 mm	100 cm	5.0 cm	AR3	
7506-11	6	0.054"/1.37 mm	100 cm	3.3 cm	AR MOD	
7508-13	6	0.054"/1.37 mm	100 cm	3.3 cm	AR1	
7508-23	6	0.054"/1.37 mm	100 cm	3.5 cm	AR2	
7508-33	6	0.054"/1.37 mm	100 cm	4.0 cm	AR3	
LEFT CORONARY BYPASS						
7709-10	4	0.042"/1.07 mm	100 cm	2.5 cm	LCB	bumper tip
7534-13	5	0.046"/1.17 mm	100 cm	2.5 cm	LCB	
7510-13	6	0.054"/1.37 mm	100 cm	2.5 cm	LCB	
RIGHT CORONARY BYPASS						
7709-20	4	0.042"/1.07 mm	100 cm	2.5 cm	RCB	bumper tip
7534-23	5	0.046"/1.17 mm	100 cm	2.5 cm	RCB	
7510-23	6	0.054"/1.37 mm	100 cm	2.5 cm	RCB	
PROGRESSIVE RIGHT CORONARY (NON TORQUE)						
7707-10	4	0.042"/1.07 mm	100 cm	2.5 cm	PRC	bumper tip
7600-11	5	0.052"/1.32 mm	100 cm	2.5 cm	PRC	non-braided
7532-13	5	0.046"/1.17 mm	100 cm	2.5 cm	PRC	
7572-11	6	0.059"/1.49 mm	100 cm	2.5 cm	PRC	non-braided
7512-13	6	0.054"/1.37 mm	100 cm	2.5 cm	PRC	
7512-23	6	0.054"/1.37 mm	125 cm	2.5 cm	PRC	
INTERNAL MAMMARY						
7706-10	4	0.042"/1.07 mm	100 cm	3.3 cm	IMA	bumper tip
7531-30	5	0.046"/1.17 mm	65 cm	3.3 cm	IMA	
7531-13	5	0.046"/1.17 mm	100 cm	3.3 cm	IMA	
7511-13	6	0.054"/1.37 mm	100 cm	3.3 cm	IMA	



CATHETERS
& EXTRUSIONS


Recommended guide wire is 0.038" (0.97 mm)

Performa Diagnostic Cardiology Catheters


Catalog Number	French Size	Internal Diameter	Length	Length of Tip	Description	# Side Ports	Small Pigtail Curve	Angle Length
STRAIGHT CARDIAC PIGTAIL								
7704-70	4	.042"/1.07 mm	80 cm	6.5 cm	braided pig	5		
7704-10	4	.042"/1.07 mm	110 cm	6.5 cm	braided pig	5		
7776-11	5	.046"/1.17 mm	90 cm	5.0 cm	braided pig	5	yes	
7776-21	5	.046"/1.17 mm	110 cm	5.0 cm	braided pig	5	yes	
7776-B1	5	.052"/1.32 mm	110 cm	4.6 cm	non-braided pig	5	yes	
7602-11	5	.052"/1.32 mm	110 cm	7.5 cm	non-braided pig	5		
7529-11	5	.046"/1.17 mm	110 cm	7.5 cm	braided pig	5		
7529-43	5	.046"/1.17 mm	125 cm	7.5 cm	braided pig	5		
7569-11	6	.059"/1.49 mm	110 cm	8.0 cm	non-braided pig	5		
7509-11	6	.054"/1.37 mm	110 cm	7.5 cm	braided pig	5		
7777-B1	6	.059"/1.49 mm	110 cm	4.8 cm	non-braided pig	5	yes	
7777-21	6	.054"/1.49 mm	110 cm	4.6 cm	braided pig	5	yes	
7509-43	6	.054"/1.49 mm	125 cm	7.5 cm	braided pig	5		
7569-43	6	.059"/1.49 mm	125 cm	8.0 cm	non-braided pig	5		
145 DEGREE CARDIAC PIGTAIL								
7704-20	4	.042"/1.07mm	110 cm	6.5 cm	braided pig	5		5.5 cm
7776-C3	5	.052"/1.32 mm	110 cm	4.6 cm	non-braided pig	5	yes	5.5 cm
7776-33	5	.046"/1.17 mm	110 cm	5.0 cm	braided pig	5	yes	5.5 cm
7529-53	5	.046"/1.17 mm	125 cm	7.5 cm	braided pig	5		7.5 cm
7602-23	5	.052"/1.32 mm	110 cm	7.5 cm	non-braided pig	5		7.5 cm
7529-23	5	.046"/1.17 mm	110 cm	7.5 cm	braided pig	5		7.5 cm
7569-23	6	.059"/1.49 mm	110 cm	8.0 cm	non-braided pig	5		5.5 cm
7509-23	6	.054"/1.37 mm	110 cm	7.5 cm	braided pig	5		5.5 cm
7777-C3	6	.059"/1.49 mm	110 cm	7.5 cm	non-braided pig	5	yes	5.5 cm
7777-33	6	.054"/1.37 mm	110 cm	4.6 cm	braided pig	5	yes	5.5 cm
7509-53	6	.054"/1.37 mm	125 cm	7.5 cm	braided pig	5		5.5 cm
7569-53	6	.059"/1.49 mm	125 cm	8.0 cm	non-braided pig	5		5.5 cm
155 DEGREE CARDIAC PIGTAIL								
7704-30	4	.042"/1.07 mm	110 cm	6.5 cm	braided pig	5		5.5 cm
7602-33	5	.052"/1.32 mm	110 cm	7.5 cm	non-braided pig	5		7.5 cm
7529-33	5	.046"/1.17 mm	110 cm	7.5 cm	braided pig	5		7.5 cm
7529-63	5	.046"/1.17 mm	125 cm	7.5 cm	braided pig	5		7.5 cm
7776-D3	5	.052"/1.32 mm	110 cm	4.6 cm	ultra high flow	5	yes	5.5 cm
7776-43	5	.046"/1.17 mm	110 cm	5.0 cm	braided pig	5	yes	5.5 cm
7569-33	6	.059"/1.49 mm	110 cm	8.0 cm	non-braided pig	5		5.5 cm
7509-33	6	.054"/1.37 mm	110 cm	7.5 cm	braided pig	5		5.5 cm
7509-63	6	.054"/1.37 mm	125 cm	7.5 cm	braided pig	5		5.5 cm
7569-63	6	.059"/1.49 mm	125 cm	8.0 cm	non-braided pig	5		5.5 cm
7777-D3	6	.059"/1.49 mm	110 cm	4.8 cm	non-braided pig	5	yes	5.5 cm
7777-43	6	.054"/1.37 mm	110 cm	4.6 cm	braided pig	5	yes	5.5 cm
SONES								
7705-10	4	.042"/1.07 mm	80 cm	7.3 cm	SONES	4	bumper tip	
7705-20	4	.042"/1.07 mm	100 cm	7.3 cm	SONES	4	bumper tip	
7530-13	5	.046"/1.17 mm	80 cm	7.3 cm	SONES	4		
7530-23	5	.046"/1.17 mm	100 cm	7.3 cm	SONES	4		
7520-13	6	.054"/1.37 mm	80 cm	7.3 cm	SONES	4		
7520-23	6	.054"/1.37 mm	100 cm	7.3 cm	SONES	4		



PIGTAIL



145° PIGTAIL



155° PIGTAIL



SONES



Recommended guide wire is 0.038" (0.97 mm)

Performa Diagnostic Cardiology Catheters

Catalog Number	French Size	Internal Diameter	Length	Length of Tip	Description	# Side Ports	
MULTIPURPOSE							
7708-30	4	.042"/1.07 mm	80 cm	5.8 cm	MULTI A2	2	bumper tip
7708-40	4	.042"/1.07 mm	80 cm	5.8 cm	MULTI B2	2	bumper tip
7708-10	4	.042"/1.07 mm	100 cm	5.8 cm	MULTI A1	-	bumper tip
7708-20	4	.042"/1.07 mm	100 cm	5.8 cm	MULTI A2	2	bumper tip
7708-50	4	.042"/1.07 mm	100 cm	5.8 cm	MULTI B1	-	bumper tip
7708-60	4	.042"/1.07 mm	100 cm	5.8 cm	MULTI B2	2	bumper tip
7533-53	5	.046"/1.17 mm	80 cm	5.8 cm	MULTI A2	2	
7533-63	5	.046"/1.17 mm	80 cm	5.8 cm	MULTI B2	2	
7533-13	5	.046"/1.17 mm	100 cm	5.8 cm	MULTI A1	-	
7533-23	5	.046"/1.17 mm	100 cm	5.8 cm	MULTI A2	2	
7533-33	5	.046"/1.17 mm	100 cm	5.8 cm	MULTI B1	-	
7533-43	5	.046"/1.17 mm	100 cm	5.8 cm	MULTI B2	2	
7513-A3	6	.054"/1.37 mm	80 cm	5.8 cm	MULTI A1	-	
7513-13	6	.054"/1.37 mm	80 cm	5.8 cm	MULTI A2	2	
7513-D3	6	.054"/1.37 mm	80 cm	5.8 cm	MULTI B1	-	
7513-33	6	.054"/1.37 mm	80 cm	5.8 cm	MULTI B2	2	
7513-B3	6	.054"/1.37 mm	100cm	5.8 cm	MULTI A1	-	
7513-23	6	.054"/1.37 mm	100 cm	5.8 cm	MULTI A2	2	
7513-E3	6	.054"/1.37 mm	100 cm	5.8 cm	MULTI B1	-	
7513-43	6	.054"/1.37 mm	100 cm	5.8 cm	MULTI B2	2	
7513-C3	6	.054"/1.37 mm	125 cm	5.8 cm	MULTI A1	-	
7513-53	6	.054"/1.37 mm	125 cm	5.8 cm	MULTI A2	2	
7513-F3	6	.054"/1.37 mm	125 cm	5.8 cm	MULTI B1	-	
7513-63	6	.054"/1.37 mm	125 cm	5.8 cm	MULTI B2	2	

Recommended guide wire is 0.038" (0.97 mm)

Catalog Number	French Size	Description	Catheters and Individual Catalog Numbers
PERFORMA MULTIPACK™			
7751-10	4	CARDIAC MULTIPACK I	PIG (Braided), JL4.0, JR4.0 (7704-10, 7700-20, 7701-20)
7560-13	5	CARDIAC MULTIPACK I	PIG (Braided), JL4.0, JR4.0 (7529-11, 7521-21, 7523-21)
7560-A3	5	CARDIAC MULTIPACK I	PIG (Braided), JL4.0, PRC (7529-11, 7521-21, 7532-13)
7560-23	5	CARDIAC MULTIPACK II	145° PIG (Braided), JL4.0, JR4.0 (7529-23, 7521-21, 7523-21)
7560-C3	5	CARDIAC MULTIPACK II	145° PIG (Braided), JL4.0, PRC (7529-23, 7521-21, 7532-13)
7560-33	5	CARDIAC MULTIPACK III	155° PIG (Braided), JL4.0, JR4.0 (7529-33, 7521-21, 7523-21)
7560-43	5	CARDIAC MULTIPACK IV	PIG (Non-Braided), JL4.0, JR4.0 (7602-11, 7521-21, 7523-21)
7560-53	5	CARDIAC MULTIPACK V	145° PIG (Non-Braided), JL4.0, JR4.0 (7602-23, 7521-21, 7523-21)
7560-G3	5	CARDIAC MULTIPACK V	PIG (Braided), JL4.0 (high flow), JR4.0 (7529-11, 7601-21, 7523-21)
7560-83	5	CARDIAC MULTIPACK V	155° PIG (Non-Braided), JL3.5, JR4.0 (7602-33, 7521-13, 7523-21)
7660-13	6	CARDIAC MULTIPACK I	PIG (Braided), JL4.0, JR4.0 (7509-11, 7501-21, 7503-21)
7660-A3	6	CARDIAC MULTIPACK I	PIG (Braided), JL4.0, PRC (7509-11, 7501-21, 7512-13)
7660-23	6	CARDIAC MULTIPACK II	145° PIG (Braided), JL4.0, JR4.0 (7509-23, 7501-21, 7503-21)
7660-B3	6	CARDIAC MULTIPACK II	145° PIG (Braided), JL4.0, PRC (7509-23, 7501-21, 7512-13)
7660-33	6	CARDIAC MULTIPACK III	155° PIG (Braided), JL4.0, JR4.0 (7509-33, 7501-21, 7503-21)
7660-43	6	CARDIAC MULTIPACK IV	PIG (Non-Braided), JL4.0, JR4.0 (7569-11, 7501-21, 7503-21)
7660-53	6	CARDIAC MULTIPACK V	145° PIG (Non-Braided), JL4.0, JR4.0 (7569-23, 7501-21, 7503-21)
7660-63	6	CARDIAC MULTIPACK VI	155° PIG (Non-Braided), JL4.0, JR4.0 (7569-33, 7501-21, 7503-21)



CATHETERS & EXTRUSIONS

Performa Diagnostic Cardiology Catheters

Catalog Number	French Size	Diagnostic Catheters 100 cm Pigtail Catheters 110 cm	Sheath (11cm) w/ 0.038 Guide Wire	Guide Wire (150 cm) 3 mm J 0.035" (0.89 mm) Fixed Core	Needle
PERFORMA MULTIPACK PLUS™					
4MPPRO-W-S	4	PRC JL4.0	Straight Pigtail	4F PRO*	Yes
4MPPRO-W-145	4	JR4.0 JL4.0	145° Pigtail	4F PRO*	Yes
5MPSW-S	5	JR 4.0 JL 4.0	Straight Pigtail	5F Prelude	Yes
5MPSW-145	5	JR 4.0 JL 4.0	145° Pigtail	5F Prelude	Yes
5MPS-S	5	JR 4.0 JL 4.0	Straight Pigtail	5F Prelude	No
5MPS-145	5	JR 4.0 JL 4.0	145° Pigtail	5F Prelude	No
5MPSW-SNB	5	JR 4.0 JL 4.0	Straight Pigtail †	5F Prelude	Yes
5MPSW-S-18G	5	JR 4.0 JL 4.0	Straight Pigtail	5F Prelude	Yes 18G
5MPPROW-S	5	JR 4.0 JL 4.0	Straight Pigtail	5F PRO*	Yes
5MPPRO-S	5	JR 4.0 JL 4.0	Straight Pigtail	5F PRO*	No
5MPPRO-145	5	JR 4.0 JL 4.0	145° Pigtail	5F PRO*	No
5MPPRO-W-145	5	JR 4.0 JL 4.0	145° Pigtail	5F PRO* **	Yes
5MPPROW-145	5	JR 4.0 JL 4.0	145° Pigtail	5F PRO*	Yes
6MPSW-S	6	JR 4.0 JL 4.0	Straight Pigtail	6F Prelude	Yes
6MPSW-145	6	JR 4.0 JL 4.0	145° PPigtail	6F Prelude	Yes
6MPS-S	6	JR 4.0 JL 4.0	Straight Pigtail	6F Prelude	No
6MPS-145	6	JR 4.0 JL 4.0	145° Pigtail	6F Prelude	No
6MPSW-S-18G	6	JR 4.0 JL 4.0	Straight Pigtail	6F Prelude	No 18G
6MPPROW-S	6	JR 4.0 JL 4.0	Straight Pigtail	6F Pro*	Yes
6MPPROW-145	6	JR 4.0 JL 4.0	145° Pigtail	6F Pro*	Yes

* Polypropylene (PRO) sheath tubing offers increased resilience and additional support during procedure.


** Contains 5F Prelude PRO™ with 0.035" guide wire.

† non braided



Catalog Number	French Size	Internal Diameter	Length	Length of Tip	Description	# Side Ports
PEDIATRIC JUDKINS RIGHT AND LEFT CATHETERS						
7700-A0	4	0.042"/1.07 mm	70 cm	2.0 cm	JL 1.5	- bumper tip
7700-B0	4	0.042"/1.07 mm	70 cm	2.0 cm	JL 2.0	- bumper tip
7700-C0	4	0.042"/1.07 mm	70 cm	2.0 cm	JL 2.5	- bumper tip
7700-D0	4	0.042"/1.07 mm	70 cm	2.0 cm	JL 3.0	- bumper tip
7701-A0	4	0.042"/1.07 mm	70 cm	2.0 cm	JR 1.5	- bumper tip
7701-B0	4	0.042"/1.07 mm	70 cm	2.0 cm	JR 2.0	- bumper tip
7701-C0	4	0.042"/1.07 mm	70 cm	2.0 cm	JR 2.5	- bumper tip
7701-D0	4	0.042"/1.07 mm	70 cm	2.0 cm	JR 3.0	- bumper tip
7521-L3	5	0.046"/1.17 mm	70 cm	2.0 cm	JL 1.5	- bumper tip
7521-M3	5	0.046"/1.17 mm	70 cm	2.0 cm	JL 2.0	- bumper tip
7521-N3	5	0.046"/1.17 mm	70 cm	2.0 cm	JL 2.5	- bumper tip
7521-P3	5	0.046"/1.17 mm	70 cm	2.0 cm	JL 3.0	- bumper tip
7523-L3	5	0.046"/1.17 mm	70 cm	2.0 cm	JR 1.5	- bumper tip
7523-M3	5	0.046"/1.17 mm	70 cm	2.0 cm	JR 2.0	- bumper tip
7523-N3	5	0.046"/1.17 mm	70 cm	2.0 cm	JR 2.5	- bumper tip
7523-P3	5	0.046"/1.17 mm	70 cm	2.0 cm	JR 3.0	- bumper tip

Recommended guide wire is 0.038" (0.97 mm)



Performa Diagnostic Cardiology Catheters

Catalog Number	French Size	Internal Diameter	Length	Length of Tip	# Side Ports	
PERFORMA STANDARD PEDIATRIC AND INFANT PIGTAIL CATHETERS - (NYLON)						
7756-10	4	0.041"/1.04 mm	65 cm	6.0 cm	5	std non-braided
7756-20	4	0.041"/1.04 mm	80 cm	6.0 cm	5	std non-braided
7757-10	4	0.041"/1.04 mm	65 cm	5.5 cm	4	infant non-braided
7757-20	4	0.041"/1.04 mm	80 cm	5.5 cm	4	infant non-braided

Recommended guide wire is 0.035" (0.89 mm)

INFANT PIG

STANDARD PIG

Catalog Number	French Size	Internal Diameter	Length	Recommended Guide Wire	# Side Ports
PEDIATRIC PIGTAIL CATHETERS - (FEP - FLUORINATED ETHYLENE PROPYLENE)					
1310-21-0043	3	0.027"/.69 mm	40 cm	0.021" (0.53 mm)	4
1310-21-0053	3	0.027"/.69 mm	50 cm	0.021" (0.53 mm)	4
1310-21-0054	4	0.040"/1.02 mm	50 cm	0.021" (0.53 mm)	4
1310-21-0064	4	0.040"/1.02 mm	65 cm	0.021" (0.53 mm)	4
1310-25-0054	4	0.040"/1.02 mm	50 cm	0.025" (0.64 mm)	4
1310-25-0064	4	0.040"/1.02 mm	65 cm	0.025" (0.64 mm)	4
1310-25-0084	4	0.040"/1.02 mm	80 cm	0.025" (0.64 mm)	4
1310-25-0065	5	0.046"/1.17 mm	60 cm	0.025" (0.64 mm)	4
1310-25-0085	5	0.046"/1.17 mm	80 cm	0.025" (0.64 mm)	4
1308-25-0085	5	0.046"/1.17 mm	80 cm	0.025" (0.64 mm)	6
1308-35-0065	5	0.046"/1.17 mm	65 cm	0.035" (0.89 mm)	6
1308-35-0085	5	0.046"/1.17 mm	80 cm	0.035" (0.89 mm)	6
1308-35-0005	5	0.046"/1.17 mm	100 cm	0.035" (0.89 mm)	6
1308-35-0086	6	0.056"/1.42 mm	80 cm	0.035" (0.89 mm)	6
1308-35-0006	6	0.056"/1.42 mm	100 cm	0.035" (0.89 mm)	6
1308-38-0006	6	0.056"/1.42 mm	100 cm	0.038" (0.97 mm)	6
1308-35-0007	7	0.065"/1.65 mm	100 cm	0.035" (0.89 mm)	6
1308-38-0007	7	0.065"/1.65 mm	100 cm	0.038" (0.97 mm)	6

PIGTAIL