

# PNEUMOCATH® / NEO-PNEUMOCATH

Thin and radiopaque drainage in slide-in cover with marker rings, three-way-stopcock with suction connection tube for drainage of air and/or fluid in the pleural cavity, as well as in all emergency situations where the application of a surgical suction drainage might come too late or is quite a risk. Due to its small size and easy handling, the Pneumocath® thoracic drain is particularly valuable in paediatrics. The NEO-PNEUMOCATH® with it enlarged lumen, offers a new and improved technique for drainage in the thoracic cavity. It is also a valuable device for simplifying the Bülau suction method and other trocar methods. A high degree of safety is ensured by compatibility of the material and an Heimlich-valve (PNEUMOVENT®) with fixed LL-mechanism ensures more security. The catheter permits sterile, fast and safe handling without any additional aids.



#### Content

#### Basic equipment consisting of:

- · Suction catheter with protective cover
- Puncture needle and three-way-stopcock with tube adapter

#### Features of the kits according to the declaration

		Stitch length		
REF	PNEUMOCATH	Needle	Catheter Length	Catheter Size
503 001	PNEUMOCATH® for adults	85 mm	500 mm	8 F (2.7 mm)
503 002	PNEUMOCATH® for children	65 mm	400 mm	8 F (2.7 mm)
503 003	PNEUMOCATH® for prenature and newborns	55 mm	300 mm	6 F (2.0 mm)
503 201	PNEUMOCATH® fix-adapter for adults	85 mm	500 mm	8 F (2.7 mm)
Box: 20 pcs.				
503 401	PNEUMOCATH® Kit consisting of • Basic equipment	85 mm	500 mm	8 F (2.7 mm)
	<ul><li>50 ml syringe</li><li>Collecting bag 2.0 ltr. with 90 cm tube</li></ul>			
503 501	PNEUMOCATH® Kit consisting of	85 mm	500 mm	8 F (2.7 mm)
	• 503 401 but with			
Box: 10 pcs.	<ul> <li>Heimlich valve (Pneumovent®)</li> </ul>			
box. To pes.		Stitch length		
REF	NEO-PNEUMOCATH®	Needle	Catheter Length	Catheter Size
503 011	NEO-PNEUMOCATH® for adults	85 mm	500 mm	10 F (3.2 mm)
503 211	NEO-PNEUMOCATH® for adults with fix-adapter	85 mm	500 mm	10 F (3.2 mm)
503 012	NEO-PNEUMOCATH® for children	65 mm	400 mm	10 F (3.2 mm)
Box: 20 pcs.				
503 014	NEO-PNEUMOCATH® emergency model with fix-Adapter in Hardblister	85 mm	500 mm	10 F (3.2 mm)
Box: 10 pcs.				



# **THORACENTESIS / Pleural-Puncture-Kits**

Kits for pleural puncture in case of:

- Pneumothorax
- Pleurodesis
- Aspiration of any effusions (i.e. ascites)

They consist of a puncture cannula with catheter (over the needle).

Besides the straight version, the combined models with self-locking-seal and sideport, are particularly well suited for pleurodesis applications. Special possibility for infusion of drugs with subsequent optional suction. Reliable handling by stopcock and a flexible silicon tube. All models have a lateral eye on the distal end of the catheter.

The models with Veress-Needle offers an additional protection against injuries.







REF	Description
503 021 503 019 503 119	Thoracentesis (14 G – 10 cm) straight version Thoracentesis (14 G – 8.5 cm) for adults with side port Thoracentesis-Ascites-Puncture Kit consisting of • Thoracentesis for adults with sideport (14 G – 8.5 cm) • Syringe 50-60 ml LL
	• Collecting bag 2.0 ltr. LL with 90 cm tube
503 061	Thoracentesis Kit with Veress-Needle (1.8 x 120 mm) and • Collecting bag 2.0 ltr. LL with 90 cm tube • Stopcock • Syringe 50-60 ml LL • 75 cm extension line
503 060	Thoracentesis Kit with Veress-Neddle (2.0 x 120 mm) and Collecting bag 2.0 ltr. LL with 90 cm tube Stopcock Syringe 50-60 ml LL To cm extension line
503 025	Pleural-Puncture-Kit consisting of • Thoracentesis (14 G – 8 cm) with puncture needle, 2-way-stopcock and syringe (5 ml)
Box: 20 pcs.	
503 022	Pleural-Puncture-Kit consisting of  • Thoracentesis (14 G – 10 cm), straight version  • Collecting bag 2.0 ltr. LL with 90 cm tube  • Stopcock  • Syringe 50-60 ml LL  • 75 cm extension line

Box: 10 pcs.

Other kits conforming customer's requests on enquiry.



### 24 Thoracic Drainage - Pleural Puncture - Accessories

# **Accessories**

intra special catheters offers suitable accessories for all thoracic catheters, such as Heimlich-type drainage valves PNEUMOVENT® (with LL-connections or also tube connections), various variable adaptors (with or without stopcock), a special fixation plate to adapt three-way stopcocks, as well as suction tubes with silicone adaptors.







REF	Description
504 001	Pneumovent® with LL-connection
504 002	Pneumovent® with tube connection
503 010	Connection tube to connect the thin drains to bigger systems
503 050	Suction tube 50 cm, male/silicone adaptor
503 180	Suction tube 180 cm, male/silicone adaptor
010 0 <b>45</b>	Collecting bag 1.5 ltr. with 90 cm tube
Box: 20 pcs.	
011 000	Variable adaptor without stopcock
011 001 Box: 50 pcs.	Variable adaptor with stopcock for NEO-PNEUMOCATH® (10 F)

All collecting bags with back valve and bottom valve.

