

M93UB



CVC Insertion Simulator II

Conduct a variety of CVC insertion and management simulations with the assistance of 3 interchangeable pads




Product Supervision

Masahiro Tanabe, M.D., Ph.D., Professor Director General Medical Education Center Chiba University School of Medicine

Kinya Sando, M.D., Ph.D., Professor Director Department of Human Dietetics Graduate School of Human Science Osaka Shoin Women's University

Masanori Hoki, M.D., Ph.D., Professor Department of Human Dietetics Graduate School of Human Science Osaka Shoin Women's University

Joho Tokumine, M.D.Ph.D., Department of Anesthesia, Kawatetsu Chiba Hospital



FEATURES

- 1 | Realistic ultrasound-guided CVC
- 2 | Landmark puncture pad with anatomically correct vein bifurcations simulates mechanical complications including pneumothorax, mislodging and artery puncture
- 3 | Transparent anatomical block for anatomical understanding and guide wire manipulation
- 4 | Introductory ultrasound training block to acquire basics of ultrasound guided venous access
- 5 | Both internal jugular and subclavian (axillary) veins are accessible

SKILLS

- | Ultrasound-guided CVC
- | Landmark guided CVC
- | Ultrasound-guided venous access
- | Prevention of mechanical complications

KEY FEATURES



Transparent anatomical block



Ultrasound training block



Ultrasound puncture pad

DESCRIPTIONS

SPECIFICATIONS

Size: W40 x D20 x H 34 cm
W15.8 x D7.9 x H 13.4 inch

Weight: 2.3kg/5lbs

MATERIALS

Soft resin
Latex free

REPLACEMENT PARTS

11347-119 2 landmark puncture pads
11347-170 2 ultrasound puncture pads
11347-210 2 introductory ultrasound training blocks
11347-320 irrigation bottle
11388-400 10 cotton swab red
11402-010 10 cotton swab blue

SET INCLUDES

1 male upper torso manikin	1 introductory ultrasound training block	1 irrigation bottle
1 landmark puncture pad	10 cotton swab red	2 syringes (2.5ml, 50ml)
1 ultrasound puncture pad	10 cotton swab blue	1 sample needle
1 skin for cannulation training	1 air bulb	
1 transparent anatomical block	2 plastic jars	

Ultrasound scanner and needles are not included.