

# Arterial Seldinger Catheters

## MICROSELD®

Radiopaque arterial catheter for measurement and diagnosis

MICROSELD® is suited for non-traumatic puncture of peripheral arteries for safe and complication-free positioning. It can be applied for any operation where a continuous blood pressure monitoring in anaesthesia and internal medicine is required and vessel damage or bleeding is to be eliminated, especially for arteria radialis, arteria brachialis and arteria femoralis. MICROSELD® is also suited for monitoring of the heart function, circulation and lung function control as well as for sampling of blood gas. The used biocompatible Teflon material avoids deposition of any kind of blood ingredients such as fibrin, erythrocytes and leucocytes. For this reason MICROSELD®-PTFE is particularly suited for long-term-monitoring of the catheter over 48 hours. But also the MICROSELD®-PEBAX has a big flexibility area, an excellent moving fortune and a good chemical resistance. It is produced without softeners.

### Content

- Micro-catheter with flexible wings
- Puncture needle with transparent handle
- Safety guide wire with flexible tip 3 cm

REF	LENGTH X SIZE	FRENCH SIZE	MATERIAL
<b>Box: 20 pcs.</b>			
302 042 / 302 042A	4 cm x 0.7 mm	2	PTFE
302 062 / 302 062A	6 cm x 0.7 mm	2	PTFE
302 063 / 302 063A	6 cm x 1.0 mm	3	PTFE
302 083 / 302 083A	8 cm x 1.0 mm	3	PTFE
302 113 / 302 113A	11 cm x 1.0 mm	3	PTFE
302 153 / 302 153A	15 cm x 1.0 mm	3	PTFE
302 084 / 302 084A	8 cm x 1.3 mm	4	PTFE
302 114 / 302 114A	11 cm x 1.3 mm	4	PTFE
302 154 / 302 154A	15 cm x 1.3 mm	4	PTFE
302 204 / 302 204A	20 cm x 1.3 mm	4	PTFE
302 115	11 cm x 1.7 mm	5	PTFE
302 205	20 cm x 1.7 mm	5	PTFE
305 043 / 305 043A	4 cm x 1.0 mm	3	Pebax
305 063 / 305 063A	6 cm x 1.0 mm	3	Pebax
305 083 / 305 083A	8 cm x 1.0 mm	3	Pebax
305 113 / 305 113A	11 cm x 1.0 mm	3	Pebax
305 114 / 305 114A	11 cm x 1.3 mm	4	Pebax

Other lengths and sizes on enquiry.

