

Endless Love 

To Guard The New Life Safe



# BM-100A

Transcutaneous Jaundice Detector



Stock code [300314]



## Brief Introduction

### ■ Effective, User-friendly, Intelligent transmission

3.0 inch color touch screen, large character indicates and read measurement value clearly and rapidly;

On-screen display of multiple measurement data, and simultaneous storage of multiple patient information, achieving data management to patients.

### ■ Long endurance, environment friendly and durable

Equipped with charging base, enduring built-in chargeable battery, measuring 2000 times after battery is full;

Xenon flash lamp with long lifespan (measurement  $\geq 150,000$  times);

Charging base with the function of integrated light checker, consumable-free, saving cost for check.


### ■ Accurate measurement and quick check

Effective check to identify patient with high-risk jaundice precisely, conform to the newest clinical guidance;

Constant measurement to multiple patients without repeating calibration;

Average measurement: this function permits 1~5 times setting of average measurement to ensure the accuracy on results;

Converting between two kinds of measurement unit mg/dL,  $\mu\text{mol/L}$ .



## Standard Configuration

Consist of main body, base and power adapter.



Ningbo David Medical Device Co.,Ltd

Add: Bldg No. 5 Ningbo Smart Park, #98 Chuangyuan Rd, Hi-tech Zone, Ningbo, China

Tel: 0086-574-8780 0008 or 8780 0007

Fax: 0086-574-8780 3714

E-mail: sales@nbdavid.com

Web.: www.nbdavid.com

P.C.: 315048