

Endless Love 

To Guard The New Life Safe



BM-100C

Transcutaneous Jaundice Detector



Stock code [300314]



Brief Introduction

Effective, User-friendly, Intelligent transmission

3.0 inch color touch screen, large character indicates and read measurement value clearly and rapidly;
On-screen display of multiple measurement data, and simultaneous storage of multiple patient information, achieving data management to patients;
Configured with bar code reader, by scanning QR and bar code, it can check or enter nurse and patient ID quickly and accurately;
Configured with micro USB data communication connector, can be connected with PC to transmit data, record and management are simple.



Long endurance, environment friendly and durable

Equipped with charging base, enduring built-in chargeable battery, measuring 2000 times after battery is full;
Xenon flash lamp with long lifespan (measurement $\geq 150,000$ times);
Charging base with the function of integrated light checker, consumable-free, saving cost for check.


Accurate measurement and quick check

Effective check to identify patient with high-risk jaundice precisely, conform to the newest clinical guidance;
Constant measurement to multiple patients without repeating calibration;
Average measurement: this function permits 1~5 times setting of average measurement to ensure the accuracy on results;
Converting between two kinds of measurement unit mg/dL, $\mu\text{mol/L}$.



Standard Configuration

Consist of main body, base, power adapter, bar code reader and USB data communication connector.



Add: Bldg No. 5 Ningbo Smart Park, #98 Chuangyuan Rd, Hi-tech Zone, Ningbo, China
Tel: 0086-574-8780 0008 or 8780 0007 Fax: 0086-574-8780 3714
E-mail: sales@nbdavid.com Web.: www.nbdavid.com P.C.: 315048