



HyLED 200 Series NEW Examination Light



mindray
healthcare within reach

Mindray Building, Keji 12th Road South,
High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R. China
Tel: +86 755 8188 8998 Fax: +86 755 26582680
E-mail: intl-market@mindray.com www.mindray.com

mindray | healthcare within reach are registered trademarks or trademarks owned by Shenzhen Mindray Bio-medical Electronics Co., LTD.
©2014 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved. Specifications subject to changes without prior notice.
P/N:ENG- HyLED200-210285x4P-20160315

mindray

HyLED 200 Series

LED Examination Light

Mindray examination light inherits the advanced technology from LED surgical light and shows excellent light performance: better light intensity, longer service life and lower power consumption.

It is ideal for patient examination, ICU, outpatient, as well as preparation and recovery areas.

- 70,000 lux illuminance for the best possible visualization
- Up to 10 hours battery(optional) duration in mobile version
- High flexibility and space-saving compact design

Flexible and Simple Handling


The following features can be operated electronically using the control panel in the light head:

- On/off switching
- 5-step brightness control



Technical Specifications

As well as mobile lights, HyLED 200 examination light is also available in ceiling, wall and even pendant mounted versions which free the floor space and facilitate cleaning.



	HyLED200 /HyLED 200 M
Central illuminance (at 1m distance)	70.000 lux
Light field diameter(at 1m distance)	170 mm
Depth of illumination (20%)	1500 mm
Color Temperature	4350 K
Color rendering index (Ra)	96
Color rendering index (R9)	96
Service life of bulb	≥60,000h
Optional battery supply duration	10 hours