

# M99

Product Code: 11351-000

## Arterial Puncture Training Wrist

An essential simulator to increase confidence and competence

Arterial blood sampling and cannulation is a common yet delicate procedure. With the innovative design of the arterial puncture simulator, a true-to-life experience can be simulated to ensure correct techniques and boost proficiency.



### FEATURES

- | Palpable arterial pulse
- | Realistic needle resistance against tissue and artery wall
- | Natural flashback

### Training skills / Applications

Radial arterial puncture / Artery blood sampling /  
Radial line placement

### SET INCLUDES

1 wrist and hand model / 1 circulation pump / 2 sets of  
puncture skins (inner and outer skin) / 4 artery tubes / 1  
silicon tape / 1 jar of coloring powder (red) / 1 spoon / 2  
plastic jars / 2 connection tubes / 1 vinyl sheet / 1 syringe  
/ 1 storage case / 1 instruction manual

World inquires and Ordering

 **KYOTO KAGAKU co.,LTD**  
rw@kyoto@kyotokagaku.co.jp

Head Office and Factories

15 KitaneKoya-cho, Fushimi-ku  
Kyoto, 612-8388, JAPAN

+81-75-605-2510

+81-75-605-2519



All American regions  
**Kyoto Kagaku America Inc.**

info@kkamerica-inc.com

+1-310-325-8860

+1-310-325-8867

3109 Lomita Boulevard,  
Torrance, CA 90505-510,  
USA

<http://www.kyotokagaku.com>