

# Pediatric FAST/Acute Abdomen Phantom

Product Supervision  
Shunsuke Nosaka, MD, PhD  
Director of Radiology  
National Medical Center for Children and Mothers



This world's first pediatric ultrasound torso phantom provides opportunities of hands-on training in ultrasound examination of young patients. Facilitates training in pediatric ultrasound that is a crucial modality, especially for radiosensitive children.

## Features

1. The phantom includes life-size 2-year-old thoracoabdominal organs, a bone structure, free fluid to learn FAST procedures and pathologies that are commonly seen in pediatrics.
2. With this phantom, trainees can acquire skills in basics of pediatric abdominal ultrasound as well as FAST/Acute sonography.

## Training Skills

- FAST procedures:  
Internal hemorrhage at perihepatic, perisplenic, pelvis, and pericardium area
- Sonography for acute pediatric patient:  
Bowel intussusception, appendicitis and biliary dilatation




Pleural Hemorrhage



Hydronephrosis



appendicitis



## Specifications

Set Includes:  
1 ultrasound phantom  
1 storage case  
1 tutorial manual (DVD)

Size:  
W41 x D15 x H15 cm  
16W x 6D x 6H inch

■ pericardial hemorrhage     ⊙ pathologies  
■ pleural hemorrhage     ● typical pediatric conditions  
■ intra-abdominal hemorrhage

\*The small intestine is also embedded.

Kyoto Kagaku America Inc.

3109 Lomita Blvd.  
Torrance, CA 90505-5108

310-325-8860  
 310-325-8867

[www.kkamerica-inc.com](http://www.kkamerica-inc.com)  
[info@kkamerica-inc.com](mailto:info@kkamerica-inc.com)