

Specifications

Compact System

| | | | |
|----------------|--|-------------------|---|
| Width | 411.9 cm/16.2 in | Height | 86.6 cm/3.4 in |
| Depth | 406.6 cm/16.0 in | Weight | 10.57 kg/23.3 lb including internal battery |
| Monitor | 15.6" 1920x1080 LCD panel | Power | 100-140 V, ~50-60 Hz, 250 VA |
| AC adapter | Input: 100-240 VAC, 50/60 HZ Output: 24 VDC @ 250 W | Power consumption | 310 VA |
| System battery | 14.4 VDC. 98 WH | AUX display port | 1920x1080 resolution |
| Connectors | Ethernet and two USB 3.0 connectors | | |

Compact Cart

| | | | |
|-------------------|--|-----------------------------|---|
| Cart base | 488.8 mm x 488.8 mm (19.2" x 19.2") | Vertical adjustment | 820-1001 mm (32.2"-39.3") |
| Wheels | Four 5" locking swivel wheels rear two wheels include steering lock mechanism | Other specification details | Integrated Ethernet connector Integrated keyboard in slide out drawer Probe holders on both sides to accommodate up to four transducers Three storage bins provided, large, small and rear handle tray Two USB connectors |
| Power consumption | 660 VA max depending on system configuration | Standard cart | Input: 100-240 VAC, 50/60 Hz Output: 24 VDC @ 250 W Integrated AC adapter in bottom compartment Integrated keyboard in slide-out drawer |
| Extended cart | Input: 100-240 VAC, 50/60 Hz Output: 24 VDC @ 250 W + 5 V @ 1 W Integrated keyboard in slide-out drawer Three batteries to provide additional | Deluxe cart | Input: 100-240 VAC, 50/60 Hz Output: 24 VDC @ 250 W + 5 V @ 1 W Integrated keyboard in slide-out drawer Multiport adapter provides ports for |

scanning time

attaching up to three
imaging transducers

Three batteries to
provide additional
scanning time

Premium cart

Input: 100-240 VAC,
50/60 Hz

Output: 24 VDC @ 250
W + 5 V @ 1 W

Integrated keyboard in
slide-out drawer

Multiport adapter
provides ports for
attaching up to three
imaging transducers

Small profile B/W video
printer

Three batteries to
provide additional
scanning time
