


# CONTROL CREATES CONFIDENCE

## RadialFLO™

Arterial Catheter

The first step in setting up a closed in-line pressure monitoring and blood sampling system.



# RadialFLO™


Arterial Catheter

Integral FloSwitch®  
to control blood flow



20G PTFE striped  
radiopaque tubing

Printed Arrow to show vent  
plug removal direction  
and level orientation




Silicone Coated  
needle for smooth  
insertion



Notched needle for  
instant blood return

THINK SAFETY AND INFECTION PREVENTION.™



## PRODUCT INFORMATION

- One handed control of blood flow
- Provides peripheral access to the arterial system
- Uses over-the-needle technique for insertion into the radial artery
- Instant blood return for safe arterial access
- Can be used for pressure monitoring and arterial blood sampling
- Includes luer-lock connection for easy tubing attachment
- Sterile, latex-free, non-pyrogenic, disposable catheter
- Length = 45mm (1.77")
- Outer diameter = 1.10mm (0.043")

RadialFlo is designed to be used with Merit's Safedraw® and Meritrans® DTXPlus® closed in-line sampling system and disposable transducers.

## ORDERING INFORMATION

Catalog Number	Description	Flow Rate
682315	20G x 45mm	49mL/min

Packaged individually; 25 units per carton

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.253.1600  
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)  
Amerikalaan 42, 6199 AE Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733

merit.com

Australia  
+61 (0) 3 8373 4943

Austria  
0800 295 374

Belgium  
0800 72 906 (Dutch)  
0800 73 172 (Français)

Brazil  
+55 11 4561 - 0062

Canada  
1.800.364.4370

China  
+86-10 8561 0788

Denmark  
80 88 00 24

Finland  
0800 770 586

France  
0800 91 60 30

Germany  
0800 182 0871

Hong Kong  
+852-2207 6888

India  
+91-80-41223376

Ireland (Republic)  
1800 553 163

Italy  
800 897 005

Luxembourg  
8002 25 22

Netherlands  
0800 022 81 84

Norway  
800 11629

Portugal  
308 801 034

Russia  
+7 495 221 89 02

Singapore  
+65 67536311

South Korea  
+8210 2205 0707

Spain  
+34 911238406

Sweden  
020 792 445

Switzerland  
(Deutsch)  
+41 22 518 02 30

(Français)  
+41 22 518 02 52

(Italiano)  
+41 22 518 00 35

UK  
0800 973 115