



S SERIES

LOW-PRESSURE ASPIRATOR

Stock code [300314]



Brief Introduction

Low pressure electromagnetism pump, no backflow when power is off;

Equipped with the liquid storage bottle;

S-2 Low pressure aspirator: Audible and Visual Alarm function will be activated when the bottle is filled over 3/4 volume with fluid;

S-3 is made by plastic mould with streamlined design: Audible and Visual Alarm function will be activated when the bottle is filled over 3/4 volume with fluid.

Specification

Power requirement: AC220V-230V/50Hz
or AC110-120V/50-60Hz, 15VA

Pressure range: 0-22Kpa

Aspirate speed is $\geq 400\text{ml}/10\text{sec}$ (water) when the negative pressure is 22Kpa

Working noise: $\leq 55\text{dB(A)}$

Capacity of liquid bottle: 1000ml (S-1/S-2/S-3)

Standard Configuration

Main body (including the electromagnetism pump, vacuum meter, filter, pressure adjustable valve), liquid tube, liquid bottle.

Product Packaging

Package: The whole unit is packed into one carton

Dimension: 385mm x 225mm x 280mm

Gross weight: 4.2kg



S-1



S-2



S-3



Ningbo David Medical Device Co.,Ltd

Add: Bldg No. 5 Ningbo Smart Park, #98 Chuangyuan Rd, Hi-tech Zone, Ningbo, China

Tel: 0086-574-8780 0008 or 8780 0007

/Fax: 0086-574-8780 3714

E-mail: sales@nbdavid.com

/Web.: www.nbdavid.com

/P.C.: 315048