



Manufacturer
B. Braun Medical
204, avenue du Maréchal Juin
BP 331
92107 Boulogne Cedex - France

Continence Care Products

Safety with simplicity for great freedom



Introduction



Genitourinary disorders (incontinence & retention) can affect men, women and children of all ages.

In all cases, the management of these disorders is critical to give autonomy to patients and to enable them to have an excellent quality of life.

In this framework and in partnership with healthcare professionals and patients, **B. Braun** has developed a wide innovative product range dedicated to continence care, to simplify on a daily basis what life complicates.

Practice	B. Braun solution
Intermittent self-catheterisation	Pre-lubricated, sterile catheters, ready to use <ul style="list-style-type: none"> • Actreen[®] Lite Cath • Actreen[®] Glys Set • Actreen[®] Glys Cath
Vesical Instillation	Luer Lock catheter, pre-lubricated, sterile, ready to use <ul style="list-style-type: none"> • Actreen[®] Glys Luer Lock
Urine Collection	Male External Catheters self-adhering <ul style="list-style-type: none"> • Urimed[®] Vision • Urimed[®] Vision Short
Urine Collection	Urine bags <ul style="list-style-type: none"> • Urimed[®] bag: 500 ml - 800 ml - 1,5 l - 2 l
Urine sampling and diuresis for paediatric use	Paediatric urine bags <ul style="list-style-type: none"> • Urinocol • Urinocol Premature
Faecal Collection	Drainable pouch with total skin protection <ul style="list-style-type: none"> • Collector I.F.




Table of contents



Catheters

Catheterisation	6
Actreen® Technology	7
Actreen® Lite Cath	8
Actreen® Glys Set	10
Actreen® Glys Cath	11
Intravesical infusion	12
Actreen® Glys Luer Lock	13



Male External Catheters

Male Incontinence	14
Urimed® Vision	15



Urine Bags

Urimed® Leg Bags	18
Urimed® Night and Day Bags	19
Urimed® Accessories	19



Paediatric and premature urine bags

Urinocol® Paediatric	20
Urinocol® Premature	21



Faecal Collector

Collector I.F.	22
----------------	----

Catheters

Catheterisation

What is Urinary Catheterisation?

Urinary Catheterisation is the act of passing a catheter into the bladder to drain urine.

Intermittent Catheterisation

What is it?

It is a technique that may be carried out by individual themselves (intermittent self-catheterisation), or by health care professionals. Once bladder drainage has ceased, the catheter is removed immediately afterwards.

What are the benefits for the patient?

- Less risk of urinary infection
- Avoid flow-back towards kidneys
- Prevent over-distension or atrophy of the bladder
- Improved quality of life and autonomy

When is this procedure indicated?

- Nerve damage from spinal injury or disease
- Before or after a surgical procedure
- Urethral stricture




Actreen® Technology



Ready to use

- Sterile PVC catheter
- No additional water required
- Hydrophilic lubricant: water and glycerine based formula maintained throughout catheterisation.



No Touch system

- Minimize risk of infection
- Catheter has no direct contact with hands and fingers



No effect on sperm mobility



An in vitro comparative study performed by an external expert demonstrates that Actreen® Glyc has no toxicological effect on sperm motility and motion.⁽¹⁾

(1) Effect of short-term exposure to two hydrophilic-coated and one gel prelubricated urinary catheters on sperm vitality, motility and kinematics in vitro (Minerva urol nefrol 2007; 59:115-24-J.Auger, R.Rihaoui, N.François, F.Eustache).

No risk of residue

When the catheter is immersed in synthetic urine, Actreen® Glyc does not present any visually detectable particules in suspension⁽²⁾.



(2) B.Braun test report 7.5.AA.034.3

Intermittent catheterisation

Actreen® Lite Cath

NEW

Simplify every day life



User friendly design

- Easy handling:
Finger grasp eyelets enabling opening with one finger only



- Universal connector



- No additional wrapping
- Handy packaging:
Two ways of opening the box
Easy accessibility to the product





No Touch and Ready to use

No Touch system :

- Maintained sterility during catheterisation



Discreet and sterile catheterisation, wherever you are



- Highly resistant film
- Additional precut opening:
 - Simplified access to catheter

Men ♂ : One 1/3 part from the catheter tip => release first part of the catheter from protective sheath

Women ♀ : Close to the connector => Release completely the catheter



Ordering information

	TIEMANN 45 cm	Size	CH08	CH10	CH12	CH14	CH16	CH18
		Article number ⁽¹⁾	228108	228110	228112	228114	228116	228118
	NELATON 45 cm	Size	CH08	CH10	CH12	CH14	CH16	CH18
		Article number ⁽¹⁾	228208	228210	228212	228214	228216	228218
	NELATON 20 cm	Size	CH06	CH08	CH10	CH12	CH14	CH16
		Article number ⁽¹⁾	228306	228308	228310	228312	228314	228316

(1) For each country a specific letter is added at the end of the article number

Actreen® Glys Set

Set for an active independent life style

Convenient in any circumstances




- Ready and easy to use:
Compact and lightweight
- Easy emptying: Tear-off corner and eyelet
- Transparent and graduated collecting bag:
Quality and volume of urine can be easily monitored
- No risk of leakage:
anti-reflux valve



Single use,
sterile,
prelubricated set



Ordering information

	TIEMANN 45 cm	Size Article number ⁽¹⁾	CH08 226108	CH10 226110	CH12 226112	CH14 226114	CH16 226116	
	NELATON 50 cm	Size Article number ⁽¹⁾	CH08 226208	CH10 226210	CH12 226212	CH14 226214	CH16 226216	CH18 226218
	NELATON 25 cm	Size Article number ⁽¹⁾	CH06 226306	CH08 226308	CH10 226310	CH12 226312	CH14 226314	CH16 226316

(1) For each country a specific letter is added at the end of the article number





Actreen® Glys Cath

Single use, sterile, prelubricated catheter

Safe and reliable system

- Ready to use
- Flexible protective sheath
- Wide range for adults and children






Adapted for hospital and homecare use

Colour coded

- Each CH size is identified by universal colour connector



Ordering information

	TIEMANN	Size	CH08	CH10	CH12	CH14	CH16	CH18
	45 cm	Article number⁽¹⁾	225108	225110	225112	225114	225116	225118
	NELATON	Size	CH08	CH10	CH12	CH14	CH16	CH18
	45 cm	Article number⁽¹⁾	225208	225210	225212	225214	225216	225218
	NELATON	Size	CH06	CH08	CH10	CH12	CH14	CH16
	20 cm	Article number⁽¹⁾	225306	225308	225310	225312	225314	225316

(1) For each country a specific letter is added at the end of the article number

Intravesical infusion

Vesical instillation

What is it?

A technique which consists of giving medication in solution through a catheter into the bladder to treat various bladder disorders or pathologies.

When is this procedure indicated?

- Detrusor hyperreflexia
- Interstitial cystitis
- Superficial bladder tumour



Actreen® Glys Luer Lock





Single use, sterile. Prelubricated catheter with luer lock connector

Reliability and safety of a secure connection

- Leak-proof
- No risk of direct contact: no splashing, no leakage
- Guaranteed volume of instilled solution



Ordering information

	TIEMANN	Size	CH08	CH10	CH12	CH14	CH16
	45 cm	Article number ⁽¹⁾	227408	227410	227412	227414	227416
	NELATON	Size	CH08	CH10	CH12	CH14	CH16
	45 cm	Article number ⁽¹⁾	227208	227210	227212	227214	227216
	NELATON	Size	CH08	CH10	CH12	CH14	CH16
	20 cm	Article number ⁽¹⁾	227308	227310	227312	227314	227316

(1) For each country a specific letter is added at the end of the article number

Male external catheters


Male Incontinence

Male urinary incontinence is a problem that affects millions of men in the world, particularly as they reach the later years of their lives.

This unintentional loss of urine can range from light, occasional dribbles to more frequent, heavy leakage.


A penile sheath or Male External Catheter (MEC) is used in urinary continence problems in men.

- This is an external device that fits over the penis, secured by an integrated adhesive band. A drainage system (Urine bag) is then connected to the sheath.



Male external catheters are worn permanently:

- They must be comfortable, reliable and must not cause any skin damage or irritation.
- For that reason the choice of the material is of utmost importance.





Urimed® Vision



Urimed Vision MEC are made of 100% silicone, the material with the highest water vapor permeability, compared to latex and elastomer*.

LATEX FREE



Biocompatibility: 100% Silicone

- Hypoallergenic
- Breathability: excellent gas exchange, less risk of maceration and irritation
- The moisture transmission of silicone is close to normal skin.



*Internal method according to the NF ISO 2528 standard

Male external catheters

Urimed® Vision

Comfortable, secure and gentle to the skin

Self-adhering Male External Catheter designed for management of male urinary incontinence

- Soft, breathable silicone, promotes healthy skin and greater patient comfort
- Built – in adhesive band allows the ease of application, improved patient confidence and longer wear time
- Reduced risk of infection: prevents from urine residues and maceration
- A kink – proof funnel helps ensure uninterrupted urine flow
- Transparency for complete monitoring of the skin
- Allows the practice of urinary catheterization through the MEC (Male External Catheter)
- Urimed® Vision fits all major types of collection bags





Adapted to all morphologies

Urimed® Vision MEC are available in five diameters to provide patients comfort, proper fit and greater reliability

Urimed® Vision Standard

- Pack of 30 male external catheters



Ordering information

Description	Ø 25 mm	Ø 29 mm	Ø 32 mm	Ø 36 mm	Ø 41 mm
Article number ⁽¹⁾	IH2525	IH2529	IH2532	IH2536	IH2541

Urimed® Vision Short

- Shorter sheath with adhesive area specifically adapted to fit short or retracted penises
- Pack of 30 male external catheters



Ordering information

Description	Ø 25 mm	Ø 29 mm	Ø 32 mm	Ø 36 mm	Ø 41 mm
Article number ⁽¹⁾	IH4525	IH4529	IH4532	IH4536	IH4541

(1) For each country a specific letter is added at the end of the article number

Urine bags

Urimed®

Leg bags

- Fuzzy back foil – avoid skin irritation for better comfort
- Free flow lever tap – easy draining with one hand
- Tube non kicking, large diameter



Urimed Bag Plus 500 ml



- 1 chamber
- Non sterile
- 60 cm tube, cut to fit



Ordering information


Article number ⁽¹⁾	28501
Box	30

Urimed Tribag Plus 800 ml

ALSO AVAILABLE
STERILE



- 3 chambers
- Non sterile
- 60 cm tube, connected connector



Optimal leg fit
Quiet

Ordering information

	Sterile	Non sterile
Article number ⁽¹⁾	28305	28306
Box	10	10

Urimed Bag Plus 1,5 l



- 1 chamber
- Non sterile
- 90 cm tube

Ordering information

Article number ⁽¹⁾	28150
Box	30

(1) For each country a specific letter is added at the end of the article number





Night & Day Bags



Medicare® 2 l

- Non sterile
- 90 cm tube

Cross outlet with stopper

Urimed® bag 2 l



- 120 cm length tube
- Large diameter tube
- Non sterile

Urimed® SP 2 l



- Needle free sample port
- Sterile
- 100 cm tube

Ordering information

	Urimed® bag 2 l	Urimed® SP 2 l	Medicare® 2 l
Article number ⁽¹⁾	28300	28201	28114
Box	30	10	10

Urimed® Accessories

Leg fixation

Urimed® Fix



- Can be worn on the thigh or lower leg

Urimed® Straps



- Adapted for leg bag

Urimed® Strip



- Fixation for tube and catheter

Bed fixation

Urimed® Handle



Ordering information

Urimed® Fix	Small	Medium	Large
Article number ⁽¹⁾	68520	68530	68540
Units / Box	2	2	2
Urimed® Straps	Article number ⁽¹⁾		68550
	Units / Box		1 pair
Urimed® Strip	Special	Small	Large
Article number ⁽¹⁾	28259	28260	28261
Units / Box	10	10	10
Urimed® Handle	Article number ⁽¹⁾		68548
	Units / Box		1 piece

(1) For each country a specific letter is added at the end of the article number

Paediatric urine bags

Urinocol® Paediatric

Baby > 2 500 g

- High skin tolerance: gentle on application and removal
- Adhesive shape adapted to boy and girl anatomy
- Transparent: permanent visual control of urine


PVC FREE - LATEX FREE

Urine sampling for biochemical, cytological or bacteriological examination

Closed bags




- For urine sampling




Collecting for diuresis management

Open bags with a double stage stopper




Sampling with syringe



Connection to urine bag

Ordering information

			Reference	units/box
	Sterile open paediatric bags	Boy	227546	50
		Girl	227566	50
	Sterile closed paediatric bags	Boy	227550	100
		Girl	227560	100



Urinocol® Premature



Bags for premature infants



- 2 sizes: Mini and Standard depending on baby weight
- With extension line: for incubator



"Parachute" shape to maintain the bag more effectively in place.

Ordering information

			Weight	Reference	units/box
	Sterile premature bags	Standard	1200–2500g	227571	50
		Mini	<1200g	227570	50
	Sterile premature bags (with extension line)	Standard	1200–2500g	227573	50
		Mini	<1200g	227572	50

Faecal Collector

Collector I.F[®]

What is faecal incontinence?

Faecal incontinence is the involuntary loss of stools from the anus. This can result in embarrassment, distress and an impaired quality of life for patient. The faecal collector can make life easier and much more comfortable for patients and for the healthcare professionals.

Description

Drainable pouch with a total skin protector specifically fitted to the patient's anatomy to collect faeces:

- 1 adhesive total skin protector
- 1 transparent drainable pouch
- 1 emptying device outlet
- 1 connector

Advantages

- Only one application per day
- Protection from odour and from infections as the pouch is waterproof
- Prevention of irritation and pressure sores by avoiding contact of the patients faeces with the skin
- Time saving: one pouch per day instead of 5 to 10 changes a day



Ordering information

Article number ⁽¹⁾	H04001
-------------------------------	--------

(1) For each country a specific letter is added at the end of the article number



