

EAR Examination Simulator II

Practice patient-focused care during an otoscopy

EAR II is an otoscopy simulator that focuses on patient safety.

With a sensitive alarm system, patient discomfort can be alerted when the otoscope is excessively inserted.



KEY FEATURES



Alarm indicating painful insertion



Quick switching between cases

FEATURES


- 1 | Alarm indicating painful insertion:
 - beep (self learning mode)
 - silent, red light/yellow light (instructor mode)
- 2 | 9 preset cases
- 3 | Quick switching between cases
- 4 | Ear canals come in 2 sizes - normal and stenosis
- 5 | Includes simulated foreign bodies and ear wax
- 6 | Anatomically accurate and soft ears with auricles that can be pulled to straighten the auditory canal
- 7 | Work with an actual otoscope
- 8 | The manikin neck can be tilted
- 9 | Training is possible with both left and right ears

SKILLS

- | Ear examination with a real otoscope
- | Ear wax and foreign body removal

CASES

- | normal
- | serious otitis media (SOM)
- | mucoid otitis media (MOM)
- | chronic otitis media with perforation
- | acute suppurative otitis media (AOM)
- | cholesteatoma
- | tympanosclerosis
- | traumatic perforations
- | cerumen block



normal

cholesteatoma

DESCRIPTIONS

For product codes of the alternate skin tone versions of this model, please refer to www.kyotokagaku.com.

SPECIFICATIONS

Size: 42 x 22 x 39 cm / 16.5 x 8.7 x 15.3 inch Weight: 1.75kg / 3.85lbs

SET INCLUDES

1 manikin head and shoulder	1 types of ear units -1 pair each
1 set of foreign body removal kit	-exam ear (normal)
-simulated ear wax, small objects and a spoon	-exam ear (stenosis)
	-FB ear (normal)
1 instruction manual	-FB ear (stenosis)

REPLACEMENT PARTS

11391-010	ear units and foreign body set
	-4 types of ear units - 1 pair each
	ear exam (normal), ear exam (stenosis), FB ear (normal), FB (stenosis)
	-1 set of foreign body removal kit
	simulated ear wax, small objects and a spoon