

MW32



# NCPR Simulator

## Manikin that facilitates Cardio-Pulmonary Resuscitation training of a newborn baby



### FEATURES

- 1 | Airway management training with observable bilateral/unilateral ventilation
- 2 | Soft and realistic umbilical cord capable of venous and arterial line placement, and blood sample collection
- 3 | Feeding tube placement confirmable through gastric fluid aspiration and auscultation.
- 4 | Chest compression

### SKILLS

#### Basic neonatal care

- | Changing body position
- | Hygiene
- | Palpation of the anterior fontanelle
- | Umbilical catheterization
- | Umbilical cord care
- | Tube feeding (nasal and oral) and gastric fluid aspiration
- | Clamping, tying and cutting of the umbilical cord
- | Suction (nasal, oral)

#### Emergency Simulation

- | Chest compression
- | Endotracheal intubation
- | BVM ventilation

\*The product is not designed for bathing.

### ANATOMY

- | ribs
- | bronchi
- | lungs
- | stomach
- | esophagus
- | umbilical cord

### DESCRIPTIONS

#### SPECIFICATIONS

Size: 50 cm / 19.6 inch Weight: 2.6kg / 5.7 lbs

#### SET INCLUDES

- |   |                    |
|---|--------------------|
| 1 | newborn manikin    |
| 1 | lubricant          |
| 1 | instruction manual |

#### MATERIALS

Soft resin  
Latex free

#### REPLACEMENT PARTS

- |           |                                  |
|-----------|----------------------------------|
| 11400-010 | 1 mask for MW21, MW32            |
| 11400-020 | 1 chest skin for MW21, MW32      |
| 11400-030 | 4 umbilical cords for MW21, MW32 |
| 11229-050 | 1 lubricant                      |



#### RECOMMENDED DEVICES

Feeding tube: 5fr.  
Suction catheter: 5fr.  
Endotracheal tube: 3.0mm ID  
Umbilical catheter: 5fr., 30cm  
Laryngoscope: #0  
BVM: neonatal

\* Instruments are not included.