



PRESCRIPTION

PATIENT ACCESS

PREPARATION

APPLICATION

DISCHARGE
MANAGEMENT

Original Perfusor[®] Line

THE CHOICE IS YOURS

Original Perfusor® Line

Color coded tubes for an individual assignment



POLYETHYLENE
(PE)



NOT MANUFACTURED
WITH DEHP



NOT MANUFACTURED
WITH PVC



2 BAR PRESSURE
RESISTANT

- Option to assign high potential drug by color
- Quick identification of specific drug applications
- Designed to reduce risk of confusion
- Low priming volume
- Excellent infusion start-up characteristics in combination with syringe pumps
- 3 different color coded tubes in blue, magenta, green

Product	Color	Material	Connection	Tube length (cm)	Diameter (mm)	Priming volume (ml)	Units per box	Code no. (REF)
Original Perfusor® Line	Blue	PE (PVC-free)	Luer-Lock	150	1.0 x 2.0	1.27	100	8723020
				200		1.67		8723021
	Magenta	PE (PVC-free)	Luer-Lock	150	1.0 x 2.0	1.27	100	8723025
				200		1.67		8723026
	Green	PE (PVC-free)	Luer-Lock	150	1.0 x 2.0	1.27	100	8723023
				200		1.67		8723024