

AUTOMATED INFUSION SYSTEMS

Original Perfusor® Syringes with enhanced performance




ALL IMPROVEMENTS AT A GLANCE:

- Enhanced gliding characteristics for continuous flow
- Improved start up of the infusion at low infusion rates
- Lower alarm triggering times after system occlusion
- Fits into B. Braun Perfusor® pumps







Comparison of flow constancy of a standard 50 ml syringe and the new Original Perfusor® Syringe 50 ml



DEDICATED PUMP SYRINGE FOR AUTOMATED INFUSION PUMPS

The new Original Perfusor® Syringe portfolio offers the right product for every scope of application with modern infusion pumps. Original Perfusor® Syringes are available without needle, with aspiration cannula, and with cannula and integrated 15 µm particle filter. For the handling of light sensitive drugs the orange colored Original Perfusor® Syringe grants a high UV light protection. The photo-resistant barrel shields medication from wave lengths with up to 520 nm.

ORDERING INFORMATION

Original Perfusor® Syringes Portfolio	Description	Color of barrel	Units per box	Code no. (REF)
	Original Perfusor® Syringe 50 ml	transparent	100	8728844F-06
	Original Perfusor® Syringe 50 ml, Aspiration Cannula	transparent	100	8728810F-06
	Original Perfusor® Syringe 50 ml, Aspiration Cannula, 15 µm Filter	transparent	100	8728852F-06
	Original Perfusor® Syringe 50 ml, UV-Protect, Aspiration Cannula, 15 µm Filter	orange	100	8728861F-06