





ACUMARK SIZING BALLOON

Designed For Precise Measurement Of
Cardiovascular Defects



ACUMARK SIZING BALLOON



Soft and tapered tip helps to smoothly enter the vessel, over a 0.035" guidewire;



The balloon is made of an ultra-thin and compliant TPU membrane to ensure a low crossing profile, minimizing the injury to the vessel;



Three radiopaque marker bands located inside the balloon for clear positioning and assistance measuring the defect.

Ordering Information

Model	Balloon Length (mm)	Maximum Defect Size (mm)	Maximum Inflation Volume (cc)
LT-SZB-24	40	27	30
LT-SZB-34	50	42	90

Model	Shaft Size (Fr)	Shaft Usable Length (mm)	Compatible Guidewire (in)
LT-SZB-24	7F	750	0.035
LT-SZB-34	8F	750	0.035



LIFETECH SCIENTIFIC (SHENZHEN) CO.,LTD.

ADD:8F, LifeTech Scientific Building, No.22, Keji 12th Road South, High-tech Industrial Park, Yuehai Subdistrict, Nanshan District, Shenzhen, 518063, P.R.China.

TEI:+86-755-86026250 **FAX:**+86-755-86026251

EMAIL: marketing@lifetechmed.com

EU authorized representative:

LifeTech Scientific (Europe) Coöperatief U.A
Ollemolenstraat 60, 6416CB Heerlen, the Netherlands

www.lifetechmed.com