



# Ankura™

AUI Stent Graft System

**Stable support  
Precise Apposition**



## PRODUCT SPECIFICATION

Code	Proximal Diameter(D1) (mm)	Distal Diameter(D2) (mm)	Total covered Length(L1) (mm)	Delivery system (Fr.)
AUI2412B160	24	12	160	21
AUI2414B140	24	14	140	21
AUI2612B160	26	12	160	21
AUI2614B140	26	14	140	21
AUI2812B160	28	12	160	21
AUI2814B140	28	14	140	21

## Lifetech Scientific (Shenzhen) co.,Ltd

ADD:Floor 1-5, Cybio Electronic Building, Langshan 2nd Street, North Area of High-tech Park,  
Nanshan District, Shenzhen 518057, P. R. China  
Tel:+86-755-86026250 Fax:86-755-86026251  
Email:marketing@lifetechmed.com  
Web:www.lifetechmed.com

EU authorized representative:  
LifeTech Scientific (Europe) Coöperatief U.A.  
Oliemolenstraat 60, 6416CB Heerlen, the Netherlands

2020-EN-1.0



# TIP CAPTURE MECHANISM



## Proximal mini wave

- enhances apposition to the aortic wall

## Connecting bar

- offers axial support
- avoids stent shortening

## No suture on the main body

- avoids pinhole leakage (Type IV endoleak)

## e-PTFE dual membrane

- excellent biocompatibility, no Type III endoleak

# Ankura™

# AUI STENT GRAFT CASES

