

ORDERING INFORMATION

Table.1 The product family of Fustar™ Steerable Introducer (550/700/800/900)

Curve Type	Specification				ID/mm	OD/mm	Deflectable Angle/°
	550	700	800	900			
Type A (30mm)	SVA5F-550	SVA5F-700	SVA5F-800	SVA5F-900	1.8	2.6	0~180
	SVA6F-550	SVA6F-700	SVA6F-800	SVA6F-900	2.1	2.9	
	SVA7F-550	SVA7F-700	SVA7F-800	SVA7F-900	2.5	3.3	
	SVA8F-550	SVA8F-700	SVA8F-800	SVA8F-900	2.8	3.6	
	SVA9F-550	SVA9F-700	SVA9F-800	SVA9F-900	3.1	3.9	
	SVA10F-550	SVA10F-700	SVA10F-800	SVA10F-900	3.5	4.3	
Type B (50mm)	SVB5F-550	SVB5F-700	SVB5F-800	SVB5F-900	1.8	2.6	0~180
	SVB6F-550	SVB6F-700	SVB6F-800	SVB6F-900	2.1	2.9	
	SVB7F-550	SVB7F-700	SVB7F-800	SVB7F-900	2.5	3.3	
	SVB8F-550	SVB8F-700	SVB8F-800	SVB8F-900	2.8	3.6	
	SVB9F-550	SVB9F-700	SVB9F-800	SVB9F-900	3.1	3.9	
	SVB10F-550	SVB10F-700	SVB10F-800	SVB10F-900	3.5	4.3	
	SVB12F-550	SVB12F-700	SVB12F-800	SVB12F-900	4.1	5.0	0~90
	SVB14F-550	SVB14F-700	SVB14F-800	SVB14F-900	4.8	5.7	

Note: 1F=1 French=1/3 mm



Navigate to Your Destination

Lifetech Scientific (Shenzhen) Co., Ltd.

Add: 8F, LifeTech Scientific Building, No.22, Keji 12th Road South, High-tech Industrial Park, Yuehai Subdistrict, Nanshan District, Shenzhen 518063, P. R. China.
 Tel: +86-755-86026250 Fax: 86-755-86026251
 Email: marketing@lifetechmed.com
 Web: www.lifetechmed.com

EU authorized representative:
 LifeTech Scientific (Europe) Coöperatief U.A
 Olievliedestraat 60, 6416CB Heerlen, the Netherlands


2022-EN-1.0

Outstanding Design


Enable to position and place occluder precisely with deflectable 0-160°



Double ASD



Structural and congenital heart disease



Heart proximal branch vessels

Vessel Protection, Minimize Surgery Risk


Reach target point with coaxial status and minimize risk of mural thrombosis or plaque abscission during advancement

Multiple Potential Applications with Only One Sheath

Apply to various interventional procedures



Abdominal aorta branch vessels




Aortic artery stents

Fustar = Delivery Sheath + Guiding Catheter


Reduce surgery time by avoiding repetitive catheter sizing selection

Safe and Ergonomic Design for Usage

Smooth transition between dilator and sheath can reduce vessel trauma
Quickly and easily access to target



Lower limb arteries



Iliac artery crossover

Distal Silicon Valve Prevents Air Embolism and Minimizes Bleeding