

Nasco *Life/form*® Simulated Burns

**LF00729U
Complete Set**



LF00724U
Chest Burn

LF00722U
Leg Burns



LF00723U
Face Burn

LF00721U
Arm Burns

About the Simulator...

Simulated Burns developed to fit any human patient, patient simulator, or trainer! Each burn is handcrafted with care. The burn that is simulated is a deep partial 2nd degree burn. Made from pliable, easy to cut materials these simulated burns can be used over and over! Use them as they are or add Moulage to create additional realism. Purchase as a complete set (LF00729U) or individually (see above for individual part numbers).

Cleaning and Maintenance

Gently clean simulated burns with soap and water and allow to air dry. Cleaning the burns will remove dirt and dust from the tacky surface of the simulated burn. Store clean, dry simulated burns on tissue paper and return them to the box.

Cautions

- Makeup or artificial blood may stain the simulated burns. The texture of the simulated burn may make blood and makeup difficult to remove.
- Makeup can be purchased separately.
- Actual product may vary slightly from the photo. Nasco reserves the right to change the product, color, materials, or function as needed.
- Prolonged contact with skin may cause irritation.

Grease Paint Makeup Wheels

Each 1-oz. wheel contains a different color palette. Sh. wt. 0.19 lbs.

Injury Shades Wheel

Contains six colors commonly used to create bruises and contusions: warm honey, maroon, dark purple, yellow, forest green, and thunder gray.
LF00760U

Primary Colors Wheel

Includes five colors: white, black, red, yellow, and blue.
LF00761U

Special Effects Shades Wheel

Six colors include: medium honey, green, light cream, ebony, red, and milk chocolate.
LF00762U

Grease Paint Makeup

Made especially for use over simulated wounds and latex. Coats the simulated wounds, latex, and skin sufficiently to allow a natural look. Offers great coverage over a wide variety of skin and simulated skin surfaces. Each jar is ½ oz. Sh. wt. 0.06 lbs.

LF00763U	Capillary Shadow	LF00774U	Forest Green
LF00764U	Blood Blister	LF00775U	Blue
LF00765U	Bliethe Spirit	LF00776U	Dark Purple
LF00766U	Rose Blush	LF00777U	Clotted Red
LF00767U	Mold Green	LF00778U	Corpse Flesh
LF00768U	White	LF00779U	Frankie Gray
LF00769U	Black	LF00780U	Light Cream
LF00770U	Red	LF00781U	Ebony
LF00771U	Green	LF00782U	Milk Chocolate
LF00772U	Deep Yellow	LF00783U	Medium Honey
LF00773U	Yellow	LF00784U	Warm Honey

Creme Makeup Wheel (Auguste)

Colors include white, Auguste light, red, and black creme foundations.

Sh. wt. 0.19 lbs.

LF00789U

Black Soot Powder

Skin-safe specialty powder for realistic burn, charring, and grease effects. 5-oz. jar. Sh. wt. 0.50 lbs.

LF00786U

Ash Dust Powder

Skin-safe specialty powder for realistic distress effects. 5-oz. jar. Sh. wt. 0.50 lbs.

LF00787U

Dirt Powder

Skin-safe powder that simulates real dirt. 5-oz. jar. Sh. wt. 0.50 lbs.

LF00788U

Other Available *Life/form*® Simulators

- LF00720U** Ultra Nursing Wound Simulation Kit
- LF00725U** Simulated Burn Rolls
- LF00793U** Basic Nursing Wound Simulation Kit
- LF00794U** Advanced Nursing Wound Simulation Kit
- LF00795U** Burn Simulation Kit

LF00795U



LF00794U



901 Janesville Avenue, P.O. Box 901
Fort Atkinson, Wisconsin 53538-0901

eNasco.com
1-800-558-9595

E-mail: lifeform@eNasco.com
COPYRIGHT © NASCO 2012 PRINTED IN U.S.A. NP064-14