




Negative-Pressure Wound Therapy

LF00948U

INSTRUCTION CARD



5-Year Warranty



Negative-Pressure Wound Therapy (NPWT) is a therapeutic technique using a vacuum dressing to promote healing in acute or chronic wounds and enhance healing of first and second degree burns. The therapy involves the controlled application of subatmospheric pressure to the wound, using sealed wound dressing and connected to a pump.

CAUTIONS

Note: The glue from the pump will create tackiness on the trainer skin material — just like the tackiness sticks on the real skin of a person. Solvents or corrosive materials will damage the simulator. Never place simulator on newsprint or any kind of printed paper or plastic; these materials will transfer indelible stains. Ballpoint pen ink will also make an indelible stain.

AVAILABLE SUPPLIES

LF09919U Nasco Cleaner

Actual product may vary slightly from photo. Nasco reserves the right to change product color, materials, supplies, or function as needed.



ONLINE: NascoHealthcare.com
EMAIL: lifeform@enasco.com
PHONE: 1.800.431.4310
ADDRESS: 901 Janesville Ave.
Fort Atkinson, WI 53538



NEGATIVE-PRESSURE WOUND THERAPY — LF00948U

Negative-Pressure Wound Therapy (NPWT) is a therapeutic technique using a vacuum dressing to promote healing in acute or chronic wounds and enhance healing of first and second degree burns. The therapy involves the controlled application of subatmospheric pressure to the wound, using sealed wound dressing and connected to a pump.

Cautions:

Note: The glue from the pump will create tackiness on the trainer skin material — just like the tackiness sticks on the real skin of a person.

Solvents or corrosive materials will damage the simulator. Never place simulator on newsprint or any kind of printed paper or plastic. These materials will transfer indelible stains. Ball-point pen ink will also make an indelible stain.

Available Supplies:

LF09919U Nasco Cleaner



901 Janesville Avenue, P.O. Box 901
Fort Atkinson, Wisconsin 53538-0901

eNasco.com
1-800-558-9595

E-mail: lifeform@eNasco.com

COPYRIGHT © NASCO 2012 PRINTED IN U.S.A. NP042-12