


SeQure™ Snare System



Multiple Use, Precise Control

FEATURES & BENEFITS

- **Titanium Nitride Coated Tungsten Loop**

Increased biocompatibility
Enhanced fracture resistance
Excellent radiopacity
Excellent shape supporting and distal shaft control ability



90°angle between Loop and Shaft

- **Flexible Nitinol shaft provides Kink Resistance and 1:1 Torque control**
- **90°angle between Loop and Shaft makes manipulation and retrieval more simple**
- **Radiopaque marker at distal tip of catheter provides excellent Fluoroscopic visualization**

ORDERING INFORMATION

- **Snare System-General**

Code	Loop Diameter	Shaft Length	Catheter Size	Catheter Length
	mm	cm	Fr	cm
LT- SG- 05	5	120	4	100
LT- SG- 10	10	120	4	100
LT- SG- 15	15	120	6	100
LT- SG- 20	20	120	6	100
LT- SG- 25	25	120	6	100
LT- SG- 30	30	120	6	100
LT- SG- 35	35	120	6	100

- **Snare System-Petite**

Code	Loop Diameter	Shaft Length	Catheter Size	Catheter Length
	mm	cm	Fr	cm
LT- SP- 10	10	65	4	50
LT- SP- 25	25	65	6	50

- **Snare System-Micro**

Code	Loop Diameter	Shaft Length	Catheter Size	Catheter Length
	mm	cm	Fr	cm
LT- SM- 02	2	200	3	175
LT- SM- 04	4	200	3	175
LT- SM- 07	7	200	3	175

For further inquiries, please contact your local sales representatives.

Lifetech Scientific (Shenzhen) Co., Ltd.

Add: 8F, LifeTech Scientific Building, No.22, Keji 12th Road South, High-tech Industrial Park, Yuehai Subdistrict, Nanshan District, Shenzhen 518063, P. R. China.

Tel: +86-755-86026250 Fax: 86-755-86026251

Email: marketing@lifetechmed.com

Web: www.lifetechmed.com

EU authorized representative:
Life Tech Scientific (Europe) Coöperatief U.A.
Oliemolenstraat 60, 6416CB Heerlen, the Netherlands