


# Surpass™

SUPER STIFF GUIDE WIRE



CROSSING BY BEST SUPPORT

AN ACCESSORY FOR PERIPHERAL INTERVENTION



# Surpass™

## SUPER STIFF GUIDE WIRE

The Super Stiff Guidewire is intended for complex diagnostic and interventional procedures where increased body, flexibility, and low surface friction of the wire guide is needed, and it is intended for use by physicians trained and experienced in diagnostic and interventional techniques.

### Easy to use


- A radiopaque gold inner coil allows precise tip visibility and location during wire-guide manipulation.

### Low friction surface

- Surpass™ providing maximum support and low friction surface with PTFE coating technology during complex diagnostic and interventional procedures.

### Soft and shapable tip

- Super long tip design stiff shaft and smooth transition at the junction to reduce transfer resistance and vascular damage.



Specification	Guidewire OD (inch)	Guidewire Length (mm)	Inner Coil Length (mm)
SSGW-35-260	0.035	2600	40

### Lifetech Scientific (Shenzhen) co.,Ltd.

Add: Floor 1-5, Cybio Electronic Building, Langshan 2nd Street, North Area of High-tech Park, Nanshan District, Shenzhen 518057, P. R. China

Tel: +86-755-86026250 Fax: +86-755-86026251

Email: [marketing@lifetechmed.com](mailto:marketing@lifetechmed.com)

Web: [www.lifetechmed.com](http://www.lifetechmed.com)

EU authorized representative:  
LifeTech Scientific (Europe) Coöperatief U.A.  
Oliemolenstraat 60, 6416CB Heerlen, the Netherlands