




## ISO80369-6 DISPOSABLE SPINAL NEEDLE "K-3 POINT"

Features :

- Smooth but moderate insertion resistance of the needle makes it possible to identify easily the tissue site through where the needle is inserted.
- Dr JAPAN K-3 point spinal needle is made of thin wall stainless steel for maximum lumen.
- Exact fit of the stylet edge point with the cannula allows close contact of both components and smooth movement of the stylet within the cannula.



Item Number	SIZE	COLOR	GUIDE NEEDLE	Item Number	SIZE	COLOR	GUIDE NEEDLE
50250	18G×90mm	PINK	—	50260	25G×90mm	ORANGE	21G×30mm
50251	19G×90mm	CREAM	—	50261	25G×90mm	ORANGE	21G×38mm
50252	20G×90mm	YELLOW	—	50262	26G×90mm	BROWN	22G×30mm
50253	21G×90mm	GREEN	—	50263	26G×90mm	BROWN	22G×38mm
50254	22G×90mm	BLACK	—	50264	27G×90mm	GRAY	22G×30mm
50255	23G×90mm	BLUE	—	50265	27G×90mm	GRAY	22G×38mm
50256	24G×90mm	PURPLE	—				
50257	25G×90mm	ORANGE	—				
50258	26G×90mm	BROWN	—				
50259	27G×90mm	GRAY	—				

**Dr. Japan Co., Ltd.**

1-1, Kagurazaka, Shinjuku-ku, Tokyo, Japan  
<http://www.drjapan-jp.com/>


SNJ-03 Rev.01  
2019.MAR



## ISO80369-6 DISPOSABLE SPINAL NEEDLE "PENCIL POINT"

Features :

- This tip configuration minimizes the risk of post-dural puncture headache.
- The lateral hole is asymmetrically shaped.  
Soft transition towards the tip facilitates smooth passage; slight ridge in the proximal region enables a "dural clock" to be well perceived penetrating.
- The distance from the tip to the proximal end of the lateral opening, which is equal to the segment segment to be introduced durally, is kept as short as possible.



Item Number	SIZE	COLOR	GUIDE NEEDLE	Item Number	SIZE	COLOR	GUIDE NEEDLE
50266	21G×90mm	GREEN	—	50273	25G×90mm	ORANGE	21G×38mm
50267	22G×90mm	BLACK	—	50274	25G×103mm	ORANGE	21G×38mm
50268	23G×90mm	BLUE	—	50275	25G×125mm	BROWN	21G×38mm
50269	24G×90mm	PURPLE	—	50276	26G×90mm	BROWN	22G×38mm
50270	25G×90mm	ORANGE	—	50277	26G×90mm	BROWN	22G×30mm
50271	26G×90mm	BROWN	—	50278	26G×103mm	BROWN	22G×38mm
50272	27G×90mm	GRAY	—	50279	26G×125mm	BROWN	22G×38mm
				50280	27G×90mm	GRAY	22G×30mm
				50281	27G×90mm	GRAY	22G×38mm
				50282	27G×103mm	GRAY	22G×38mm
				50283	27G×125mm	GRAY	22G×38mm

**Dr. Japan Co., Ltd.**

1-1, Kagurazaka, Shinjuku-ku, Tokyo, Japan  
<http://www.drjapan-jp.com/>

SNJ-04 Rev.01  
2019.MAR