

# DANi Teen & Young Adult Advanced Training Solution

Product Name & Overview	<b>DANi Teen &amp; Young Adult Advanced Training Solution</b>
Addresses Training Needs/Curricula	<ul style="list-style-type: none"> <li>• Basic Life Support (BLS)</li> <li>• Advanced Cardiac Life Support (ACLS)</li> <li>• Certification in Sexual Abuse Trauma (CSAT)</li> <li>• Sexual Assault Nurse Examiner (SANE)</li> <li>• Sexual Assault Nurse Examiner-Pediatric/Adolescent (SANE-P)</li> </ul>
Configuration	<ul style="list-style-type: none"> <li>• Configurable to fit unique training needs and budgets</li> <li>• A range of accessories available to maximize versatility</li> <li>• Provides simulation scenarios not previously available in this patient group</li> </ul>
Consumables	<p><b>IV arm:</b></p> <ul style="list-style-type: none"> <li>• Skins</li> <li>• Veins kit</li> <li>• IM pads</li> <li>• Fluid supply bag</li> <li>• Simulated blood</li> </ul> <p><b>Urinary catheterization:</b></p> <ul style="list-style-type: none"> <li>• Male genitalia foreskins</li> <li>• Foley catheter</li> <li>• Urinary catheterization bladders</li> <li>• Simulated urine</li> </ul> <p><b>Additional consumables:</b></p> <ul style="list-style-type: none"> <li>• Lung bags</li> <li>• IM pads for shoulder, thighs and buttocks</li> <li>• Stomach/stoma/gastronomy bags</li> <li>• Simulated blood</li> <li>• Lubricant</li> <li>• Stomach bags</li> <li>• Enema reservoir</li> </ul>
Product Code/SKUs	<p><b>DANi Teen &amp; Young Adult Advanced Patient Care Trainer</b> LF04500 (Light skin tone) LF04500M (Medium skin tone)</p> <p><b>DANi Teen &amp; Young Adult Advanced Patient Skills Trainer</b> LF04502 (Light skin tone) LF04502M (Medium skin tone)</p> <p><b>DANi Teen &amp; Young Adult Advanced Patient Simulator</b> LF04504 (Light skin tone) LF04504M (Medium skin tone) Skillqube hardware configurations varies.</p>
Dimensions/Weight	<p>Dimensions: 22.25" H x 44.5" L x 25.5" W Ship Weight: 80 lbs.</p>



Be  
**READY**



**Nasco Healthcare Inc.**  
16 Simulaid Drive, Saugerties, NY, 12477 USA  
Phone: 1-833-NASCOHC (627-2642) Email: info@nascohealthcare.com  
www.nascohealthcare.com



# Be **READY** to make a difference in a young person's life.

Evolving from our full body KERi™ manikin, DANi is an advanced training solution designed to support training on the fundamentals of nursing and EMS procedures for issues relating to drug abuse, sexual assault, mental health, or self-harm.

DANi accurately reproduces scenarios and challenging dynamics for learners to administer critical care for this important patient demographic.

A realistic, hands-on learning experience ideally suited for classroom, in-situ, or hospital instruction.

**Each version is affordable and durable with upgradable items.**

### Adult Advanced Patient Care Trainer

General patient care, patient positioning, bandaging and wound dressing, oxygen administration, injection site and sexual assault, drug use and self-harm identification.

### Adult Advanced Patient Skills Trainer

Respiratory support, injection/IV with IM pads, intubation, CPR with quality performance metrics, non-invasive blood pressure, enema administration (female), catheterization, vaginal, pelvic and prostate examinations.

### Adult Advanced Patient Simulator

Comprehensive and complete simulation platform including patient monitoring capability and auscultation (heart and lung sounds).



## DANi Teen & Young Adult Advanced Training Solution

Accurate physical assessments (body map) for signs of physical trauma are essential for legal purposes. DANi allows learners to practice some forensic evidence gathering from blood, semen samples, and more.

### Trauma Simulation<sup>1-3</sup>

- Physical assessment
- Airway, breathing, circulation
- Body trauma scan
- Wound care
- Self-harm/drug use moulage
- Asymmetric pupils



### Realistic Features<sup>1-3</sup>

- Pupils
- Airway
- Acne
- Smooth skin
- Tattoo sleeve
- Piercings & earrings
- Dyed and black wigs
- Foot wounds

### Female & Male Genitalia Inserts<sup>1-3</sup>

Includes female sexual assault genitalia insert

### Patient Positioning<sup>1-3</sup>

Superior range of motion, seating, laying, finger and toe manipulation



### Nursing Care<sup>1-3</sup>

Grooming, daily living assistance, eye/ear/hair care, catheter, injection, oxygen, blood pressure, pericare, ostomy care, urinary catheter care, oxygen administration

### Respiratory Support<sup>2,3</sup>

Manual ventilation with chest rise, respiratory support, intubation, ventilation, tracheostomy care & suctioning

### Auscultation<sup>3</sup>

- 6 anterior heart sites with 12 heart conditions
- 5 anterior, 6 upper posterior, 4 lower posterior, 2 midaxillary lung sites with 12 lung conditions

### Basic Life Support (BLS) & Advanced Cardiac Life Support (ACLS)<sup>2,3</sup>

CPR with quality performance metrics, carotid pulse, intubation

### Advanced Skills<sup>2,3</sup>

Gastrostomy, nasogastric procedures, airway/trach suctioning, enema, pap smear, pelvic exam, intravenous injections



### SKILLQUBE® Patient Vital Signs & Monitoring<sup>3</sup>

- Patient Monitor/Defibrillator
- Instructor controlled scenario editor with library
- Auscultation, glucometer, thermometer, checklists, operating procedures
- Realistic laboratory values
- Optional CPR feedback with real-time evaluation (add-on)

### REFERENCE GUIDE:

1. Available on Advanced Patient Care Trainer
2. Available on Advanced Patient Skills Trainer
3. Available on Advanced Patient Simulator

This manikin is intended to be used as an instructional aide for teaching purposes. Actual product may vary slightly from photo. Nasco Healthcare reserves the right to change product color, materials, supplies, or function as needed.